## FLASHPAK® STERILIZATION CONTAINER SYSTEM

### INSTRUCTIONS FOR USE

#### **Description**

The FlashPak Sterilization Container System consists of a family of rigid reusable containers that provide an effective sterilization packaging method for surgical instruments requiring immediate-use (flash) steam sterilization. Each FlashPak system is comprised of an upper lid that seals by means of a silicone gasket and latches to a lower base creating a totally enclosed environment. A stainless steel wire basket is utilized inside the container to facilitate handling of sterilized items. The lid and base incorporate pressure actuated valves that open to allow sterilant to enter the container during the pressurization portion of the sterilization cycle and then close to seal the container once the cycle is over so the container and its contents can be removed from the sterilizer and immediately transported to the point of use without the risk of recontamination.

#### **Intended Use**

FlashPak is a reusable rigid container system to be used during flash sterilization by hospitals and healthcare facilities. It is intended to enable flash sterilization of the enclosed devices and prevent recontamination during immediate transport to the point of use. The container is compatible with gravity-displacement flash sterilization using a 10 minute cycle at 132° C and with pre-vacuum flash sterilization using a 4 minute cycle at 132° C. FlashPak is recommended for surface sterilization of stainless steel instruments and for lumens with the following limits: gravity-displacement (5.5mm inner diameter or larger and up to 184mm in length), pre-vacuum (1 mm inner diameter or larger and up to 203mm in length).

### **WARNINGS AND PRECAUTIONS**

#### **Warnings**

- Read these instructions completely before using the device.
- FlashPak is designed and intended for immediate-use (flash) steam sterilization only.
- Care should be taken when removing FlashPak after sterilization because it will be hot.

#### **Precautions**

- This product is packaged and shipped non-sterile. It must be cleaned with a neutral pH detergent prior to use.
- Do not use a FlashPak container that is chipped, cracked or has a damaged lid gasket, valve seal or vent button.
- Do not stack the FlashPak containers in the sterilizer as this may inhibit sterilant access.
- FlashPak is recommended for use with an approved multi-parameter chemical integrator (CI). Always verify proper multi-parameter chemical integrator (CI) exposure after every sterilization cycle.
- After sterilization the valves must seat against silicone seals to prevent possible recontamination. Always visually verify full seal contact after every sterilization cycle.
- Immediate-use (flash) steam sterilization does not include dry times so condensation in the container is expected. Care should be taken when transporting the container and when removing the container lid and contents.
- Unless otherwise recommended by the device manufacturer's instructions, do not use FlashPak for sterilization of batteries, motorized, electric or electronic devices.
- The FlashPak valves and pressure balancing element (i.e. metal microfilter) must be periodically replaced in order to ensure continued effectiveness. Replacement of parts should only be performed by trained personnel.

### **CLEANING**

#### **Manual Cleaning**

1. Wash all surfaces (inside and outside) using a neutral 6.0 - 8.0 pH detergent solution. FlashPak may be immersed for cleaning. Soft sponges or cloth wipes are recommended for cleaning.
2. Rinse thoroughly with warm water.
3. Dry the FlashPak thoroughly using lint free absorbent towels or clean pressurized air.

#### **Automated Cleaning**

1. FlashPak may be cleaned with automated washing equipment using a low-sudsing neutral 6.0 - 8.0 pH detergent solution. Follow the washer manufacturers recommended directions for this activity and ensure the container is dry before using or storing it.

## MAINTENANCE

### Daily Activities

1. Clean FlashPak before first use when visibly soiled, and daily according to the recommended directions above. FlashPak may be washed manually or in automated washers.
2. Visually inspect the entire container with each use and daily to assure there are no chips or cracks. Visually inspect the silicone lid gasket, silicone valve seals in the lid and base and the silicone lid vent to assure there are no cuts, tears or other damage.

Valve seals under mounting plate (not visible)  
See Picture for Sterilization Step 1  
Lid gasket

**CAUTION:** Do not use a FlashPak that is chipped, cracked or has damaged gaskets, seals or lid vent button. Remove the affected unit from service and replace damaged or worn parts before re-using.

Lid vent

### Periodic Activities

- Replace both pressure actuated valves, in the lid and base, at least once per year.
- For FlashPak containers used exclusively for gravity-displacement sterilization replace the valves after approximately 1500 uses.
- For FlashPak containers used exclusively for dynamic-air-removal (pre-vacuum) sterilization replace the valves after approximately 500 uses.
- Replace the lid vent button if it becomes damaged.
- Replace the pressure balancing element (or metal microfilter) at least once per year or when it becomes discolored. The microfilter allows the internal container pressure to reach equilibrium with the outside ambient pressure while the valves are closed. Over time it will become clogged with particles from the steam and will appear discolored (e.g. rust spots). Quality of the steam will affect the element life.

**CAUTION:** Replacement of parts should only be performed by trained personnel. The local FlashPak sales representative can provide in-service training to hospital personnel or perform a variety of maintenance services or non-warranty repairs.

### Sterilization Process

1. Visually verify that the valve plate fully engages the silicone valve seals around its entire circumference of both the lid and base valve assemblies.
2. Read and follow the device manufacturer's Instructions for Use (IFU) regarding immediate-use (flash) steam sterilization. Clean and decontaminate instruments and devices per the manufacturer's instructions prior to sterilization. Disassemble instruments and devices as required per the manufacturer's instructions prior to sterilization.
3. Do not attempt to sterilize lumens of smaller diameter, longer length or in greater numbers than shown in the table.


| Lumen Size Recommendations |        |     |            |        |     |
|----------------------------|--------|-----|------------|--------|-----|
| Gravity                    |        |     | Pre-Vacuum |        |     |
| ID                         | Length | Qty | ID         | Length | Qty |
| 5.5 mm                     | 184 mm | 3   | 1 mm       | 203 mm | 5   |

4. See chart below for recommended instrument weight per FlashPak model. Recommended weights are based upon limiting the total loaded container weight to approximately 25 lbs. (11.4 kg) per AAMI ST77:2006/(R)2010, Containment Devices for Reusable Medical Devices. Instruments should be uniformly distributed across the wire basket. Ensure that instruments are properly loaded so that they do not contact the valve or interfere with closing the lid. Hinged devices should be in an open position when placed in the FlashPak.

| FlashPak Model                                                                                                                                                       | Maximum Instrument Load Weight Recommendation |                     |
|----------------------------------------------------------------------------------------------------------------------------------------------------------------------|-----------------------------------------------|---------------------|
|                                                                                                                                                                      | Gravity                                       | Pre-Vacuum          |
| 9020<br>w/ 9020-08 Basket                                                                                                                                            | 3 pounds (1.36 Kg)                            | 3 pounds (1.36 Kg)  |
| 9030<br>w/ 9030-08 Basket                                                                                                                                            | 10 pounds (4.55 Kg)                           | 10 pounds (4.55 Kg) |
| 9040<br>w/ 9040-08 Basket                                                                                                                                            | 14 pounds (6.35 Kg)                           | 14 pounds (6.35 Kg) |
| 9050<br>w/ 9050-08 Basket                                                                                                                                            | 16 pounds (7.27 Kg)                           | 16 pounds (7.27 Kg) |
| Note: Symmetry Surgical has validated the FlashPak model 9050 with a combined instrument and container weight of 35 pounds for Pre-vacuum and 28 pounds for Gravity. |                                               |                     |


5. Place the wire basket containing the instruments inside FlashPak. Instruments may be placed directly into the FlashPak without using the wire basket but care must be taken to ensure the instruments do not contact the valve plate. Contact with the valve plate may interfere with or damage the valve and sterilization efficacy may be affected.

6. Place an approved multi-parameter chemical integrator (CI) in the center of the basket for process monitoring. Additional multi-parameter chemical integrators (CI) may be added as required by hospital procedures and recommended practices.

7. Place the lid on top of the base and secure the four latches.

8. Routine biological testing should be performed in accordance with hospital policy and procedures.
9. Place FlashPak in sterilizer chamber ensuring that its feet sit level and flat on the rack.

**CAUTION: Do not stack on top of another container or another FlashPak.**

10. Select the appropriate sterilization parameters from the table below:


| Sterilization Mode / Device Type                        | Minimum Cycle Time (minutes) | Minimum Temperature |
|---------------------------------------------------------|------------------------------|---------------------|
| Gravity-displacement /<br>(non-porous & porous devices) | 10                           | 132°C (270°F)       |
| Pre-Vacuum /<br>(non-porous & porous devices)           | 4                            | 132°C (270°F)       |

**CAUTION: Read and follow the device manufacturer's immediate-use (flash) steam sterilization instructions for specialty items such as complex devices, power equipment, and batteries before sterilizing them in FlashPak. Failure to follow the manufacturer's instructions may result in irreparable damage to the device or in the device not being sterilized.**

11. Initiate the sterilization cycle.

12. After the sterilization cycle ends, open the sterilizer door and remove FlashPak.

**WARNING: The container will be hot. FlashPak Gloves 9000-14 or towels may be necessary to prevent minor burns. Do not handle or carry FlashPak by the lid.**

13. Transfer the container to the FlashPak Delivery Cart Product Code 9111 or other suitable cart for transport to the point of use. As a result of its closed system design and the raised integral feet FlashPak may be placed on non-sterile surfaces.

**CAUTION: Condensate may have formed inside the FlashPak so use care when transporting.**


**CAUTION: Never touch or probe the valve plate after sterilization as this may breach the sterile barrier and render the contents non-sterile.**

**CAUTION: Never touch or probe the pressure relief button prior to opening the case at the point of use as this may breach the sterile barrier and render the contents non-sterile.**


14. Place FlashPak on a table near the point of use and visually verify that the valve plate in the lid fully engages the silicone seal around its entire circumference. See picture for sterilization step 1 lid gasket.

**CAUTION: If the valve appears to be open or is not fully engaged with the silicone seal assume the instruments are contaminated and reprocess accordingly.**

15. Unlatch the lid. Depress the pressure relief button and remove the lid.

16. Grasping the lid on diagonally opposite corners while lifting may aid lid removal.


17. Verify that the multi-parameter chemical integrator(CI) indicates acceptable exposure. As the terminology used varies by CI brand, follow the CI instructions to verify that sufficient exposure for sterilization has occurred.

**CAUTION: Do not use the instruments if the multi-parameter chemical integrator (CI) does not verify acceptable exposure. Incomplete change in the multi-parameter integrator (CI) indicates that the sterilization process was incomplete or interrupted. Troubleshoot the process according to established procedures and after resolving the problem, repeat these instructions beginning at step 6.**

18. Aseptically remove the basket and contents in accordance with hospital policy and procedures.

19. Visually verify that the base valve plate fully engages the silicone seal around its entire circumference before the basket or its contents are placed on the sterile field. See picture for sterilization step 1 lid gasket.


**CAUTION: If the valve appears to be open or is not fully engaged with the silicone seal assume the instruments are contaminated and reprocess accordingly. Do not touch the valve during this examination.**

20. Residual water condensate is typical after immediate-use (flash) steam sterilization cycles and may be present in the bottom of the container. Remove any liquid water from the FlashPak and dry the container before using it again.


### **Valve Replacement Instructions:**

1. Examine packaging containing replacement valves. Packages damaged during shipping should not be accepted and a claim should be filed with the carrier. Remove the replacement valve from its packaging and examine it for any obvious signs of damage. Valves with obvious damage such as breakage, holes or crushed, bent or distorted parts should not be used.
2. Remove the existing valves from FlashPak base and lid by removing the four (4) retaining nuts for each valve. **Note: Both valves should be replaced at the same time.**
3. If replacing an existing FlashPak Evolution valve with another Evolution valve, skip to step 6.
4. If replacing a legacy diaphragm valve (#9020-09) with an Evolution Valve Conversion Kit (#9090) then the legacy mounting hardware must be discarded. The legacy mounting hardware to discard includes the sheet metal FlashPak Valve Mount bracket, the four (4) plastic spacer nuts, the four (4) plastic valve retaining nuts and the four (4) stainless steel carriage bolts. **Keep the perforated top plate and silicone valve gasket.**

**Perforated Top Plate**

**Silicone Valve Gasket**

5. Use the FlashPak Valve Conversion Kit hardware to install the new Evolution valve. The conversion kit hardware is required to install the Evolution valve in all FlashPak models (Le. 9020, 9030, 9040 and 9050).
  - a. Insert the four (4) stainless steel carriage bolts (#9090-93) from the kit through the perforated top plate first, then the silicone valve gasket and finally through the lid or base. **Be sure that the silicone valve gasket flange is directed toward the inside of the lid or base.**
  - b. Thread and firmly hand-tighten four plastic retaining nuts (#9090-92) onto the carriage bolts; do not use any tools to tighten.
  - c. The new valve is now ready for assembly. Go to step 6 to complete the valve replacement.
6. Align the four holes in the replacement valve with the four stainless steel carriage bolts and push into place. Thread and hand-tighten the four (4) plastic retaining nuts.

## Replacement Parts:

Numerous replacement parts are available to maximize FlashPak in-service time. Refer to the replacement parts chart below. To use the chart, first locate the item that requires replacement and then the model number of the FlashPak itself. The catalog number at the intersection of the appropriate row and column is that of the required part.

Please contact your local sales representative for assistance selecting replacement parts or repairing FlashPak.

| Replacement Part                | Fits FlashPak Model Number |          |          |          |
|---------------------------------|----------------------------|----------|----------|----------|
|                                 | 9020                       | 9030     | 9040     | 9050     |
| Evolution Valve Conversion Kit  | 9090                       | 9090     | 9090     | 9090     |
| Evolution Valve Replacement Kit | 9090-90                    | 9090-90  | 9090-90  | 9090-90  |
| Valve Gasket                    | 9020-10                    | 9020-10  | 9020-10  | 9020-10  |
| Retaining Nuts (set of 8)       | 9090-92                    | 9090-92  | 9090-92  | 9090-92  |
| Carriage Bolts (set of 4)       | 9090-93                    | 9090-93  | 9090-93  | 9090-93  |
| Pressure Balance Element (PBE)  | 9020-11                    | 9020-11  | 9020-11  | 9020-11  |
| PBE Gasket                      | 9020-12                    | 9020-12  | 9020-12  | 9020-12  |
| Pressure Relief Button          | 9020-16M                   | 9020-16M | 9020-16M | 9020-16L |
| Base Only. No valve of PBE      | 9020-01                    | 9030-01  | 9040-01  | 9050-01  |
| Base w/PBE & Gasket No Valve    | 9020-13                    | 9030-13  | 9040-13  | 9050-13  |
| Lid Only w/Gasket. No Valve     | 9020-02                    | 9030-02  | 9040-02  | 9050-02  |
| Wire Basket                     | 9020-08                    | 9030-08  | 9040-08  | 9050-08  |
| Silicone Pin Mat for Basket     | 9020-04                    | 9030-04  | 9040-04  | 9050-04  |
| FlashPak Gloves (Pair)          | 9000-14                    | 9000-14  | 9000-14  | 9000-14  |
| FlashPak Delivery Cart          | 9111                       | 9111     | 9111     | 9111     |

## Symbol Definitions:

Symbol for Manufacturer; the name and address of the Manufacturer is next to this symbol

Caution—Consult accompanying documents

Lot number or batch number

Catalog number

Consult instructions for use

|                                                                                                                                                                                                                                                            |                                                                                                                                                                                                                                              |                                                                                                                                                                         |
|------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|-------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| <br><b>In USA Contact:</b><br>Symmetry Surgical<br>3034 Owen Drive<br>Antioch, TN 37013 USA<br>Phone: 800-251-3000<br>Fax: + 1-615-964-5566<br>www.symmetrysurgical.com | <br><b>In Europe contact:</b><br>Symmetry Surgical GmbH<br>Maybachstraße 10, 78532 Tuttlingen,<br>Germany<br>Phone: +49 7461 96490<br>Fax: +49 7461 77921 | <b>Notified Body</b><br>BSI Product Services<br>Maylands Avenue<br>Hemel Hempstead<br>Hertfordshire<br>HP2 4SQ UK<br>Tel: +44 (0)1442 278549<br>Fax: +44 (0)1442 278575 |
|------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|-------------------------------------------------------------------------------------------------------------------------------------------------------------------------|

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