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## INTRODUCTION

### Indications for Use

The ULTRA SYSTEM SP™ is a reusable sterilization container system designed specifically for use in STERRAD™ 50 and 100S Sterilization Systems. The container system is intended to allow sterilization of surgical instruments by health care providers and to maintain sterility of the enclosed devices for up to 180 days. Stainless steel baskets are provided to facilitate placement of the enclosed surgical instruments.

ULTRA SYSTEM SP™ containers are readily identified because of their silver colored finish. Standard ULTRA SYSTEM™ containers have a dark gray finish.

### Summary of Cautionary Statements

Please read the following cautionary statements prior to use.

- **CAUTION:** Sterility is the responsibility of the hospital. Please read and understand this manual before using the product.
- **CAUTION:** Use only ULTRA SYSTEM SP™ containers in STERRAD™ Sterilizers.
- **CAUTION:** Follow all instructions provided with your STERRAD™ Sterilizer.
- **CAUTION:** Complex instruments (e.g., air powered, endoscopes, and instruments with lumens and channels) should be prepared and sterilized according to the instrument manufacturer's instructions.
- **CAUTION:** Use only the stainless steel basket provided with the container. Do not use an aluminum basket.

- **CAUTION:** Only sterilize lumened instruments where the lumens meet the following criteria:
  - Must be stainless steel
  - Diameter is 3 mm or greater
  - Length is 400 mm or less
- **CAUTION:** Do not sterilize more than 8.35 pounds of instruments in containers used in the STRERRAD™ 50 or more than 16.70 pounds of instruments in containers used in the STERRAD™ 100S. A maximum of 5 lumens may be sterilized in each load.
- **CAUTION:** Clean and decontaminate the container system after each use.
- **CAUTION:** Inspect the container system and gasket carefully before and after each use. If the container or gasket is damaged, DO NOT use the container until it is repaired.
- **CAUTION:** If the filter is not visible at the inspection points prior to opening, the contents must be re-sterilized.
- **CAUTION:** If any problems are encountered at any time during the use of the container, the contents should not be used. The container should be inspected and repaired if necessary and the contents re-sterilized.
- **CAUTION:** Do not clean the container and basket with solutions that contain CAUSTIC SODA, TRISODIUM PHOSPHATE, or ACID SOLUTIONS stronger than TWO-TENTHS OF ONE PERCENT.
- **CAUTION:** The gasket must be replaced every 12 months. Read instructions carefully prior to replacement. Failure to fully seat the gasket may result in a compromised seal.
- **CAUTION:** The contents of the container will remain sterile up to 180 days after sterilization.

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## LOADING

**CAUTION:** Inspect the container system and gasket carefully before each use. If the container or gasket is damaged, DO NOT use the container until it is repaired.

**CAUTION:** Use only the stainless steel basket provided with the container. DO NOT use an aluminum basket.

1. Center a new filter over the raised vent area on the container top. Ensure that the filter overlaps the entire perimeter of the raised vent.


2. Hold the flexible retainer plate with both hands, one on the handle and the other on the side opposite the handle.

3. Guide the side of the plate opposite the handle snugly (hugging) against the raised vent area.

4. While pressing down firmly on the side opposite the handle, press down on the handle to seat the retainer plate.


5. Press down on the retainer plate with flat palms ensuring that all the locking tabs are fully engaged.

6. Repeat steps 1-5 to place a filter and filter retainer on the raised vent area on the container bottom.

7. Inspect the top and bottom to verify that the filters are visible at all inspection points.

8. Load the stainless steel basket with instruments according to hospital protocol. Make certain that instruments do not project above the top of the basket so that they would touch the container top and interfere with proper container closure. Place the filled stainless steel basket into the container bottom ensuring that the basket handles are down.

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

9. Place a process indicator inside the container according to hospital protocol.
10. Place the top onto the container bottom. Make sure the top is properly seated to ensure an effective seal of the gasket.
11. Close the latch on each end by placing the latch plates on top of the lock bar and pushing down on the latch release levers.

12. Insert a tamper-resistant locking key through all three openings in the latch at each end of the container. Insert from left to right with the process indicator dot facing toward you. Make sure each tamper-resistant locking key is fully inserted as shown and the arrow is completely visible beyond the second opening.

**CAUTION:** Complex instruments (e.g., air powered, endoscopes, and instruments with lumens and channels) should be prepared and sterilized according to the instrument manufacturer's instructions.

**CAUTION:** Only sterilize lumened instruments where the lumens meet the following criteria:

- Must be stainless steel
- Diameter is 3 mm or greater
- Length is 400 mm or less

**CAUTION:** Do not sterilize more than 8.35 pounds of instruments in containers used in the STERRAD™ 50 or more than 16.70 pounds of instruments in containers used in the STERRAD™ 100S. A maximum of 5 lumens may be sterilized in each load.

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## PROCESSING

1. Place the prepared containers into the STERRAD™ sterilizer according to instructions provided with the sterilizer.
2. Process the containers according to the sterilizer manufacturer's recommendations and/or hospital protocol.
3. After processing, inspect the dot on the tamper-resistant lock to ensure that sterilization conditions were met and the color changed from yellow to blue or red to blue.

**CAUTION:** Use only ULTRA SYSTEM SP™ containers in STERRAD™ sterilizers.

**CAUTION:** Follow all instructions provided with your STERRAD™ sterilizer.

**CAUTION:** The contents of the container will remain sterile up to 180 days after sterilization.


## OPENING

**CAUTION: If the filter is not visible at the inspection points prior to opening, the contents must be re-sterilized.**


1. Place the container on a surface that facilitates an aseptic opening.
2. Inspect the dot on each tamper-resistant lock to confirm that the color has changed from yellow or red to blue. Inspect each tamper-resistant lock to ensure the latch has not been opened.

3. Check the filter inspection points to ensure that the filters are properly in place.

Inspection Points


4. Unlock the top by simultaneously snapping both latch release levers open, breaking the tamper-resistant locking keys.

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
5. Remove the top by placing fingers under the lock bar and lifting up.

6. Remove the basket by grasping the handles and lifting straight up and out of the container.

7. Place the basket on a sterile surface.

8. After the surgical procedure, place instruments in the basket, place the basket in the container, and close the container.

**CAUTION: If any problems are encountered at any time during the inspection or use of the container, the contents should not be used. The container should be inspected and repaired if necessary and the contents should be re-sterilized.**

## CLEANING/DECONTAMINATION

**CAUTION: Clean and decontaminate the container system after each use.**

1. Unlock the top by quickly snapping both latch release levers open.

2. Remove the top by placing fingers under the lock bar and lifting up.

3. Remove instruments from the basket and process them according to hospital protocol.

4. Remove each retainer plate by pulling up on the retainer plate handle as shown.

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
5. Remove and discard filters.
  
6. Remove and discard any parts of the tamper-resistant locking keys remaining in the latches.
  
7. Inspect the container for damage such as:
  - ✓ broken latches
  - ✓ dented edge of top or bottom
  - ✓ damaged gasket
  - ✓ damaged filter retainer plate
  - ✓ any significant damage
  
8. Clean container tops, bottoms, filter retainer plates and baskets separately according to hospital protocol. Clean the gasket when cleaning the container top. It is not necessary to remove the gasket for cleaning. Use an instrument cleaning detergent suitable for use with anodized aluminum. A detergent with excellent rinsing and controlled sudsing is recommended. Please consult your cleaning/detergent supplier for more information.
  
9. If you use an automatic washer, ensure that container components are secured properly to avoid damage to the container during the cleaning cycle. Completely rinse the container after cleaning.

**CAUTION: Inspect the container system and gasket carefully after each use. If the container or gasket is damaged, DO NOT use the container until it is repaired.**

**CAUTION: Do not use solutions that contain the following: CAUSTIC SODA, TRISODIUM PHOSPHATE, or ACID SOLUTIONS stronger than TWO-TENTHS OF ONE PERCENT.**

## GASKET REPLACEMENT

**CAUTION:** The gasket must be replaced every 12 months.

1. Check the date code on the container top to determine replacement schedule.
2. Grasp the old gasket with a hemostat or similar instrument.

3. Pull up on the gasket and remove it from the groove in the container top.

4. Install the new gasket beginning at either end where the handles are located. Using your finger or thumb, press the new gasket into the groove of the container cover, stopping just before the corners.

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5. Press the gasket into the groove at the opposite end, stopping just before the corners.

6. Press the gasket into the grooves on both sides, stopping just before the corners.

7. Press the gasket into each corner, stretching as needed.

8. Using a smooth instrument, such as the ring handle of a hemostat, fully seat the gasket into the groove.

**CAUTION: Failure to fully seat the gasket may result in a compromised seal.**

**In USA Contact**

Symmetry Surgical  
3034 Owen Drive  
Antioch, TN 37013 USA  
Tel: 800-251-3000  
Fax: +1-615-964-5566  
www.symmetrysurgical.com

**In Europe Contact**

Symmetry Surgical GmbH  
Maybachstraße 10, 78532 Tuttlingen, Germany  
Phone: +49746196490  
Fax: +49746177921

**Notified Body**

BSI Product Services  
Maylands Avenue  
Hemel Hempstead  
Hertfordshire  
HP2 4SQ UK  
Phone: +44 (0)1442 278549  
Fax: +44 (0)1442 278575

The Ultra System is a trademark of Symmetry Surgical, Inc.  
11002-99 Rev. E

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