RHOTON® Titanium Micro Dissector Kit

Symmetry Surgical Inc.
3034 Owen Drive,
Antioch, TN 37013 USA
☎ 1-800-251-3000
Fax: 1-615-964-5566
wwwu u

Symmetry Surgical GmbH
Maybachstraße 10
78532 Tuttlingen, Germany
☎ +49 7461 96490
Fax: +49 7461 77921

LCN 800-1695-85/J
© 1999–2015 Symmetry Surgical
♻ Revised 06/1

ENGLISH

IMPORTANT INFORMATION

Please Read Before Use

RHOTON® Titanium Micro Dissector Kit

Description

The RHOTON® Micro Dissectors are made of titanium alloy and are used during aneurysm or tumor dissections, particularly acoustic neuromas.

The kit (catalog no. 80-1660) is provided nonsterile and must be sterilized before use.

Instruments (see Figure 1):

Round Dissector, small, catalog no. 80-1680
Round Dissector, medium, catalog no. 80-1681
Round Dissector, large, catalog no. 80-1682
General Purpose Elevator, small, catalog no. 80-1683
General Purpose Elevator, large, catalog no. 80-1684
Spatula Dissector, small, catalog no. 80-1685
Spatula Dissector, medium, catalog no. 80-1686
Spatula Dissector, large, catalog no. 80-1687
Semi-sharp Hook, 90°, catalog no. 80-1688
Blunt Hook, 90°, catalog no. 80-1689
Semi-sharp Hook, 45°, catalog no. 80-1690
Straight Point Needle, catalog no. 80-1691
Micro Curette, small, catalog no. 80-1692
Micro Curette, large, catalog no. 80-1693
Dissector Str. 3 mm, catalog no. 80-1679
Dissector 90° 3 mm, catalog no. 80-1696
Dissector 90° 5 mm, catalog no. 80-1697
Dissector 40° 4 mm, catalog no. 80-1698
Dissector 40° 8 mm, catalog no. 80-1699
Dissector case, catalog no. 80-1661

The above instruments are supplied with a convenient sterilizing-carrying case. Each instrument has an identification number etched on the base that corresponds with a numbered slot in the case.

Reprocessing Instructions

These reprocessing instructions apply to:

- Reusable surgical instruments supplied by Symmetry Surgical
- Instruments intended for reprocessing in a health care facility setting

These instructions **do not apply** to single-use devices.

These reprocessing instructions have been validated as being capable of preparing reusable Symmetry Surgical instruments for reuse. It is the responsibility of the reprocessor to ensure that the reprocessing is actually performed using appropriate equipment, materials, and personnel to achieve the desired result. This normally requires validation and routine monitoring of the process. Any deviation by the reprocessor from these instructions should be evaluated for effectiveness and potential adverse consequences.

Please telephone Symmetry Surgical Customer Service at 1 800 251 3000 or fax: 1 615 964 5566 if you need additional assistance.

Warnings for Reprocessing

- Follow the instructions and warnings issued by the suppliers of any cleaning and disinfection agents and equipment used.
- Do not exceed 140°C during reprocessing steps.
- Highly alkaline conditions can damage products with aluminum parts.
- Complex devices, such as those with tubes, hinges, retractable features, mated surfaces, and textured surface finishes, require special attention during cleaning. Manual pre-cleaning of such device features is required before automated cleaning processing.
- Avoid exposure to hypochlorite solutions, as these will promote corrosion.

Limitations on Reprocessing

- Repeated processing has minimal effects on instrument life and function.

- End of useful life is generally determined by wear or damage in surgical use.
- Carefully inspect instruments between uses to verify proper functioning.
- Send damaged instruments to a supplier of authorized repair or refurbishment services.

Decontamination Considerations – Creutzfeldt-Jakob Disease (CJD)

Under certain classifications of risk, the World Health Organization (WHO), or local regulatory authorities recommend special CJD inactivation processing procedures. Consult WHO and local regulations for further information.

Care at the Point of Use

- Clean instruments as soon as possible after use. If cleaning must be delayed, immerse instruments in a compatible detergent solution to prevent drying and encrustation of surgical soil.
- Avoid prolonged exposure to saline to minimize the chance of corrosion.
- Remove excessive soil with a disposable wipe.

Containment and Transportation

Reprocess instruments as soon as is reasonably possible after use.

Preparation for Cleaning

For instruments that require disassembly for cleaning, perform disassembly as instructed in the instructions for use provided with the instrument.

Manual Cleaning

A. ALL INSTRUMENTS

- Clean delicate instruments separately from other instruments.
- Disassemble instruments that are intended to come apart for cleaning.
- Prepare an enzymatic cleaning solution in accordance to the manufacturer's instructions.
- Soak soiled instruments for 5 minutes.
- Follow the additional steps in Cleaning sections *B.*, or *C.*, which follow. Select the category that most closely matches the type of instrument.
- Use a soft bristle brush to remove all traces of blood and debris; pay close attention to any hard to reach areas, textured surfaces, or crevices.
- Rinse the instrument thoroughly with warm tap water.
- Ultrasonically clean the instrument for 10 minutes in neutral pH detergent, prepared in accordance with the manufacturer's instructions.
- Rinse the instrument thoroughly with warm tap water.
- Dry the instrument immediately after final rinse.

B. Instruments with cannulations or lumens (i.e. tubes), or holes

- Follow the steps in the section *Manual Cleaning – A. ALL INSTRUMENTS.*
- When cleaning, use a tight-fitting, soft, non-metallic cleaning brush or pipe cleaner to scrub the cannula, lumen, or hole. Push in and out, using a twisting motion to remove debris. Use a syringe filled with enzymatic cleaning solution to flush hard to reach internal areas.
- When rinsing, pay particular attention to flush the cannulations, lumens, or holes with warm tap water.
- Dry internal areas with filtered compressed air.

C. Articulating Instruments (those with moveable parts)

- Follow the steps in the section *Manual Cleaning – A. ALL INSTRUMENTS.*
- When cleaning, fully immerse the instrument in the cleaning solution to avoid aerosol generation. Brush with a soft non-metallic bristle brush to remove all traces of blood and debris. Pay close attention to threads, crevices, seams, and any hard to reach areas. Actuate any moveable mechanisms, such as hinged joints, box locks, or spring-loaded features, to free trapped blood and debris. If the components of the instrument can be retracted, retract or open the part while cleaning the area. For instruments with flexible shafts, flex the instrument under the cleaning solution.
- When rinsing; pay particular attention to internal areas and moveable parts. Actuate moveable parts while rinsing. If the components of the instrument can be retracted, retract or open the part while rinsing the area. For instruments with flexible shafts, flex the instrument under the rinse solution.
- Dry internal areas with filtered compressed air.

Automated Cleaning

- For instrument types with complex design features, such as those described in *Manual Cleaning* sections B. and C. above, it is important to manually clean prior to automated processing to improve the removal of adherent soil. Brush the instrument, actuate mechanisms, agitate and/or irrigate under the surface of the cleaning solution to prevent the creation of aerosols.
- Clean, using the "INSTRUMENTS" cycle in a validated washer disinfectant and a pH neutral cleaning agent intended for use in automated cleaning. The cleaning cycle should incorporate enzyme pre-wash, wash, rinse, thermal rinse, and drying steps.
- Load instruments so that hinges are open and cannulations and holes can drain.
- Place heavier instruments on the bottom of containers. Do not place heavy instruments on top of delicate instruments.
- For instruments with concave surfaces, such as curettes, place instrument with the concave surface facing downward to facilitate draining.

Cleaning Inspection

- Inspect all instruments before sterilization or storage to ensure the complete removal of soil from surfaces, tubes and holes, moveable parts.
- If areas are difficult to inspect visually, check for blood by immersing or flushing the instrument in a 3% hydrogen peroxide solution. If bubbling is observed, blood is present. **Rinse instruments thoroughly after using hydrogen peroxide solution.**
- If soil is still present, reclean the instrument.

Disinfection

Symmetry Surgical instruments must be terminally sterilized prior to surgical use. See *Sterilization* instructions.

Maintenance

Between uses, lubricate moving parts with a water-soluble lubricant in accordance with the manufacturer's instructions. PRESERVE® Lubricant is recommended for instruments to be steam sterilized (REF 43-1033).

Inspection and Functional Testing

- Visually inspect the instrument and check for damage and wear.
- Cutting edges should be free of nicks and have a continuous edge.
- Jaws and teeth should align properly.
- Moveable parts should have smooth movement without excessive play.
- Locking mechanisms should fasten securely and close easily.
- Long, thin instruments should be free of bending and distortion.

Packaging

- If desired, use instrument trays to contain instruments that are provided in sets.
- Double wrap instruments in accordance with local procedures, using standard wrapping techniques such as those described in ANSI/AAMI ST46-1993.

Sterilization

- Use a validated, properly maintained and calibrated steam sterilizer.
- Effective steam sterilization can be achieved using the following cycles:

Cycle Type	Temperature	Exposure Time
Prevacuum	132–134°C	4 minutes
Prevacuum	134–137°C	3 minutes

Storage

Store sterile packaged instruments in a manner that provides protection from dust, moisture, insects, vermin, and extremes of temperature and humidity.

Service and Repair

For service and repair outside the United States, contact your local Symmetry Surgical representative.

Inside the United States, send instruments for service or repair to:

Symmetry Surgical
3034 Owen Drive,
Antioch, TN 37013 USA

Warranty

Symmetry Surgical warrants that this medical device is free from defects in both materials and workmanship. **Any other express or implied warranties, including warranties of merchantability or fitness, are hereby disclaimed. Suitability for use of this medical device for any particular surgical procedure should be determined by the user in conformance with the manufacturer's instructions for use. There are no warranties that extend beyond the description on the face hereof.**

PRESERVE® is a registered trademark of Symmetry Surgical
RHOTON® is a registered trademark of Al Rhoton, M.D.

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