# Bovie®

## Electrosurgical Electrodes; Sterile

REF

Dermal Tip Electrodes: A805, A807, A807DE  
H-Type Dermal Tip Electrodes: H10008, H10108

EN	Dermal Tip Electrodes; Sterile
DA	Dermal spids-elektroder; Sterile
NL	Dermale puntelektroden; Steriel
FI	Iholle käytettävät kärkielektrodit; Steriili
FR	Électrodes à pointe dermique ; Stériles
DE	Dermale Punktelektroden; Steril
IT	Elettrodi a punta dermica; Sterile
PR	Eléctrodos de Ponta Dérmica; Esterilizados
ES	Puntas dérmicas de electrodos Estériles
SV	Dermala spetsselektroder; Sterila
EL	Δερματικά ηλεκτρόδια άκρου, Στείρα
CN	皮肤端头电极；无菌
JA	経皮先端電極; 滅菌

R<sub>x</sub> ONLY


STERILE EO

STORAGE  
10 to 75 % non-condensing

OPERATE  
30 to 75 % non-condensing

TRANSPORT  
10 to 100 % including condensation

EN	Not made with natural rubber latex
DA	Ikke fremstillet med naturlig gummitatex
NL	Niet vervaardigd uit latex van natuurlijk rubber
FI	Valmistuksessa ei ole käytetty luonnonkumilateksia
FR	Exempt de latex de caoutchouc naturel
DE	Ohne Naturkautschuklatex hergestellt
IT	Non realizzato con lattice di gomma naturale
NO	Ikke laget med naturgummie latex
PT	Não fabricado com látex de borracha natural
ES	No está fabricado con látex de caucho natural
SV	Ej tillverkad med naturgummi
EL	Δεν είναι κατασκευασμένο από φυσικό καουτσούκ λάτεξ
CN	非天然橡膠乳膠製成
JA	天然ゴムラテックスで作られていません

## **EN** DERMAL TIP ELECTRODES; STERILE

### DESCRIPTION

Bovie Medical Corporation electrodes have a standard 2.36 mm ( $\frac{3}{32}$  in.) diameter stainless steel shaft and plastic insulator, and fit Bovie® and most other electrosurgical pencils.

### WARNINGS:

**Danger – Explosion Hazard:** Do not use electrosurgery in the presence of flammable anesthetics.

**Fire Hazard:** Both oxygen (O<sub>2</sub>) and Nitrous Oxide (N<sub>2</sub>O) support combustion. Avoid O<sub>2</sub> and N<sub>2</sub>O enriched atmospheres. Enriched atmospheres may result in fires and burns to patients or surgical personnel.

The metal shaft of the electrode must fit completely and securely into the pencil. If metal is visible, do not use the electrode. Exposed metal may result in electric shock or burn to the patient or surgical personnel.

Confirm proper electrosurgical settings prior to and during a procedure. Use the lowest power settings to achieve the desired effect. If increased power settings are requested, check the patient return electrode and all accessory connections before major power setting adjustments.

**Fire Hazard:** Always place the active Electrode in a clean, dry, insulated safety holster when not in use.

Electrosurgical accessories that are activated or hot from use can cause unintended burns to the patient or surgical personnel.

Electrosurgical accessories may cause fire or burn if placed close to or in contact with flammable materials such as gauze or surgical drapes.

**Fire/Explosion Hazard:** The following substances contribute to increased fire and explosion hazards in the operating room:

- Oxygen enriched environments
- Oxidizing agents such as nitrous oxide (N<sub>2</sub>O) atmospheres
  - \* Verify all anesthesia circuit connections are leak free before and during use of electrosurgery.
  - \* Verify endotracheal tubes are leak free and that the cuff seals properly to prevent oxygen leaks.
  - \* If an uncuffed tube is in use, pack the throat with wet sponges around the uncuffed tube.
  - \* If possible, stop supplemental oxygen at least one minute before and during use of electrosurgery.
  - \* Alcohol-based skin prepping agents and tinctures
  - \* Activate the electrosurgical unit only after vapors from skin prep solutions and tinctures have dissipated.
  - \* Naturally occurring flammable gases (such as methane) that may accumulate in body cavities.

**Fire Hazard:** The sparking and heating associated with electrosurgery can provide an ignition source.

**Observe fire precautions at all times:**

- When using electrosurgery in the presence of gases or flammable substances, prevent pooling of fluids and the accumulation of gases under surgical drapes. Tent drapes to allow vapors to mix with room air before using electrosurgery.
- Tissue buildup (eschar) on the tip of an active electrode poses a fire hazard, especially in oxygen enriched environments. With sufficient heating, eschar can become a glowing ember and pose a fire hazard both as an ignition source and as a fuel. Keep the electrode clean and free of all debris.
- Facial and other body hair is flammable. Water soluble surgical lubricating jelly may be used to cover hair close to the surgical site to decrease flammability.

---

This medical device cannot be effectively cleaned and/or sterilized by the user and therefore cannot be safely reused. It is intended for single use only. Any attempt by the user to clean and resterilize this device may result in bioincompatibility, infection, or other risks of device failure to the patient.

---

Conductive fluids (e.g., blood or saline) in direct contact with an active electrode or in close proximity to any active accessory may disperse electrical current and cause unintended burns to the patient. This can happen as a result of either direct coupling with the active electrode or capacitive coupling between the active electrode and the external surface of the electrode insulation. Therefore, to prevent unintended burns in the presence of conductive fluids:

- Always keep the external surface of the active electrode away from adjacent tissue while activating the electrosurgical generator.
- Clear conductive fluid from the electrode before activating the electrosurgical pencil.

**WARNING:** This product is sterilized utilizing Ethylene Oxide. These products can expose you to chemicals including Ethylene Oxide (EtO), which is known to the state of California to cause cancer and/or birth defects or other reproductive harm. For more information go to [www.P65Warnings.ca.gov](http://www.P65Warnings.ca.gov).

---

.....

**CAUTIONS:**

**R<sub>x</sub> ONLY**

Federal (USA) law restricts this device to sale by or on the order of a physician.

Use aseptic techniques to open packaging.

Before beginning the procedure, verify compatibility of all instruments and accessories. Generators and electrical accessories complying with IEC60601-1-1 and IEC60601-2-2 standards are deemed to be compatible.

Always refer to the instruction manuals of a high frequency generator.

During electrosurgery, follow general instructions for each surgical procedure.

.....

Device is rated for 4k Vpeak.

Always use the lowest power setting that achieves the desired surgical effect. Use the active electrode for the minimum time necessary in order to reduce the possibility of unintended burn injury.

Do not contact or arc the activated electrode to metal surfaces. Blade or needle breakage may result.

Confirm that the active electrode is properly connected to the active accessory, and that the accessory is firmly connected to the correct power output receptacle on the electrosurgical generator.

Before use, examine the electrosurgical unit and accessories for defects. Do not use cables or accessories with damaged (cracked, burned, or taped) insulation or connectors. Localized burns to the patient or physician may result.

Do not modify or add to the insulation of active electrodes.

Activate the electrosurgical unit only when you are ready to deliver electrosurgical current and the active tip is in view and near target tissue. Deactivate the electrosurgical unit before the tip leaves the surgical site.

The electrodes are intended for single use only. Discard after use in a puncture resistant container to prevent injury. These electrodes are not designed to withstand reesterilization. Do not reesterilize.

For monopolar electrosurgical procedures, ensure that the patient return electrode is appropriately selected, properly applied to the patient, and connected to the electrosurgical generator.

When applying the patient return electrode, ensure that the entire pad adheres to the patient. Do not reuse single use patient return electrodes. Do not use the return electrode if the packaging is damaged, or if the gel / adhesive material has dried.

If the patient has moved or is repositioned after the patient return electrode is applied, ensure that the return electrode firmly contacts the skin and that all connections are intact.

Avoid spilling any solutions on the return electrode. If this occurs, remove the return electrode, dry the skin, and apply a new return electrode. A moist environment may cause the return electrode to peel off.

**NOTICE:**

*Do not modify the electrode tip. Modifications to the tip may result in tip breakage or other damage.*

**IMPORTANT:**

*Clean the electrode often with moist gauze or other material appropriate for the electrode type.*

**Medical Purpose/ Indication**

Dermal Tip is used for coagulation of tissues

**Operating Principle**

Electrodes provide a point source transfer of RF energy to tissue

**Patient Population**

No restriction

**Site of Use and Site Conditions**

Site of use – No restriction

Site conditions – Aseptic

**Intended User Profile**

This device is intended to be used by trained physicians or individuals familiar with electrosurgery.

**Intended Conditions for Use**

Environment

- Doctor Office, Surgery Center, or Hospital
- Clear and unobstructed view of point of use

**Frequency of use**

Single patient multiple activation

**Disposal**

Biohazard – Dispose of electrode according to facility procedures

**Maximum Power****NOTICE:**

*Exceeding maximum power settings may result in patient injury or product damage.*

Bovie Medical Corporation recommends the following maximum power settings:

<b>Electrode</b>	<b>Mode</b>	<b>Power Setting</b>
Dermal Tip	Coag or Blend	55W

## Install the Electrode

1. Ensure the pencil is not connected to the generator, or the generator is Off or in the Standby mode, if available.
2. Grasp the insulating sleeve on the electrode. Remove tip protector. Insert the electrode into the pencil.
3. Ensure the electrode is fully inserted into the pencil. The shank and insulating sleeve should fit securely into the pencil. If the shank and/or insulating sleeve does not fit, or the insulation will not insert 3.2 mm ( $\frac{1}{8}$  in.), do not use this electrode/pencil combination.
4. Refer to the generator instruction manual for proper procedures for connecting the active accessory to the generator.

### STORAGE AND HANDLING:

Device must be stored in a clean, cool and dry place.

It must be handled with care to avoid damage to the packaging and its content during transportation and storage.

Instrument must not be in contact with acids or other corrosive liquids.

**Sterility is guaranteed unless the package is opened or damaged.**

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