


Our **RapidClean® Convertible Kerrisons** are designed for optimal performance with easy cleaning and reassembly.

- Convenient hinge-design provides access for optimal cleaning and assembly and reduces the possibility of misplaced parts during cleaning and sterilization
- Ceramic coating provides strong, hard surface for enhanced sharpness and durability with reduced glare
- One step thumb lever gives access to inner shaft for easy cleaning
- Available with and without a bone ejector



RapidClean®
Convertible Kerrisons



Visit symmetrysurgical.com or contact your local Sales Representative for more information.

CODE	DESCRIPTION
58-1180-BD	Convertible Kerrison Rongeur, 40° Up, 1mm, Thin Footplate, 7 Inches
58-2180-BD	Convertible Kerrison Rongeur, 40° Up, 2mm, Thin Footplate, 7 Inches
58-3180J-BD	Convertible Kerrison Rongeur, 40° Up, 3mm, Thin Footplate with Bone Ejector, 7 Inches
58-4180J-BD	Convertible Kerrison Rongeur, 40° Up, 4mm, Thin Footplate with Bone Ejector, 7 Inches
58-5180J-BD	Convertible Kerrison Rongeur, 40° Up, 5mm, Thin Footplate with Bone Ejector, 7 Inches
58-1200-BD	Convertible Kerrison Rongeur, 40° Up, 1mm, Thin Footplate, 8 Inches
58-2200-BD	Convertible Kerrison Rongeur, 40° Up, 2mm, Thin Footplate, 8 Inches
58-3200J-BD	Convertible Kerrison Rongeur, 40° Up, 3mm, Thin Footplate with Bone Ejector, 8 Inches
58-4200J-BD	Convertible Kerrison Rongeur, 40° Up, 4mm, Thin Footplate with Bone Ejector, 8 Inches
58-5200J-BD	Convertible Kerrison Rongeur, 40° Up, 5mm, Thin Footplate with Bone Ejector, 8 Inches
58-1230-BD	Convertible Kerrison Rongeur, 40° Up, 1mm, Thin Footplate, 9 Inches
58-2230-BD	Convertible Kerrison Rongeur, 40° Up, 2mm, Thin Footplate, 9 Inches
58-3230J-BD	Convertible Kerrison Rongeur, 40° Up, 3mm, Thin Footplate with Bone Ejector, 9 Inches
58-4230J-BD	Convertible Kerrison Rongeur, 40° Up, 4mm, Thin Footplate with Bone Ejector, 9 Inches
58-5230J-BD	Convertible Kerrison Rongeur, 40° Up, 5mm, Thin Footplate with Bone Ejector, 9 Inches

