



Symmetry Surgical's **RapidClean® Modular Carbon Fiber Handle (CFH) Laparoscopic** instrument line has more than 70 patterns of graspers, dissectors and scissors.

- All instruments have comfortable grip handles.
- All instruments include a monopolar electrocautery connection.
- Shafts are constructed of durable, high dielectric strength polymer.



RapidClean® Modular CFH Laparoscopic Instruments

RapidClean® Modular CFH Laparoscopic Instruments

- Can be purchased as complete instruments or individual components.
- May be disassembled for thorough cleaning.
- Are validated for reprocessing either assembled or disassembled.
- May reduce overall maintenance cost.

Inserts can be replaced if jaws become dull or damaged.

Flush port allows for cleaning without disassembly.

Non-indexing 360° rotation

Ergonomic handle design provides a comfortable grip.

Ratchet on demand included on all graspers.

Visit [symmetrysurgical.com](https://www.symmetrysurgical.com) or contact your local Sales Representative for more information on our entire line of RapidClean® Laparoscopic Instruments.

 Laparoscopic Instrument Sheath With Flush Port	Shaft Dimensions				Product Code
	Length		Diameter		
	33cm	45cm	5mm	10mm	
	x		x		96-0526
		x	x		96-0526-45
	x			x	96-0126
	x		x	96-0126-45	
	Non-Locking Handle				96-0002NL
	Locking Handle with Ratchet On Demand				96-0003R
	Traditional Ratchet Locking Handle				96-0004TR

	RapidClean® Modular CFH Laparoscopic Scissors	Shaft Dimensions				Product Code	
		Length		Diameter		Complete Instrument	Insert Only
	With 96-0002NL Non-Locking Handle	33cm	45cm	5mm	10mm		
	Metzenbaum Curved Scissors, Double Action	x		x		96-2003NL	96-2003Q
			x	x		96-2003NL-45	96-2003Q-45
	Metzenbaum Curved Serrated Scissors, Double Action	x		x		96-1145NL	96-1145Q
	Metzenbaum Curved Slimline Scissors, Double Action	x		x		96-2006NL	96-2006Q
			x	x		96-2006NL-45	96-2006Q-45
	Metzenbaum Straight Scissors, Double Action	x		x		96-2045NL	96-2045Q
	Hook Scissors, Single Action	x		x		96-2020NL	96-2020Q
			x	x		96-2020NL-45	96-2020Q-45
	Micro Scissors, Extremely Delicate Tips, Single Action	x		x		96-2044NL	96-2044Q
			x	x		96-2044NL-45	96-2044Q-45

	RapidClean® Modular CFH Laparoscopic Graspers	Shaft Dimensions				Product Code	
		Length		Diameter			
		With 96-0003R Locking Handle	33cm	45cm	5mm	10mm	Complete Instrument
	Wave Grasper With Atraumatic Serrations, Single Action	x		x		96-1051R	96-1051Q
			x	x		96-1051R-45	96-1051Q-45
	Wave Grasper With Atraumatic Serrations, Double Action	x		x		96-1053R	96-1053Q
			x	x		96-1053R-45	96-1053Q-45
	Fundus Grasper, Double Action	x		x		96-1050R	96-1050Q
			x	x		96-1050R-45	96-1050Q-45
	Babcock Grasper, Double Action	x		x		96-1089R	96-1089Q
			x	x		96-1089R-45	96-1089Q-45
	Babcock Grasper, Double Action	x			x	96-1079R	96-1079Q
			x		x	96-1079R-45	96-1079Q-45
	Allis Atraumatic Grasper, 3x4 Teeth, Double Action	x		x		96-1078R	96-1078Q
			x	x		96-1078R-45	96-1078Q-45
	Allis Grasper, Double Action	x		x		96-1087R	96-1087Q
	Allis Grasper, Double Row of Teeth, Single Action	x		x		96-1088R	96-1088Q
	Alligator Grasper, Double Action	x		x		96-1016R	96-1016Q
			x	x		96-1016R-45	96-1016Q-45
	Retraction Grasper, Single Action	x		x		96-1056R	96-1056Q
			x	x		96-1056R-45	96-1056Q-45

	RapidClean® Modular CFH Laparoscopic Graspers	Shaft Dimensions				Product Code	
		Length		Diameter			
		With 96-0003R Locking Handle	33cm	45cm	5mm	10mm	Complete Instrument
	Claw Grasper With 2x3 Teeth, Single Action	x		x		96-1005R	96-1005Q
		x			x	96-1015R	96-1015Q
			x		x	96-1015R-45	96-1015Q-45
	Snake Tooth (Cobra) Grasper, Double Action, 2x2 Teeth	x		x		96-1069R	96-1069Q
			x	x		96-1069R-45	96-1069Q-45
	Hunter Type Grasper, Double Action	x		x		96-1093R	96-1093Q
	DeBakey Straight Grasper, 40mm Jaw Length, Double Action	x		x		96-1092R	96-1092Q
			x	x		96-1092R-45	96-1092Q-45
	DeBakey Curved Grasper, Double Action	x		x		96-2070R	96-2070Q
	Bullet Nose Grasper, With Cup, Double Action	x		x		96-1033R	96-1033Q
			x	x		96-1033R-45	96-1033Q-45
	Tenaculum Forceps, 1x1 Sharp Teeth, Double Action	x		x		96-1008R	96-1008Q
			x	x		96-1008R-45	96-1008Q-45
		x			x	96-1009R	96-1009Q
			x		x	96-1009R-45	96-1009Q-45
	Strong Grasper, Double Action	x		x		96-1020R	96-1020Q
			x	x		96-1020R-45	96-1020Q-45
	Serrated Grasper, Rounded Jaws, Double Action	x		x		96-1024R	96-1024Q
			x	x		96-1024R-45	96-1024Q-45
	Cup Tooth Grasper, Double Action	x		x		96-1035R	96-1035Q
			x	x		96-1035R-45	96-1035Q-45

	RapidClean® Modular CFH Laparoscopic Graspers	Shaft Dimensions				Product Code	
		Length		Diameter			
		With 96-0003R Locking Handle	33cm	45cm	5mm	10mm	Complete Instrument
	Fenestrated Grasper, Slightly Flat Jaws, Double Action	x		x		96-1042R	96-1042Q
			x	x		96-1042R-45	96-1042Q-45
	Fenestrated Grasper, Micro Jaw, Double Action	x		x		96-1018R	96-1018Q
			x	x		96-1018R-45	96-1018Q-45
	Spoon Forceps, Rounded Jaw, Single Action	x			x	96-1072R	96-1072Q
	Biopsy Spoon, Single Action	x		x		96-1070R	96-1070Q
			x	x		96-1070R-45	96-1070Q-45
	Biopsy Spoon, With Spike, Double Action	x		x		96-1071R	96-1071Q
			x	x		96-1071R-45	96-1071Q-45

	RapidClean® Modular CFH Laparoscopic Dissectors	Shaft Dimensions				Product Code	
		Length		Diameter			
	With 96-0002NL Non-Locking Handle	33cm	45cm	5mm	10mm	Complete Instrument	Insert Only
	Micro Fenestrated Dissector, Double Action	x		x		96-1044NL	96-1044Q
	Bottle Nose (Dolphin) Dissector, Short Jaw, Double Action	x		x		96-1030NL	96-1030Q
			x	x		96-1030NL-45	96-1030Q-45
	Bullet Nose Dissector, With Cup, Double Action	x		x		96-1033NL	96-1033Q
			x	x		96-1033NL-45	96-1033Q-45
	Duckbill Forceps, Serrated, Double Action	x		x		96-1022NL	96-1022Q
			x	x		96-1022NL-45	96-1022Q-45
	Needle Nose Dissector, Double Action	x		x		96-1045NL	96-1045Q
	Maryland Dissector, Double Action	x		x		96-1067NL	96-1067Q
			x	x		96-1067NL-45	96-1067Q-45
	Maryland Dissector, Diamond Serrations, Double Action	x			x	96-1167NL	96-1167Q
	Mixer Dissector, 90° Angle, Double Action	x		x		96-1066NL	96-1066Q
			x	x		96-1066NL-45	96-1066Q-45
		x			x	96-1068NL	96-1068Q
			x		x	96-1068NL-45	96-1068Q-45
	Right Angle Dissector, Double Action	x		x		96-1063NL	96-1063Q

 Laparoscopic Instrument Sheath With Flush Port	Shaft Dimensions				Product Code
	Length		Diameter		
	33cm	45cm	5mm	10mm	
	x		x		96-0526
		x	x		96-0526-45
	x			x	96-0126
		x		x	96-0126-45

	Non-Locking Handle	96-0002NL
	Locking Handle with Ratchet On Demand	96-0003R
	Traditional Ratchet Locking Handle	96-0004TR

Visit symmetrysurgical.com or contact your local Sales Representative for more information on our entire line of RapidClean Laparoscopic Instruments.

