


# Symmetry Surgical<sup>®</sup> Secto<sup>®</sup> Dissectors & Tonsil Sponges

Symmetry Surgical's **Secto<sup>®</sup> Dissectors** are designed to provide a gentle, atraumatic way to separate tissue during surgical procedures. These sterile, cotton-constructed dissectors are x-ray detectable and available in four patterns: Cherry, Triangle, Peanut and Kittner.

- Available in a dispenser box including 96 packages of 5 each; 4 dispenser boxes in each case

Symmetry Surgical's **Secto<sup>®</sup> Tonsil Sponges** are sterile, cotton, non-linting, x-ray detectable gauzes with identification tape and string for use during tonsillectomies.

- Available 5 each per sterile pack, 25 packs per dispenser carton, 2 cartons per case, for a total of 250 per case.



Secto<sup>®</sup> Cherry  
Dissectors  
81-1002




Secto<sup>®</sup> Triangle  
Dissectors  
81-1006



Secto<sup>®</sup> Peanut  
Dissectors  
81-1000



Secto<sup>®</sup> Kittner  
Dissectors  
81-1001



Secto<sup>®</sup> Tonsil  
Sponges  
81-1012

Visit [symmetrysurgical.com](http://symmetrysurgical.com) or contact your local Sales Representative for more information.

### Secto® Cherry Dissectors

Dispenser box includes 96 packages of 5 each; 4 dispenser boxes to a case.

CODE	DESCRIPTION	DIAMETER	
81-1002	Sterile, X-Ray Detectable	3/8 in	10 mm

### Secto® Kittner Dissectors

Dispenser box includes 96 packages of 5 each; 4 dispenser boxes to a case.

CODE	DESCRIPTION	DIAMETER	
81-1001	Sterile, X-Ray Detectable	5/16 in	8 mm

### Secto® Peanut Dissectors

Dispenser box includes 96 packages of 5 each; 4 dispenser boxes to a case.

CODE	DESCRIPTION	DIAMETER	
81-1000	Sterile, X-Ray Detectable	3/8 in	10 mm

### Secto® Triangle Dissectors

Dispenser box includes 96 packages of 5 each; 4 dispenser boxes to a case.

CODE	DESCRIPTION	DIAMETER	
81-1006	Sterile, X-Ray Detectable	1/2 in	13 mm

### Secto® Tonsil Sponges

X-ray detectable, sterile, 5 per card/pouch, 25 pouches per dispenser, 2 dispenser boxes per case.

CODE	CORD LENGTH		DIAMETER	
81-1010	10in	254mm	5/8 in	16 mm
81-1011	10in	254mm	1 in	25 mm
81-1012	10in	254mm	1-1/4 in	32 mm

