



Symmetry Surgical can help you deliver exceptional care with easy-to-clean laparoscopic instruments. Our **RapidClean® Glide Lock instruments** offer an innovative two-piece, modular design allowing for easy, efficient cleaning and decontamination.

- Ergonomic **Glide Lock handle** locks on demand and remains out of the way when not needed
- High-volume flush port passes liquid through the shaft – aiding in the effective removal of debris before sterilization
- Devices validated to be terminally sterilized assembled or disassembled
- Available in standard and Opti-Length® Extended Length versions



**RapidClean® Glide Lock
Instruments**



For use with 99-21302 and 99-21302L



For use with 99-21323 and 99-21323L



For use with 99-21333 and 99-21333L



For use with 99-21345 and 99-21345L



For use with 99-21346 and 99-21346L



For use with 99-21364 and 99-21364L



For use with 99-21369 and 99-21369L

The RapidClean® Glide Lock Instruments below represent the complete instrument inclusive of the Glide Lock handle, insulated shaft and distal tip.

Visit symmetrysurgical.com/catalog for the full RapidClean® Glide Lock Portfolio.

RapidClean® Glide Lock Devices (Complete Device)

CODE	DESCRIPTION	Tip Size	Length
99-21302	5mm Traumatic Grasper, 2x3 Teeth, Insert, Handle & Lumen, Standard Length	22mm	37cm
99-21302L	5mm Traumatic Grasper, 2x3 Teeth, Insert, Handle & Lumen, Opti-Length®	22mm	49cm
99-21323	5mm Dorsey Fenestrated Grasper Insert, Handle & Lumen, Standard Length	30mm	38cm
99-21323L	5mm Dorsey Fenestrated Grasper Insert, Handle & Lumen, Opti-Length®	30mm	49cm
99-21333	5mm Atraumatic Grasper Insert, Handle & Lumen, Standard Length	22mm	37cm
99-21333L	5mm Atraumatic Grasper Insert, Handle & Lumen, Opti-Length®	22mm	49cm
99-21345	5mm Single Action Traumatic Grasper, 2x3 Teeth, Insert, Handle & Lumen, Standard Length	22mm	37cm
99-21345L	5mm Single Action Traumatic Grasper, 2x3 Teeth, Insert, Handle & Lumen, Opti-Length®	22mm	48cm
99-21346	5mm Atraumatic Grasper Single Action Insert, Handle & Lumen, Standard Length	22mm	37cm
99-21346L	5mm Atraumatic Grasper Single Action Insert, Handle & Lumen, Opti-Length®	22mm	48cm
99-21364	5mm Tenaculum Insert, Handle & Lumen, Standard Length	23mm	37cm
99-21364L	5mm Tenaculum Insert, Handle & Lumen, Opti-Length®	23mm	49cm
99-21369	5mm Fenestrated Bowel Grasper Insert, Handle & Lumen, Standard Length	20mm	37cm
99-21369L	5mm Fenestrated Bowel Grasper Insert, Handle & Lumen, Opti-Length®	20mm	48cm

Visit symmetrysurgical.com/rapidclean or contact your local Sales Representative for more information.

