




Scan the code above
to learn more about
our Energy portfolio.

Energy Portfolio



Energy



Electrosurgery Reference	ENG-2
Olsen® Single-Use Bayonet Bipolar Forceps	ENG-3
Olsen® Single-Use Bipolar Forceps	ENG-8
Olsen®/Ultra™ Reusable Bayonet Bipolar Forceps	ENG-12
Olsen®/Ultra™ Reusable Bipolar Forceps	ENG-22
Olsen® Single-Use Monopolar Forceps	ENG-26
Olsen® Reusable Bayonet Monopolar Forceps	ENG-27
Olsen® Reusable Monopolar Forceps	ENG-29
Olsen® Monopolar Electrodes	ENG-31
Monopolar Electrodes	ENG-34
Ultra™ Monopolar Laparoscopic Forceps	ENG-35

Olsen® Monopolar Blade Electrodes	ENG-50
Olsen® Electrode Extenders	ENG-51
Ultra™ Laparoscopic Electrodes	ENG-51
Olsen® Monopolar Loop Electrodes	ENG-52
Olsen® Monopolar Needle Electrodes	ENG-53
Olsen® Suction Coagulators	ENG-57
Olsen® Electrosurgical Cords	ENG-58
Olsen® Exposed Top Forceps	ENG-60
Olsen® Monopolar Adapters	ENG-61
Olsen® Monopolar Pencils	ENG-62
Olsen® Energy Accessories	ENG-64

Questions? Contact us at customerservice@symmetrysurgical.com or learn more at symmetrysurgical.com



The Surgical Instrument Desk Reference

How to Use Our Catalog

The Symmetry Surgical catalog is designed with you in mind, to be your trusted resource for surgical instrumentation. We took great care in creating an easy-to-use, robust catalog and we think you'll be pleased with the design as you become more familiar with our broad portfolio.

Layout

The catalog is organized by surgical specialty **and** product function.

And we know that it is difficult to pick up an instrument and know what specialty it's listed in. To make your job easier, we've sorted the instruments in our portfolio by surgical function, too, including Energy, Laparoscopic, Retractors and Sterilization Solutions.

U.S. Customer Service

Your Symmetry Surgical Customer Service team wants to make ordering surgical instruments the easiest part of your day. We accept orders via phone, fax, email, our website, Global Healthcare Exchange (GHX) and through the mail.

PHONE

800.251.3000

(Available Monday through Friday, 8:30 AM–6:30 PM EST)

FAX

800.342.3272

EMAIL

customerservice@symmetrysurgical.com

ONLINE

symmetrysurgical.com

MAILING ADDRESS

Symmetry Surgical
Attn: Customer Service
3034 Owen Drive
Antioch, TN 37013

The Surgical Instrument Desk Reference

International Customer Service

PHONE

Symmetry Surgical Switzerland GmbH
+41 (0) 52 557 59 00
(Available Monday through Friday, 7:00-16:00 GMT)

Symmetry Surgical, Inc., USA
615.964.5672
(Available Monday through Friday, 8:30 AM-6:30 PM EST)

FAX

Symmetry Surgical Switzerland GmbH
+41 (0) 52 557 59 99

EMAIL

customerservice@symmetrysurgical.com
internationalorders@symmetrysurgical.com

ONLINE

symmetrysurgical.com

MAILING ADDRESS

Symmetry Surgical Switzerland GmbH
Bühlstrasse 27
CH-8200 Schaffhausen
Switzerland

Symmetry Surgical, Inc.
3034 Owen Drive
Antioch, TN 37013

INTERNATIONAL PAYMENT TERMS/RETURN GOODS POLICY/INSTRUMENT REPAIR & RESTORE PROGRAM/LITERATURE REQUESTS

Please contact your local distributor. Locate a local distributor contact, International Customer Service at +41 (0) 52 557 59 00 or 615.964.5672 (US).

Cardiovascular/Thoracic

Energy

Ear, Nose & Throat

General

Laparoscopic

Neurology

Obesity

OB/GYN

Ophthalmology

Orthopedics & Arthroscopy

Pediatrics

Rectal

Retractors


Spine

Urology

Sterilization Solutions

Alphabetical & Numeric Indexes

Electrosurgery Reference






Abbreviation Key

“L” Total Length


“WL” Working Length

Forceps Measuring Key

To determine the correct length of desired forceps, measure per corresponding key as illustrated. Use the rulers to the left to help determine the forcep length.



Forceps Tip Sizes



Single-use products are sold sterile.
Kits Include Forceps, 12 ft (3.6 m) Cord and Surgi-Wipe.

Olsen® Single-Use Bayonet Bipolar Forceps



20-1320

Single-Use Bayonet Bipolar Forceps

Smooth Tips, 10 per Box

Length (L)	6 ½ in 165 mm	6 ½ in 165 mm	6 ½ in 165 mm	6 ½ in 165 mm
Working (WL)	3 ¾ in 95 mm	3 ¾ in 95 mm	3 ¾ in 95 mm	3 ¾ in 95 mm
TIP SIZE	FORCEPS, NON-INSULATED	FORCEPS, INSULATED	KIT, NON-INSULATED	KIT, INSULATED
0.5 mm	20-1315	20-1315I	20-1315K	20-1315KI
0.75 mm	20-1320-OL	20-1320I	20-1320K	20-1320KI
1.0 mm	20-1325	20-1325I	20-1325K	20-1325KI
1.5 mm	20-1330	20-1330I	20-1330K	20-1330KI



20-1370

Single-Use Bayonet Bipolar Forceps

Smooth Tips, 10 per Box

Length (L)	7 ¼ in 184 mm	7 ¼ in 184 mm	7 ¼ in 184 mm	7 ¼ in 184 mm
Working (WL)	2 ¾ in 70 mm	2 ¾ in 70 mm	2 ¾ in 70 mm	2 ¾ in 70 mm
TIP SIZE	FORCEPS, NON-INSULATED	FORCEPS, INSULATED	KIT, NON-INSULATED	KIT, INSULATED
0.5 mm	20-1355	20-1355I	20-1355K	20-1355KI
0.75 mm	20-1360	20-1360I	20-1360K	20-1360KI
1.0 mm	20-1340	20-1340I	20-1340K	20-1340KI
1.5 mm	20-1370	20-1370I	20-1370K	20-1370KI
2.0 mm	20-1375	20-1375I	20-1375K	20-1375KI

Single-use products are sold sterile.
Kits Include Forceps, 12 ft (3.6 m) Cord and Surgi-Wipe.

Questions? Contact us at customerservice@symmetrysurgical.com
or learn more at symmetrysurgical.com

Olsen® Single-Use Bayonet Bipolar Forceps



Single-Use Bayonet Bipolar Forceps

Smooth Tips, 10 per Box

Length (L)	7 ¾ in 197 mm	7 ¾ in 197 mm	7 ¾ in 197 mm	7 ¾ in 197 mm
Working (WL)	3 in 76 mm	3 in 76 mm	3 in 76 mm	3 in 76 mm

TIP SIZE	FORCEPS, NON-INSULATED	FORCEPS, INSULATED	KIT, NON-INSULATED	KIT, INSULATED
0.5 mm	20-1386	20-1386I	20-1386K	20-1386KI
0.75 mm	20-1384	20-1384I	20-1384K	20-1384KI
1.0 mm	20-1382	20-1382I	20-1382K	20-1382KI
1.5 mm	20-1380	20-1380I	20-1380K	20-1380KI




Single-Use Bayonet Bipolar Forceps

Smooth Tips, 10 per Box

Length (L)	8 ¼ in 210 mm	8 ¼ in 210 mm	8 ¼ in 210 mm	8 ¼ in 210 mm
Working (WL)	4 ⅞ in 105 mm	4 ⅞ in 105 mm	4 ⅞ in 105 mm	4 ⅞ in 105 mm

TIP SIZE	FORCEPS, NON-INSULATED	FORCEPS, INSULATED	KIT, NON-INSULATED	KIT, INSULATED
0.5 mm	20-1420	20-1420I	20-1420K	20-1420KI
0.75 mm	20-1400	20-1400I	20-1400K	20-1400KI
1.0 mm	20-1440	20-1440I	20-1440K	20-1440KI
1.5 mm	20-1410	20-1410I	20-1410K	20-1410KI
2.0 mm	20-1415	20-1415I	20-1415K	20-1415KI



Single-Use Bayonet Bipolar Forceps, Angled Up

Smooth Tips, 10 per Box

Length (L)	8 ¼ in 210 mm	8 ¼ in 210 mm	8 ¼ in 210 mm	8 ¼ in 210 mm
Working (WL)	4 ⅞ in 105 mm	4 ⅞ in 105 mm	4 ⅞ in 105 mm	4 ⅞ in 105 mm

TIP SIZE	FORCEPS, NON-INSULATED	FORCEPS, INSULATED	KIT, NON-INSULATED	KIT, INSULATED
1.0 mm	20-1445	20-1445I	20-1445K	20-1445KI

Single-use products are sold sterile.
Kits Include Forceps, 12 ft (3.6 m) Cord and Surgi-Wipe.

Olsen® Single-Use Bayonet Bipolar Forceps



Single-Use Bayonet Bipolar Forceps

Smooth Tips, 10 per Box

Length (L)	10 ½ in 267 mm	10 ½ in 267 mm	10 ½ in 267 mm	10 ½ in 267 mm
Working (WL)	6 in 152 mm	6 in 152 mm	6 in 152 mm	6 in 152 mm

TIP SIZE	FORCEPS, NON-INSULATED	FORCEPS, INSULATED	KIT, NON-INSULATED	KIT, INSULATED
0.5 mm	20-1484	20-1484I	20-1484K	20-1484KI
0.75 mm	20-1486	20-1486I	20-1486K	20-1486KI
1.0 mm	20-1480	20-1480I	20-1480K	20-1480KI
1.5 mm	20-1481	20-1481I	20-1481K	20-1481KI



Single-Use Delashaw* Bayonet Bipolar Forceps

Smooth Tips, 10 per Box

CODE	DESCRIPTION	TIP SIZE	LENGTH
20-1350	Forceps, Non-Insulated	Micro	7 ¼ in 184 mm
20-1350I	Forceps, Insulated	Micro	7 ¼ in 184 mm
20-1350K	Kit, Non-Insulated	Micro	7 ¼ in 184 mm
20-1350KI	Kit, Insulated	Micro	7 ¼ in 184 mm

*Developed by Johnny Delashaw, Jr., M.D.

Single-Use Bayonet Bipolar Forceps, Acute Angle

Smooth Tips, 10 per Box

Length (L)	6 ¼ in 159 mm	6 ¼ in 159 mm	6 ¼ in 159 mm	6 ¼ in 159 mm
Working (WL)	3 ½ in 89 mm	3 ½ in 89 mm	3 ½ in 89 mm	3 ½ in 89 mm

TIP SIZE	FORCEPS, NON-INSULATED	FORCEPS, INSULATED	KIT, NON-INSULATED	KIT, INSULATED
0.75 mm	20-1500	20-1500I	20-1500K	20-1500KI
1.0 mm	20-1510	20-1510I	20-1510K	20-1510KI
1.5 mm	20-1520	20-1520I	20-1520K	20-1520KI

Single-use products are sold sterile.

Kits Include Forceps, 12 ft (3.6 m) Cord and Surgi-Wipe.

Questions? Contact us at customerservice@symmetrysurgical.com
or learn more at symmetrysurgical.com

Olsen® Single-Use Bayonet Bipolar Forceps

Single-Use Bayonet Bipolar Forceps, Acute Angle

Smooth Tips, 10 per Box

Length (L)	7½ in 191 mm	7½ in 191 mm	7½ in 191 mm	7½ in 191 mm
Working (WL)	4⅛ in 105 mm	4⅛ in 105 mm	4⅛ in 105 mm	4⅛ in 105 mm

TIP SIZE	FORCEPS, NON-INSULATED	FORCEPS, INSULATED	KIT, NON-INSULATED	KIT, INSULATED
0.75 mm	20-1530	20-1530I	20-1530K	20-1530KI
1.0 mm	20-1540	20-1540I	20-1540K	20-1540KI
1.5 mm	20-1550	20-1550I	20-1550K	20-1550KI

Length (L)	8⅜ in 213 mm	8⅜ in 213 mm	8⅜ in 213 mm	8⅜ in 213 mm
Working (WL)	5⅞ in 130 mm	5⅞ in 130 mm	5⅞ in 130 mm	5⅞ in 130 mm

TIP SIZE	FORCEPS, NON-INSULATED	FORCEPS, INSULATED	KIT, NON-INSULATED	KIT, INSULATED
0.75 mm	20-1560	20-1560I	20-1560K	20-1560KI
1.0 mm	20-1570	20-1570I	20-1570K	20-1570KI
1.5 mm	20-1580	20-1580I	20-1580K	20-1580KI



Single-Use Bayonet Bipolar Forceps, Angled Up

Smooth Tips, 10 per Box

Length (L)	10½ in 267 mm	10½ in 267 mm	10½ in 267 mm	10½ in 267 mm
Working (WL)	6 in 152 mm	6 in 152 mm	6 in 152 mm	6 in 152 mm

TIP SIZE	FORCEPS, NON-INSULATED	FORCEPS, INSULATED	KIT, NON-INSULATED	KIT, INSULATED
1.0 mm		20-1682I		20-1682KI



Single-Use Bayonet Bipolar Forceps, with Midas Touch® Non-Stick Tips

Smooth Tips, 5 per Box

Length (L)	6½ in 165 mm	7½ in 191 mm	8¼ in 210 mm	10½ in 267 mm
Working (WL)	3¾ in 95 mm	3 in 76 mm	4⅛ in 105 mm	6 in 152 mm

TIP SIZE	INSULATED	INSULATED	INSULATED	INSULATED
0.5 mm	20-1315IMT	20-1386IMT	20-1420IMT	
0.75 mm	20-1320IMT	20-1384IMT	20-1400IMT	
1.0 mm	20-1325IMT	20-1382IMT	20-1440IMT	20-1480IMT
1.5 mm	20-1330IMT	20-1380IMT	20-1410IMT	
2.0 mm			20-1415IMT	


Single-use products are sold sterile.
Kits Include Forceps, 12 ft (3.6 m) Cord and Surgi-Wipe.

Olsen® Single-Use Bayonet Bipolar Forceps

Single-Use Bayonet Tonsil Bipolar Forceps

Smooth Tips, Angled Down, Insulated, 10 per Box

Length (L)	7 ½ in	191 mm	7 ½ in	191 mm
Working (WL)	3 in	76 mm	3 in	76 mm
TIP SIZE				
FORCEPS, INSULATED		KIT, INSULATED		
1.5 mm	20-1381I	20-1381KI		



Single-Use Irrigating Bayonet Bipolar Forceps

Smooth Tips, 5 per Box

Length (L)	6 ¼ in 159 mm	7 ¼ in 184 mm	8 ¼ in 210 mm	10 ½ in 267 mm
Working (WL)	2 ¾ in 70 mm	2 ¾ in 70 mm	4 ¼ in 108 mm	6 ¼ in 159 mm
TIP SIZE				
INSULATED		INSULATED	INSULATED	INSULATED
0.5 mm	20-0311	20-0361	20-0411	
1.0 mm	20-0321	20-0371	20-0421	20-0481
1.5 mm	20-0331	20-0381	20-0431	20-0491
Micro		20-0351		



Single-Use Irrigating Bayonet Bipolar Forceps, Single-Use Kits for Rotary Pumps

Smooth Tips, 5 per Box, with 12 ft Integrated Tubing/Cord Set

Length (L)	6 ¼ in 159 mm	7 ¼ in 184 mm	8 ¼ in 210 mm	10 ½ in 267 mm
Working (WL)	2 ¾ in 70 mm	2 ¾ in 70 mm	4 ¼ in 108 mm	6 ¼ in 159 mm
TIP SIZE				
INSULATED		INSULATED	INSULATED	INSULATED
0.5 mm	20-0311R	20-0361R	20-0411R	
1.0 mm	20-0321R	20-0371R	20-0421R	20-0481R
1.5 mm	20-0331R	20-0381R	20-0431R	20-0491R

Single-use products are sold sterile.

Kits Include Forceps, 12 ft (3.6 m) Cord and Surgi-Wipe.

Questions? Contact us at customerservice@symmetrysurgical.com or learn more at symmetrysurgical.com

Olsen® Single-Use Bipolar Forceps

Single-Use Adson Bipolar Forceps, Serrated Tip

10 per Box

CODE	DESCRIPTION	TIP SIZE	LENGTH
20-1080	Forceps, Non-Insulated	1.0 mm	5 in 127 mm
20-1080I	Forceps, Insulated	1.0 mm	5 in 127 mm
20-1080K	Kit, Non-Insulated	1.0 mm	5 in 127 mm
20-1080KI	Kit, Insulated	1.0 mm	5 in 127 mm



Single-Use Adson Bipolar Forceps, Smooth Tip

10 per Box

CODE	DESCRIPTION	TIP SIZE	LENGTH
20-1082	Forceps, Non-Insulated	1.0 mm	5 in 127 mm
20-1082I	Forceps, Insulated	1.0 mm	5 in 127 mm
20-1082K	Kit, Non-Insulated	1.0 mm	5 in 127 mm
20-1082KI	Kit, Insulated	1.0 mm	5 in 127 mm



Single-Use Adson Bipolar Forceps, with Midas Touch® Non-Stick Tips

5 per Box

CODE	DESCRIPTION	TIP SIZE	LENGTH
20-1082IMT	Forceps, Non-Stick	1.0 mm	5 in 127 mm



Single-Use Coaptation Bipolar Forceps

Smooth Tips, 10 per Box

CODE	DESCRIPTION	TIP SIZE	LENGTH
20-1010-OL	Forceps, Non-Insulated	Micro	4 ½ in 114 mm
20-1010I	Forceps, Insulated	Micro	4 ½ in 114 mm
20-1010K	Kit, Non-Insulated	Micro	4 ½ in 114 mm
20-1010KI	Kit, Insulated	Micro	4 ½ in 114 mm



Single-use products are sold sterile.
Kits Include Forceps, 12 ft (3.6 m) Cord and Surgi-Wipe.

Olsen® Single-Use Bipolar Forceps

Single-Use Cushing Bipolar Forceps

Smooth Tips, 10 per Box

CODE	DESCRIPTION	TIP SIZE	LENGTH
20-1160	Forceps, Non-Insulated	1.0 mm	7 in 178 mm
20-1160I	Forceps, Insulated	1.0 mm	7 in 178 mm
20-1160K	Kit, Non-Insulated	1.0 mm	7 in 178 mm
20-1160KI	Kit, Insulated	1.0 mm	7 in 178 mm
20-1162	Forceps, Non-Insulated	1.5 mm	7 in 178 mm
20-1162I	Forceps, Insulated	1.5 mm	7 in 178 mm
20-1162K	Kit, Non-Insulated	1.5 mm	7 in 178 mm
20-1162KI	Kit, Insulated	1.5 mm	7 in 178 mm



Single-Use Gerald Bipolar Forceps

Smooth Tips, 10 per Box

CODE	DESCRIPTION	TIP SIZE	LENGTH
20-1180	Forceps, Non-Insulated	1.0 mm	7 ¼ in 184 mm
20-1180I	Forceps, Insulated	1.0 mm	7 ¼ in 184 mm
20-1180K	Kit, Non-Insulated	1.0 mm	7 ¼ in 184 mm
20-1180KI	Kit, Insulated	1.0 mm	7 ¼ in 184 mm



Single-Use Iris Bipolar Forceps, Straight

Serrated Tips, 10 per Box

CODE	DESCRIPTION	TIP SIZE	LENGTH
20-1030	Forceps, Non-Insulated	1.0 mm	4 ½ in 114 mm
20-1030I	Forceps, Insulated	1.0 mm	4 ½ in 114 mm
20-1030K	Kit, Non-Insulated	1.0 mm	4 ½ in 114 mm
20-1030KI	Kit, Insulated	1.0 mm	4 ½ in 114 mm



Single-Use Iris Bipolar Forceps, Angled

Serrated Tips, 10 per Box

CODE	DESCRIPTION	TIP SIZE	LENGTH
20-1040	Forceps, Non-Insulated	1.0 mm	4 ½ in 114 mm
20-1040I	Forceps, Insulated	1.0 mm	4 ½ in 114 mm
20-1040K	Kit, Non-Insulated	1.0 mm	4 ½ in 114 mm
20-1040KI	Kit, Insulated	1.0 mm	4 ½ in 114 mm



Single-use products are sold sterile.

Kits Include Forceps, 12 ft (3.6 m) Cord and Surgi-Wipe.

Questions? Contact us at customerservice@symmetrysurgical.com or learn more at symmetrysurgical.com

Olsen® Single-Use Bipolar Forceps

Single-Use Jewelers Bipolar Forceps, Straight

Smooth Tips, 10 per Box

CODE	DESCRIPTION	TIP SIZE	LENGTH
20-1060	Forceps, Non-Insulated	Micro	4 ¾ in 121 mm
20-1060I	Forceps, Insulated	Micro	4 ¾ in 121 mm
20-1060K	Kit, Non-Insulated	Micro	4 ¾ in 121 mm
20-1060KI	Kit, Insulated	Micro	4 ¾ in 121 mm



Single-Use Jewelers Bipolar Forceps, Curved

Smooth Tips, 10 per Box

CODE	DESCRIPTION	TIP SIZE	LENGTH
20-1070	Forceps, Non-Insulated	Micro	4 ¾ in 121 mm
20-1070I	Forceps, Insulated	Micro	4 ¾ in 121 mm
20-1070K	Kit, Non-Insulated	Micro	4 ¾ in 121 mm
20-1070KI	Kit, Insulated	Micro	4 ¾ in 121 mm



Single-Use Jewelers Bipolar Forceps, Straight, Ridged Handle

Smooth Tips, 10 per Box

CODE	DESCRIPTION	TIP SIZE	LENGTH
20-1203	Forceps, Non-Insulated	Micro	3 ⅞ in 98 mm
20-1203I	Forceps, Insulated	Micro	3 ⅞ in 98 mm
20-1203K	Kit, Non-Insulated	Micro	3 ⅞ in 98 mm
20-1203KI	Kit, Insulated	Micro	3 ⅞ in 98 mm



Single-Use Jewelers Bipolar Forceps, with Midas Touch® Non-Stick Tips

Insulated, Smooth Tips, 5 per Box

CODE	DESCRIPTION	TIP SIZE	LENGTH
20-1060IMT	Forceps, Non-Stick	Micro	4 ¾ in 121 mm



Single-use products are sold sterile.
Kits Include Forceps, 12 ft (3.6 m) Cord and Surgi-Wipe.

Olsen® Single-Use Bipolar Forceps

Single-Use McPherson Bipolar Forceps, Angled

Smooth Tips, 10 per Box


CODE	DESCRIPTION	TIP SIZE	LENGTH
20-1200	Forceps, Non-Insulated	0.5 mm	3 ⁷ / ₈ in 98 mm
20-1200I	Forceps, Insulated	0.5 mm	3 ⁷ / ₈ in 98 mm
20-1200K	Kit, Non-Insulated	0.5 mm	3 ⁷ / ₈ in 98 mm
20-1200KI	Kit, Insulated	0.5 mm	3 ⁷ / ₈ in 98 mm



Single-Use McPherson Bipolar Forceps, Straight

Smooth Tips, 10 per Box

CODE	DESCRIPTION	TIP SIZE	LENGTH
20-1201	Forceps, Non-Insulated	0.5 mm	3 ⁷ / ₈ in 98 mm
20-1201I	Forceps, Insulated	0.5 mm	3 ⁷ / ₈ in 98 mm
20-1201K	Kit, Non-Insulated	0.5 mm	3 ⁷ / ₈ in 98 mm
20-1201KI	Kit, Insulated	0.5 mm	3 ⁷ / ₈ in 98 mm
20-1020	Forceps, Non-Insulated	0.5 mm	4 ¹ / ₄ in 108 mm
20-1020I	Forceps, Insulated	0.5 mm	4 ¹ / ₄ in 108 mm
20-1020K	Kit, Non-Insulated	0.5 mm	4 ¹ / ₄ in 108 mm
20-1020KI	Kit, Insulated	0.5 mm	4 ¹ / ₄ in 108 mm



Single-Use Semkin Bipolar Forceps

Smooth Tips, 10 per Box

CODE	DESCRIPTION	TIP SIZE	LENGTH
20-1190	Forceps, Non-Insulated	1.0 mm	5 ¹ / ₂ in 140 mm
20-1190I	Forceps, Insulated	1.0 mm	5 ¹ / ₂ in 140 mm
20-1190K	Kit, Non-Insulated	1.0 mm	5 ¹ / ₂ in 140 mm
20-1190KI	Kit, Insulated	1.0 mm	5 ¹ / ₂ in 140 mm
20-1090	Forceps, Non-Insulated	0.75 mm	5 ³ / ₄ in 146 mm
20-1090I	Forceps, Insulated	0.75 mm	5 ³ / ₄ in 146 mm
20-1090K	Kit, Non-Insulated	0.75 mm	5 ³ / ₄ in 146 mm
20-1090KI	Kit, Insulated	0.75 mm	5 ³ / ₄ in 146 mm



Single-Use Semkin Bipolar Forceps, with Midas Touch® Non-Stick Tips

Insulated, Smooth Tips, 5 per Box

CODE	DESCRIPTION	TIP SIZE	LENGTH
20-1090IMT	Forceps, Non-Stick	0.75 mm	5 ³ / ₄ in 146 mm



Single-use products are sold sterile.
Kits Include Forceps, 12 ft (3.6 m) Cord and Surgi-Wipe.

Questions? Contact us at customerservice@symmetrysurgical.com
or learn more at symmetrysurgical.com

Olsen®/Ultra™ Reusable Bayonet Bipolar Forceps



Bayonet Bipolar Forceps

Reusable, Smooth Tips, 1 Each

Length (L)	6½ in 165 mm	7¼ in 184 mm	8¼ in 210 mm	10½ in 267 mm
Working (WL)	2¾ in 70 mm	2¾ in 70 mm	4¼ in 108 mm	6¼ in 159 mm

TIP SIZE	NON-INSULATED	NON-INSULATED	NON-INSULATED	NON-INSULATED
0.5 mm	10-1311	10-1360	10-1421	10-1483
0.75 mm	10-1321	10-1361	10-1401	10-1485
1.0 mm	10-1322	10-1362	10-1431	10-1481
1.5 mm	10-1331	10-1371	10-1411	10-1484
2.0 mm		10-1372	10-1412	10-1487

Length (L)	6¼ in 159 mm	7⅞ in 181 mm	8¼ in 210 mm	10½ in 267 mm
Working (WL)	2¾ in 70 mm	2¾ in 70 mm	4¼ in 108 mm	6¼ in 159 mm

TIP SIZE	INSULATED	INSULATED	INSULATED	INSULATED
0.5 mm	10-1311I	10-1360I	10-1421I	10-1483I
0.75 mm	10-1321I	10-1361I	10-1401I	10-1485I
1.0 mm	10-1322I	10-1362I	10-1431I	10-1481I
1.5 mm	10-1331I	10-1371I	10-1411I	10-1484I
2.0 mm		10-1372I	10-1412I	10-1487I



Bayonet Bipolar Forceps

Reusable, Smooth Tips, 1 Each

Length (L)	6¾ in 171 mm	7¾ in 197 mm	7¾ in 197 mm	8⅝ in 219 mm	12 in 305 mm
Working (WL)	2¾ in 70 mm	3 in 76 mm	3⅝ in 92 mm	4½ in 114 mm	7⅞ in 200 mm

TIP SIZE	INSULATED	INSULATED	INSULATED	INSULATED	INSULATED
Micro	80205	80304			
0.5 mm	80204	80303	80434	80574	
0.75 mm			80433	80573	
1.0 mm	80202	80302			80582
1.2 mm			80432	80572	
1.5 mm	80201	80301	80431	80571	80581
2.0 mm	80200	80300	80430	80570	80580

Single-use products are sold sterile.
Kits Include Forceps, 12 ft (3.6 m) Cord and Surgi-Wipe.

Olsen®/Ultra™ Reusable Bayonet Bipolar Forceps



Bayonet Bipolar Forceps, Acute Angle

Reusable, Smooth Tips, 1 Each

Length (L)	6 1/8 in 156 mm	7 1/2 in 191 mm	8 3/8 in 213 mm	6 1/8 in 156 mm	7 1/2 in 191 mm	8 3/8 in 213 mm
Working (WL)	2 5/8 in 67 mm	2 3/4 in 70 mm	8 3/8 in 213 mm	2 5/8 in 67 mm	2 3/4 in 70 mm	8 3/8 in 213 mm

TIP SIZE	NON-INSULATED	NON-INSULATED	NON-INSULATED	INSULATED	INSULATED	INSULATED
1.0 mm	10-1501	10-1531	10-1561	10-1501I	10-1531I	10-1561I
1.5 mm	10-1511	10-1541	10-1571	10-1511I	10-1541I	10-1571I
2.0 mm	10-1521	10-1551	10-1581	10-1521I	10-1551I	10-1581I



Bayonet Bipolar Forceps, Angled Up

Reusable, Smooth Tips, 1 Each

Length (L)	8 1/4 in 210 mm	10 1/2 in 267 mm	8 1/4 in 210 mm	10 1/2 in 267 mm
Working (WL)	4 1/4 in 108 mm	6 1/4 in 159 mm	4 1/4 in 108 mm	6 1/4 in 159 mm

TIP SIZE	NON-INSULATED	NON-INSULATED	INSULATED	INSULATED
1.0 mm	10-1445	10-1482	10-1445I	10-1482I
1.5 mm	10-1447	10-1486	10-1447I	10-1486I
2.0 mm	10-1449	10-1488	10-1449I	10-1488I



Bayonet Bipolar Forceps, Irrigating

Reusable, Smooth Tips, Insulated with Stylet, 1 Each

Length (L)	6 1/2 in 165 mm	7 3/4 in 197 mm	8 5/8 in 219 mm	10 1/2 in 267 mm
Working (WL)	2 3/4 in 70 mm	3 5/8 in 92 mm	4 3/8 in 111 mm	6 1/4 in 159 mm


TIP SIZE	INSULATED	INSULATED	INSULATED	INSULATED
0.5 mm		10-0351I	10-0401I	10-0461I
0.75 mm	10-0321I	10-0361I	10-0411I	10-0471I
1.0 mm	10-0331I	10-0371I	10-0421I	10-0481I
1.5 mm	10-0341I	10-0381I	10-0431I	10-0491I
2.0 mm		10-0391I	10-0441I	10-0493I

Single-use products are sold sterile.

Kits Include Forceps, 12 ft (3.6 m) Cord and Surgi-Wipe.

Questions? Contact us at customerservice@symmetrysurgical.com or learn more at symmetrysurgical.com

Olsen®/Ultra™ Reusable Bayonet Bipolar Forceps



Bayonet Bipolar Forceps, with Midas Touch® Non-Stick Tips

Reusable, Smooth Tips, 1 Each

Length (L)	7 ¾ in 197 mm	8 ⅝ in 219 mm	10 ½ in 267 mm
Working (WL)	3 ⅝ in 92 mm	4 ⅜ in 111 mm	6 ¼ in 159 mm
TIP SIZE	INSULATED	INSULATED	INSULATED
0.5 mm	90-7004	90-8012	90-9008
0.75 mm	90-7000	90-8000	90-9007
1.2 mm	90-7005	90-8013	90-9006
1.5 mm	90-7001	90-8001	90-9005
2.0 mm	90-7006	90-8003	90-9009



Cushing Bayonet Bipolar Forceps


Reusable, Smooth Tips, 1 Each

CODE	DESCRIPTION	TIP SIZE	LENGTH
08-0000	Forceps, Non-Insulated	0.7 mm	7 ½ in 191 mm
08-0026	Forceps, Non-Insulated	1.0 mm	7 ½ in 191 mm
08-0027	Forceps, Non-Insulated	1.5 mm	7 ¾ in 197 mm
08-0000I	Forceps, Insulated	0.7 mm	7 ¾ in 197 mm
08-0026I	Forceps, Insulated	1.0 mm	7 ¾ in 197 mm
08-0027I	Forceps, Insulated	1.5 mm	7 ¾ in 197 mm

Cushing Bayonet Bipolar Forceps, Irrigating

Reusable, Smooth Tips, 1 Each

CODE	DESCRIPTION	TIP SIZE	LENGTH
08-0026IG	Irrigating, Insulated	1.0 mm	7 ½ in 191 mm



Gerald Bayonet Bipolar Forceps, Irrigating

Reusable, 1 Each

CODE	DESCRIPTION	TIP SIZE	LENGTH
08-0041	Non-Insulated	1.0 mm	7 ½ in 191 mm
08-0013IG	Insulated	1.0 mm	7 ¼ in 184 mm

Single-use products are sold sterile.
Kits Include Forceps, 12 ft (3.6 m) Cord and Surgi-Wipe.

Olsen®/Ultra™ Reusable Bayonet Bipolar Forceps



Hardy™ Bayonet Bipolar Forceps

Reusable, Smooth Tips, 1 Each

CODE	DESCRIPTION	TIP SIZE	LENGTH
08-0043	Forceps, Non-Insulated	0.5 mm	8 ½ in 216 mm
08-0043I	Forceps, Insulated	0.5 mm	8 ½ in 216 mm

K-Style Standard Acute Angle Bayonet Bipolar Forceps, Straight

Reusable, Smooth Tip, Midas Touch® Non-Stick Tips, 1 Each

Length (L)	8 ⅝ in 219 mm	9 ½ in 241 mm
Working (WL)	4 ½ in 114 mm	5 ⅛ in 130 mm
TIP SIZE	INSULATED	INSULATED
1.0 mm	90-8008	90-9504



K-Style Standard Angle Bayonet Bipolar Forceps, Straight

Reusable, Smooth Tip, Midas Touch® Non-Stick Tips, 1 Each


Length (L)	8 ⅝ in 219 mm	9 ½ in 241 mm	10 ½ in 267 mm
Working (WL)	4 ½ in 114 mm	5 ⅛ in 130 mm	6 ¼ in 159 mm
TIP SIZE	INSULATED	INSULATED	INSULATED
0.5 mm	90-8015	90-9507	
0.75 mm	90-8016	90-9508	
1.0 mm	90-8006	90-9501	90-9001
1.5 mm		90-9500	90-9000

Single-use products are sold sterile.

Kits Include Forceps, 12 ft (3.6 m) Cord and Surgi-Wipe.

Questions? Contact us at customerservice@symmetrysurgical.com
or learn more at symmetrysurgical.com


Olsen®/Ultra™ Reusable Bayonet Bipolar Forceps



K-Style Standard Angle Bayonet Bipolar Forceps, Curved Up

Reusable, Smooth Tip, Midas Touch® Non-Stick Tips, 1 Each


Length (L)	8 ⁵ / ₈ in 219 mm	9 ¹ / ₂ in 241 mm
Working (WL)	4 ¹ / ₂ in 114 mm	5 ¹ / ₈ in 130 mm
TIP SIZE	INSULATED	INSULATED
0.5 mm	90-8017	90-9509
0.75 mm	90-8018	90-9510
1.0 mm	90-8007	90-9502



K-Style Standard Angle Bayonet Bipolar Forceps, Curved Down

Reusable, Smooth Tip, Midas Touch® Non-Stick Tips, 1 Each

Length (L)	8 ⁵ / ₈ in 219 mm	9 ¹ / ₂ in 241 mm
Working (WL)	4 ¹ / ₂ in 114 mm	5 ¹ / ₈ in 130 mm
TIP SIZE	INSULATED	INSULATED
1.0 mm	90-9503	90-9506



K-Style Standard Angle Bayonet Bipolar Forceps, Curved Up

Reusable, Smooth Tip, Midas Touch® Non-Stick Tips, 1 Each

Length (L)	8 ⁵ / ₈ in 219 mm	9 ¹ / ₂ in 241 mm
Working (WL)	4 ¹ / ₂ in 114 mm	5 ¹ / ₈ in 130 mm
TIP SIZE	INSULATED	INSULATED
1.0 mm	90-8009	90-9505

Single-use products are sold sterile.
Kits Include Forceps, 12 ft (3.6 m) Cord and Surgi-Wipe.

Olsen®/Ultra™ Reusable Bayonet Bipolar Forceps



Olsen® Bayonet Bipolar Forceps

Reusable, 1 Each

CODE	DESCRIPTION	TIP SIZE	LENGTH
08-0090	Forceps, Non-Insulated	0.5 mm	5 ½ in 140 mm
08-0091	Forceps, Non-Insulated	1.5 mm	7 in 178 mm
08-0089	Forceps, Non-Insulated	1.5 mm	7 ½ in 191 mm

Olsen® Bayonet Bipolar Forceps, Irrigating

Reusable, 1 Each

CODE	DESCRIPTION	TIP SIZE	LENGTH
08-0223	Forceps, Irrigating	0.25 mm	7 ½ in 191 mm
08-0143	Forceps, Irrigating	0.5 mm	7 ¾ in 197 mm
08-0149	Forceps, Irrigating	1.5 mm	8 ¾ in 222 mm



Rhoton® Bayonet Bipolar Forceps

Reusable, Round Handle, 1 Each

Length (L)	7 ½ in 181 mm	7 ⅞ in 200 mm	8 ⅝ in 219 mm	10 in 254 mm	7 ⅞ in 200 mm	8 ⅝ in 219 mm
------------	---------------	---------------	---------------	--------------	---------------	---------------

TIP SIZE	NON-INSULATED	NON-INSULATED	NON-INSULATED	NON-INSULATED	INSULATED	INSULATED
0.3 mm			08-0308			
0.5 mm		08-0303	08-0309			
0.7 mm			08-0310	08-0319	08-0304I	08-0310I
1.0 mm	08-0300		08-0311			08-0311I
1.5 mm		08-0306	08-0312		08-0306I	08-0312I
2.0 mm	08-0301	08-0307	08-0313		08-0301I	08-0313I

Rhoton® Bayonet Bipolar Forceps, Irrigating

Reusable, Round Handle, 1 Each

CODE	DESCRIPTION	TIP SIZE	LENGTH
08-0313IG	Irrigating Forceps, Insulated	2.0 mm	8 ⅝ in 219 mm

Single-use products are sold sterile.

Kits Include Forceps, 12 ft (3.6 m) Cord and Surgi-Wipe.

Questions? Contact us at customerservice@symmetrysurgical.com
or learn more at symmetrysurgical.com

Olsen®/Ultra™ Reusable Bayonet Bipolar Forceps



Scoville-Greenwood Bayonet Bipolar Forceps


Reusable, 1 Each

CODE	DESCRIPTION	TIP SIZE	LENGTH	
08-0030	Forceps, Non-Insulated	1.5 mm	7 ½ in	191 mm
08-0092	Forceps, Non-Insulated	1.5 mm	7 in	178 mm
08-0028I	Forceps, Insulated	0.25 mm	7 ¾ in	197 mm
08-0030I	Forceps, Insulated	1.5 mm	7 ½ in	191 mm

Scoville-Greenwood Bayonet Bipolar Forceps, Irrigating

Reusable, 1 Each

CODE	DESCRIPTION	TIP SIZE	LENGTH	
08-0030IG	Irrigating, Insulated	1.5 mm	7 ¾ in	197 mm



Titanium Bayonet Bipolar Forceps, Straight Tips

Reusable, Flat Handle, 1 Each

CODE	DESCRIPTION	TIP SIZE	LENGTH	
08-0020	Forceps, Insulated	1.0 mm	8 ½ in	216 mm

Titanium Bayonet Bipolar Forceps, Tips Angled Down, Straight Tips

Reusable, Flat Handle, 1 Each

CODE	DESCRIPTION	TIP SIZE	LENGTH	
08-0104	Forceps, Non-Insulated	0.5 mm	8 ¾ in	222 mm
08-0105I	Forceps, Insulated	0.5 mm	8 ¾ in	222 mm

Titanium Bayonet Bipolar Forceps, Irrigating

Reusable, Flat Handle, 1 Each

CODE	DESCRIPTION	TIP SIZE	LENGTH	
08-0075	Forceps, Insulated	0.25 mm	8 ¾ in	222 mm
08-0100IG	Forceps, Insulated	0.5 mm	8 ¾ in	222 mm
08-0101IG	Forceps, Insulated	1.5 mm	8 ¾ in	222 mm

Single-use products are sold sterile.
Kits Include Forceps, 12 ft (3.6 m) Cord and Surgi-Wipe.

Olsen®/Ultra™ Reusable Bayonet Bipolar Forceps

Titanium Bayonet Bipolar Forceps, Straight Tips

Reusable, 1 Each

CODE	DESCRIPTION	TIP SIZE	LENGTH	
08-0100	Forceps, Non-Insulated	0.5 mm	8 ¾ in	222 mm
08-0101	Forceps, Non-Insulated	1.5 mm	8 ¾ in	222 mm
08-0100I	Forceps, Insulated	0.5 mm	8 ¾ in	222 mm
08-0101I	Forceps, Insulated	1.5 mm	8 ¾ in	222 mm

Titanium Bayonet Bipolar Forceps, with Midas Touch® Non-Stick Tips

Reusable, Straight, 1 Each

CODE	DESCRIPTION	TIP SIZE	LENGTH	
08-0062NS	Forceps, Non-Insulated	0.25 mm	8 ¾ in	222 mm



Titanium Rhoton® Bayonet Bipolar Forceps

Reusable, Round Handle, 1 Each

CODE	DESCRIPTION	TIP SIZE	LENGTH	
08-0122	Forceps, Non-Insulated	0.5 mm	9 ½ in	241 mm
08-0123	Forceps, Non-Insulated	1.0 mm	9 ½ in	241 mm
08-0122I	Forceps, Insulated	0.5 mm	9 ½ in	241 mm
08-0120I	Forceps, Insulated	1.5 mm	8 ½ in	216 mm

Titanium Rhoton® Bayonet Bipolar Forceps, Straight

Reusable, Round Handle, 1 Each

CODE	DESCRIPTION	TIP SIZE	LENGTH	
08-0058	Titanium Forceps	0.5 mm	6 in	152 mm
08-0130	Titanium Forceps	0.5 mm	7 in	178 mm
08-0509	Titanium Forceps	1.5 mm	8 7/8 in	225 mm

Titanium Rhoton® Bayonet Bipolar Forceps, Curved

Reusable, Round Handle, 1 Each

CODE	DESCRIPTION	TIP SIZE	LENGTH	
08-0061	Titanium Forceps	0.5 mm	7 1/8 in	181 mm

Single-use products are sold sterile.

Kits Include Forceps, 12 ft (3.6 m) Cord and Surgi-Wipe.

Questions? Contact us at customerservice@symmetrysurgical.com
or learn more at symmetrysurgical.com

Olsen®/Ultra™ Reusable Bayonet Bipolar Forceps



Yasargil Bayonet Bipolar Forceps

Reusable, Flat Handle, 1 Each

CODE	DESCRIPTION	TIP SIZE	LENGTH
08-0015	Forceps, Non-Insulated	0.7 mm	6¼ in 159 mm
08-0016	Forceps, Non-Insulated	0.7 mm	6¼ in 159 mm
08-0017	Forceps, Non-Insulated	0.5 mm	7⅛ in 181 mm
08-0018	Forceps, Non-Insulated	0.7 mm	7⅛ in 181 mm
08-0019	Forceps, Non-Insulated	1.0 mm	7⅛ in 181 mm
08-0024	Forceps, Non-Insulated	1.3 mm	7⅛ in 181 mm
08-0022	Forceps, Non-Insulated	0.7 mm	7⅞ in 200 mm
08-0025	Forceps, Non-Insulated	1.3 mm	7⅞ in 200 mm
08-0023I	Forceps, Insulated	0.7 mm	7⅞ in 200 mm
08-0025I	Forceps, Insulated	1.3 mm	7⅞ in 200 mm



Yasargil Bayonet Bipolar Forceps, Irrigating

Reusable, Flat Handle, 1 Each

CODE	DESCRIPTION	TIP SIZE	LENGTH
08-0022IG	Forceps, Insulated	0.7 mm	7⅞ in 181 mm
08-0023IG	Forceps, Insulated	1.0 mm	7⅞ in 200 mm


Single-use products are sold sterile.
Kits Include Forceps, 12 ft (3.6 m) Cord and Surgi-Wipe.

Olsen®/Ultra™ Reusable Bayonet Bipolar Forceps

Adson Bipolar Forceps, Smooth Tip

Reusable, Smooth Tips, 1 Each

CODE	DESCRIPTION	TIP SIZE	LENGTH
08-0004	Forceps, Non-Insulated	1.0 mm	4 ¾ in 121 mm
10-1082	Forceps, Non-Insulated	1.0 mm	5 in 127 mm
08-0004I	Forceps, Insulated	1.0 mm	4 ¾ in 121 mm
10-1082I	Forceps, Insulated	1.0 mm	5 in 127 mm



Adson Bipolar Forceps, Smooth Tip, with Midas Touch® Non-Stick Tips

Reusable, 1 Each


CODE	DESCRIPTION	TIP SIZE	LENGTH
90-2001	Forceps, Insulated	0.5 mm	5 in 127 mm
90-2000	Forceps, Insulated	1.0 mm	5 in 127 mm



Adson Bipolar Forceps, Serrated Tip

Reusable, 1 Each


CODE	DESCRIPTION	TIP SIZE	LENGTH
10-1081	Forceps, Non-Insulated	1.0 mm	5 in 127 mm
10-1081I	Forceps, Insulated	1.0 mm	5 in 127 mm



Adson Bipolar Forceps, Irrigating

Reusable, 1 Each

CODE	DESCRIPTION	TIP SIZE	LENGTH
10-0201I	Insulated, with Stylet	0.5 mm	5 in 127 mm
10-0200I	Insulated, with Stylet	1.0 mm	5 in 127 mm



Single-use products are sold sterile.
Kits Include Forceps, 12 ft (3.6 m) Cord and Surgi-Wipe.


Questions? Contact us at customerservice@symmetrysurgical.com
or learn more at symmetrysurgical.com

Olsen®/Ultra™ Reusable Bipolar Forceps

Coaption Bipolar Forceps

Reusable, 1 Each

CODE	DESCRIPTION	TIP SIZE	LENGTH
08-0008	Forceps, Non-Insulated	Micro	4¼ in 108 mm
08-0048	Forceps, Non-Insulated	0.4 mm	4¼ in 108 mm
10-1011	Forceps, Non-Insulated	Micro	4½ in 114 mm
10-1011I	Forceps, Insulated	Micro	4½ in 114 mm



Cushing Bipolar Forceps

Reusable, Smooth Tips, 1 Each

CODE	DESCRIPTION	TIP SIZE	LENGTH
10-1160	Forceps, Non-Insulated	1.0 mm	7¼ in 184 mm
10-1162	Forceps, Non-Insulated	1.5 mm	7¼ in 184 mm
08-0010	Forceps, Non-Insulated	1.5mm	7⅛ in 181 mm
10-1160I	Forceps, Insulated	1.0 mm	7¼ in 184 mm
10-1162I	Forceps, Insulated	1.5 mm	7¼ in 184 mm



Gerald Bipolar Forceps


Reusable, Smooth Tips, 1 Each

CODE	DESCRIPTION	TIP SIZE	LENGTH
08-0012	Forceps, Non-Insulated	1.0 mm	7 in 178 mm
10-1180	Forceps, Non-Insulated	1.0 mm	7¼ in 184 mm
08-0013	Forceps, Non-Insulated	1.0 mm	7¼ in 184 mm
10-1180I	Forceps, Insulated	1.0 mm	7¼ in 184 mm



Replacement Stylets for Irrigating Bipolar Forceps

CODE	DESCRIPTION	LENGTH	QTY
08-1130	Accessory		12/Box
60-0300	Accessory	10¼ in 260 mm	5/Box



Single-use products are sold sterile.
Kits Include Forceps, 12 ft (3.6 m) Cord and Surgi-Wipe.

Olsen®/Ultra™ Reusable Bipolar Forceps

Jewelers Bipolar Forceps

Reusable, Insulated, Smooth Tips, 1 Each


CODE	DESCRIPTION	TIP SIZE	LENGTH
80004	Forceps, Straight	Micro	4 ¼ in 108 mm
80003	Forceps, Straight	0.5 mm	4 ¼ in 108 mm
80002	Forceps, Straight	1.0 mm	4 ¼ in 108 mm
80001	Forceps, Straight	1.5 mm	4 ¼ in 108 mm
80000	Forceps, Straight	2.0 mm	4 ¼ in 108 mm
80104	Forceps, Curved	Micro	5 in 127 mm
80103	Forceps, Curved	0.5 mm	5 in 127 mm
80102	Forceps, Curved	1.0 mm	5 in 127 mm
80101	Forceps, Curved	1.5 mm	5 in 127 mm
80100	Forceps, Curved	2.0 mm	5 in 127 mm



Jewelers Bipolar Forceps

Reusable, Non-Insulated, Smooth Tips, 1 Each

CODE	DESCRIPTION	TIP SIZE	LENGTH
10-1061	Forceps, Straight	Micro	4 ½ in 114 mm
10-1062	Forceps, Straight	Blunt	4 ½ in 114 mm
08-0002	Forceps, Curved	0.4 mm	4 in 102 mm
08-0005	Forceps, Curved	0.4 mm	4 in 102 mm



Jewelers Bipolar Forceps

Reusable, Insulated, Smooth Tips, 1 Each

CODE	DESCRIPTION	TIP SIZE	LENGTH
10-1061I	Forceps, Straight	Micro	4 ½ in 114 mm
10-1062I	Forceps, Straight	Blunt	4 ½ in 114 mm
08-0001I	Forceps, Straight	0.4 mm	4 ½ in 114 mm



Jewelers Bipolar Forceps, with Midas Touch® Non-Stick Tips

Reusable, Smooth Tips, 1 Each

CODE	DESCRIPTION	TIP SIZE	LENGTH
90-1000	Forceps, Insulated	Micro	4 ⅝ in 117 mm



Single-use products are sold sterile.
Kits Include Forceps, 12 ft (3.6 m) Cord and Surgi-Wipe.

Questions? Contact us at customerservice@symmetrysurgical.com
or learn more at symmetrysurgical.com

Olsen®/Ultra™ Reusable Bipolar Forceps

Jewelers Bipolar Forceps, Ridged Handle, Angled

Reusable, Smooth Tips, 1 Each

CODE	DESCRIPTION	TIP SIZE	LENGTH
10-1202	Forceps, Non-Insulated	0.75 mm	4 ½ in 114 mm
10-1202I	Forceps, Insulated	0.75 mm	4 ½ in 114 mm

Jewelers Bipolar Forceps, Ridged Handle, Straight

Reusable, Smooth Tips, 1 Each

CODE	DESCRIPTION	TIP SIZE	LENGTH
10-1203	Forceps, Non-Insulated	Micro	3 ¾ in 99 mm
10-1203I	Forceps, Insulated	Micro	3 ¾ in 99 mm

McPherson Bipolar Forceps, Angled

Reusable, Smooth Tips, 1 Each

CODE	DESCRIPTION	TIP SIZE	LENGTH
08-0007	Forceps, Non-Insulated	0.5 mm	3 ¾ in 95 mm
10-1200	Forceps, Non-Insulated	0.5 mm	3 ¾ in 99 mm
10-1021	Forceps, Non-Insulated	0.5 mm	4 in 102 mm
08-0007I	Forceps, Insulated	0.5 mm	3 ¾ in 95 mm
10-1200I	Forceps, Insulated	0.5 mm	3 ¾ in 99 mm
10-1021I	Forceps, Insulated	0.5 mm	4 in 102 mm

McPherson Bipolar Forceps, Straight


Reusable, Smooth Tips, 1 Each

CODE	DESCRIPTION	TIP SIZE	LENGTH
10-1201	Forceps, Non-Insulated	0.5 mm	3 ¾ in 99 mm
10-1201I	Forceps, Insulated	0.5 mm	3 ¾ in 99 mm

Potts-Smith Bipolar Forceps

Reusable, Smooth Tips, 1 Each

CODE	DESCRIPTION	TIP SIZE	LENGTH
80404	Forceps, Insulated	Micro	7 ½ in 191 mm
80403	Forceps, Insulated	0.5 mm	7 ½ in 191 mm
80402	Forceps, Insulated	1.0 mm	7 ½ in 191 mm
80401	Forceps, Insulated	1.5 mm	7 ½ in 191 mm
80400	Forceps, Insulated	2.0 mm	7 ½ in 191 mm



Single-use products are sold sterile.
Kits Include Forceps, 12 ft (3.6 m) Cord and Surgi-Wipe.

Olsen®/Ultra™ Reusable Bipolar Forceps

Semkin Bipolar Forceps

Reusable, Smooth Tips, 1 Each

CODE	DESCRIPTION	TIP SIZE	LENGTH
10-1091	Forceps, Non-Insulated	0.75 mm	5 ¾ in 146 mm
10-1191	Forceps, Non-Insulated	1.0 mm	5 ¾ in 146 mm
10-1091I	Forceps, Insulated	0.75 mm	5 ¾ in 146 mm
10-1191I	Forceps, Insulated	1.0 mm	5 ¾ in 146 mm



Semkin Bipolar Forceps, with Midas Touch® Non-Stick Tips

Reusable, Smooth Tips, 1 Each

CODE	DESCRIPTION	TIP SIZE	LENGTH
90-3000	Forceps, Insulated	0.75 mm	5 ¾ in 146 mm



Titanium Jewelers Bipolar Forceps

Reusable, 1 Each

CODE	DESCRIPTION	TIP SIZE	LENGTH
08-0052	30° Tip	0.4 mm	3 ½ in 89 mm
08-0056	90° Tip	0.4 mm	4 ¼ in 108 mm

Titanium Jewelers Bipolar Forceps, Curved

Reusable, 1 Each

CODE	DESCRIPTION	TIP SIZE	LENGTH
08-0053	Titanium Forceps	0.4 mm	3 ½ in 89 mm

Titanium Jewelers Bipolar Forceps, Straight

Reusable, 1 Each


CODE	DESCRIPTION	TIP SIZE	LENGTH
08-0054	Titanium Forceps	0.4 mm	3 ½ in 89 mm
08-0057	Titanium Forceps	0.4 mm	4 ¼ in 108 mm

Single-use products are sold sterile.

Kits Include Forceps, 12 ft (3.6 m) Cord and Surgi-Wipe.

Questions? Contact us at customerservice@symmetrysurgical.com
or learn more at symmetrysurgical.com

Olsen® Single-Use Monopolar Forceps



Monopolar Forceps


Single-Use, Insulated, Serrated Tips, Uses 960 Series Adapter, 10 per Box

Length (L)	8 in 203 mm	8 ½ in 216 mm	8 ¼ in 210 mm	7 ½ in 191 mm
Working (WL)		4 ⅛ in 105 mm	4 ¼ in 108 mm	3 in 76 mm
Tip Style	Straight	Straight	Straight	Cathedral
Forceps Style	Potts-Smith	Bayonet	Bayonet	Bayonet

TIP SIZE	FOOT-ACTIVATED ¹	FOOT-ACTIVATED ¹	HAND-CONTROLLED ²	HAND-CONTROLLED ²
1.5 mm	20-2650K	20-2920K	20-2920HK	20-2960HK

¹Foot-activated forceps include a 10 ft (3.0m) attached cord with crimped wire generator end.
²Hand-controlled forceps include 10 ft (3.0m) attached cord with molded generator plug.


Olsen® Reusable Bayonet Monopolar Forceps



Bayonet Monopolar Forceps, Foot-Activated, Straight

Reusable, Threaded Connector, Serrated Tips, 1 Each

Length (L)	6 ³ / ₄ in 171 mm	7 ³ / ₄ in 197 mm	8 in 203 mm	8 ¹ / ₂ in 216 mm
Working (WL)	3 ¹ / ₈ in 79 mm	3 ⁷ / ₈ in 99 mm	3 in 76 mm	4 ³ / ₄ in 121 mm
TIP SIZE	INSULATED	INSULATED	INSULATED	INSULATED
0.5 mm	41103	41203	41403	41605
1.0 mm	41102	41202	41402	41604
1.5 mm	41101	41201	41401	41602
2.0 mm	41100	41200	41400	41601



Bayonet Monopolar Forceps, Hand-Controlled

Reusable, Serrated Tips, Insulated with 10 ft (3.0 m) Attached Cord, 1 Each


Length (L)	6 ¹ / ₂ in 165 mm	8 ¹ / ₄ in 210 mm	7 ¹ / ₂ in 191 mm	7 ¹ / ₂ in 191 mm
Working (WL)	3 ¹ / ₄ in 83 mm	3 ¹ / ₂ in 89 mm	3 ³ / ₄ in 95 mm	3 ³ / ₄ in 95 mm
TIP SIZE	SERRATED TIP	SERRATED TIP	STRAIGHT	CATHEDRAL
0.5 mm	11104	11605	11206	11213
1.0 mm	11102	11604	11205	11212
1.5 mm	11101	11602	11204	11210
2.0 mm	11100	11600	11200	11208

Single-use products are sold sterile.

Kits Include Forceps, 12 ft (3.6 m) Cord and Surgi-Wipe.

Questions? Contact us at customerservice@symmetrysurgical.com or learn more at symmetrysurgical.com

Olsen® Reusable Bayonet Monopolar Forceps



Bayonet Monopolar Forceps, Foot-Activated, Curved Up

Reusable, Threaded Connector, Serrated Tips, 1 Each

Length (L)	7 ¾ in 197 mm	8 ½ in 216 mm
Working (WL)	3 7/8 in 99 mm	4 ¾ in 121 mm
TIP SIZE	INSULATED	INSULATED
0.5 mm	41302	41703
1.0 mm	41304	41702
1.5 mm	41301	41701
2.0 mm	41300	41700

Product Complaint? Call 800-251-3000

Patient safety is our #1 priority. We are committed to delivering exceptional quality and service. Compliance with domestic and international medical device regulations is strictly maintained at all times.

Product complaints are addressed promptly by our Complaint Handling Team. Contact Customer Service and request a Return Material Authorization (RMA) to begin the investigation process. Please return the instrument (with the original package and label – if available) along with the RMA form.

A representative will contact you within 48 hours of receiving the returned instrument with details regarding the steps we will take to ensure you and your patients are satisfied with the instruments you purchase from us.



Olsen® Reusable Monopolar Forceps

Adson Monopolar Forceps, Hand-Controlled

Reusable, Serrated Tips, Insulated, 1 Each

CODE	DESCRIPTION	TIP SIZE	LENGTH
12204	10 ft (3.0 m) Attached Cord	Micro	4 ¾ in 121 mm
12203	10 ft (3.0 m) Attached Cord	0.5 mm	4 ¾ in 121 mm
12202	10 ft (3.0 m) Attached Cord	1.0 mm	4 ¾ in 121 mm
12200	10 ft (3.0 m) Attached Cord	2.0 mm	4 ¾ in 121 mm



Adson Monopolar Forceps, Foot-Controlled

Reusable, Threaded Connector, Serrated Tips, Insulated, 1 Each

CODE	DESCRIPTION	TIP SIZE	LENGTH
42505	Monopolar Forceps, Insulated	0.5 mm	5 ¼ in 133 mm
42503	Monopolar Forceps, Insulated	1.0 mm	5 ¼ in 133 mm
42502	Monopolar Forceps, Insulated	1.5 mm	5 ¼ in 133 mm
42500	Monopolar Forceps, Insulated	2.0 mm	5 ¼ in 133 mm



Cushing Monopolar Forceps, Hand-Controlled

Reusable, Serrated Tips, Insulated, 1 Each

CODE	DESCRIPTION	TIP SIZE	LENGTH
12103	10 ft (3.0 m) Attached Cord	0.5 mm	7 in 178 mm
12102	10 ft (3.0 m) Attached Cord	1.0 mm	7 in 178 mm
12101	10 ft (3.0 m) Attached Cord	1.5 mm	7 in 178 mm
12100	10 ft (3.0 m) Attached Cord	2.0 mm	7 in 178 mm



Jewelers Monopolar Forceps, Foot-Activated

Reusable, Threaded Connector, Serrated Tips, Insulated, 1 Each

CODE	DESCRIPTION	TIP SIZE	LENGTH
43104	Insulated, Straight	Micro	4 ¾ in 121 mm
43110	Insulated, Straight	Micro	4 ¾ in 121 mm
43103	Insulated, Straight	0.5 mm	4 ¾ in 121 mm
43109	Insulated, Straight	0.5 mm	4 ¾ in 121 mm
43204	Insulated, Curved	Micro	4 ¾ in 121 mm
43203	Insulated, Curved	0.5 mm	4 ¾ in 121 mm



Single-use products are sold sterile.
Kits Include Forceps, 12 ft (3.6 m) Cord and Surgi-Wipe.

Questions? Contact us at customerservice@symmetrysurgical.com
or learn more at symmetrysurgical.com

Olsen® Reusable Monopolar Forceps

Potts-Smith Monopolar Forceps, Hand-Controlled

Reusable, Serrated Tips, Insulated, 1 Each

CODE	DESCRIPTION	TIP SIZE	LENGTH
11804	10 ft (3.0 m) Attached Cord	0.5 mm	7 in 178 mm
11803	10 ft (3.0 m) Attached Cord	1.0 mm	7 in 178 mm
11802	10 ft (3.0 m) Attached Cord	1.5 mm	7 in 178 mm
11800	10 ft (3.0 m) Attached Cord	2.0 mm	7 in 178 mm
11906	10 ft (3.0 m) Attached Cord	0.5 mm	8 ¼ in 210 mm
11904	10 ft (3.0 m) Attached Cord	1.0 mm	8 ¼ in 210 mm
11903	10 ft (3.0 m) Attached Cord	1.5 mm	8 ¼ in 210 mm
11901	10 ft (3.0 m) Attached Cord	2.0 mm	8 ¼ in 210 mm
12000	10 ft (3.0 m) Attached Cord	2.0 mm	9 ½ in 241 mm
12003	10 ft (3.0 m) Attached Cord	1.5 mm	9 ½ in 241 mm
12004	10 ft (3.0 m) Attached Cord	1.0 mm	9 ½ in 241 mm
12005	10 ft (3.0 m) Attached Cord	0.5 mm	9 ½ in 241 mm



Potts-Smith Monopolar Forceps, Foot-Controlled

Reusable, Serrated Tips, Insulated, 1 Each

CODE	TIP SIZE	LENGTH
42203	0.5 mm	7 ½ in 191 mm
42202	1.0 mm	7 ½ in 191 mm
42201	1.5 mm	7 ½ in 191 mm
42200	2.0 mm	7 ½ in 191 mm
42304	0.5 mm	8 ½ in 216 mm
42302	1.0 mm	8 ½ in 216 mm
42301	1.5 mm	8 ½ in 216 mm
42300	2.0 mm	8 ½ in 216 mm
42405	0.5 mm	9 ¾ in 248 mm
42404	1.0 mm	9 ¾ in 248 mm
42402	1.5 mm	9 ¾ in 248 mm
42400	2.0 mm	9 ¾ in 248 mm



Single-use products are sold sterile.
Kits Include Forceps, 12 ft (3.6 m) Cord and Surgi-Wipe.

Olsen® Monopolar Electrodes



1/8 in (3.2 mm) Ball



3/16 in (4.8 mm) Ball

Ball Electrodes, 1/8 in (3.2 mm) Ball

Reusable, Straight, 1 Each

CODE	DESCRIPTION	LENGTH	
66200	Straight Ball Electrode	2 in	51 mm
66214	Straight Ball Electrode	3 in	76 mm
66202	Straight Ball Electrode	4 in	102 mm
66204	Straight Ball Electrode	5 in	127 mm
66212	Straight Ball Electrode	5 ½ in	140 mm
66206	Straight Ball Electrode	6 in	152 mm
66208	Straight Ball Electrode	7 in	178 mm
66210	Straight Ball Electrode	8 in	203 mm

Ball Electrodes, 3/16 in (4.8 mm) Ball

Reusable, Straight, 1 Each

CODE	DESCRIPTION	LENGTH	
66201	Straight Ball Electrode	2 in	51 mm
66215	Straight Ball Electrode	3 in	76 mm
66203	Straight Ball Electrode	4 in	102 mm
66205	Straight Ball Electrode	5 in	127 mm
66213	Straight Ball Electrode	5 ½ in	140 mm
66207	Straight Ball Electrode	6 in	152 mm
66209	Straight Ball Electrode	7 in	178 mm
66211	Straight Ball Electrode	8 in	203 mm

Ball Electrodes

Single-Use, Straight

CODE	DESCRIPTION	LENGTH		QTY
97200	Straight Ball Electrode	2 ¾ in	70 mm	50
97201	Straight Ball Electrode	5 in	127 mm	25

Ball Electrodes, with Midas Touch® Non-Stick

Single-Use, Straight

CODE	DESCRIPTION	LENGTH		QTY
30-2001	Straight Ball Electrode	2 ¾ in	70 mm	10
30-2003	Straight Ball Electrode	5 in	127 mm	10

Blade Electrodes, Straight

Reusable, Standard 10 mm Exposure, 1 Each

CODE	DESCRIPTION	LENGTH	
66000	Straight Electrode	2 ¾ in	70 mm
66001	Straight Electrode	4 in	102 mm
66002	Straight Electrode	5 in	127 mm
66003	Straight Electrode	6 in	152 mm
66004	Straight Electrode	7 in	178 mm
66005	Straight Electrode	8 in	203 mm

Single-use products are sold sterile.

Kits Include Forceps, 12 ft (3.6 m) Cord and Surgi-Wipe.

Questions? Contact us at customerservice@symmetrysurgical.com or learn more at symmetrysurgical.com

Olsen® Monopolar Electrodes

Blade Electrodes, Angled

Reusable, Standard 10 mm Exposure, 1 Each

CODE	DESCRIPTION	LENGTH	
66006	Angled Electrode	2 ³ / ₄ in	70 mm
66007	Angled Electrode	4 in	102 mm
66008	Angled Electrode	5 in	127 mm
66009	Angled Electrode	6 in	152 mm
66010	Angled Electrode	7 in	178 mm
66011	Angled Electrode	8 in	203 mm
66014	Angled Electrode	11 in	279 mm



Blade Electrodes, Short/Close-Up

Reusable, Standard 10 mm Exposure, 1 Each

CODE	DESCRIPTION	LENGTH	
66700	Stainless Steel Standard Blade	2 in	51 mm
66704	Stainless Steel Blunt Blade	2 in	51 mm
66703	Stainless Steel Modified Blade	2 ¹ / ₂ in	64 mm
66705	Titanium Standard Blade	2 ¹ / ₂ in	64 mm



Blade Electrodes, Straight

Reusable, 1 Each

CODE	DESCRIPTION	LENGTH		EXPOSURE
68000	Straight Electrode	2 ³ / ₄ in	70 mm	2.0 mm
68001	Straight Electrode	2 ³ / ₄ in	70 mm	3.0 mm
68003	Straight Electrode	2 ³ / ₄ in	70 mm	5.0 mm
68002	Straight Electrode	4 in	102 mm	2.0 mm
68004	Straight Electrode	4 in	102 mm	6.0 mm
68006	Straight Electrode	4 ¹ / ₂ in	114 mm	2.0 mm
68005	Straight Electrode	6 in	152 mm	2.0 mm



Single-use products are sold sterile.
Kits Include Forceps, 12 ft (3.6 m) Cord and Surgi-Wipe.

Olsen® Monopolar Electrodes

Nasal Needle Electrode

Reusable, 1 Each

CODE	DESCRIPTION	LENGTH
68603	6.0 mm Exposure	8 in 203 mm
68601	15 mm Exposure	8 in 203 mm
68602	50 mm Exposure	8 in 203 mm
68600	102 mm Exposure	8 in 203 mm



Laparoscopic Monopolar Electrode

Reusable

CODE	DESCRIPTION	LENGTH
90-9522	L-Hook	12½in 320 mm
90-9523	Spatula	12½in 320 mm
90-9525	J-Hook	12½in 320 mm
90-9526	Needle	12½in 320 mm
90-9527	Button Tip	12½in 320 mm



Laparoscopic Monopolar Suction/Irrigation Electrode

Reusable, 5.0 mm Shaft, with Valve, Fully Insulated

CODE	DESCRIPTION	LENGTH
90-4005	J-Hook	12½in 320 mm
90-4010	Button Tip	12½in 320 mm
90-4012	L-Hook	12½in 320 mm
90-4015	Spatula	12½in 320 mm

Single-use products are sold sterile.
Kits Include Forceps, 12 ft (3.6 m) Cord and Surgi-Wipe.

Questions? Contact us at customerservice@symmetrysurgical.com
or learn more at symmetrysurgical.com

Monopolar Electrodes

Rectal Snare

Insulated

CODE	DESCRIPTION	LENGTH
50-9731	Monopolar Instrument, Bayonet Shape	11 ³ / ₄ in 298 mm



Wire Loop for Rectal Snare

CODE	DESCRIPTION	LENGTH
50-9749	Loop Electrode	11 ³ / ₄ in 298 mm



Claw Grasper

Direct Drive, Monopolar Laparoscopy Forceps, 5.0 mm , Insulated

CODE	DESCRIPTION	LENGTH
94-1005L-26	Locking Handle	10 ¹ / ₄ in 26 cm
94-1005L	Locking Handle	13 in 33 cm
94-1005L-45	Locking Handle	17 ³ / ₄ in 45 cm

Dolphin Nose Dissector

Direct Drive, Monopolar Laparoscopy Forceps, 5.0 mm, Insulated

CODE	DESCRIPTION	LENGTH
94-1030-26	Non-Locking Handle	10 ¹ / ₄ in 26 cm
94-1030	Non-Locking Handle	13 in 33 cm
94-1030-45	Non-Locking Handle	17 ³ / ₄ in 45 cm

Single-use products are sold sterile.
Kits Include Forceps, 12 ft (3.6 m) Cord and Surgi-Wipe.

Ultra™ Monopolar Laparoscopic Forceps

Grasper

Direct Drive, Monopolar Laparoscopy Forceps, 5.0 mm, Insulated

CODE	DESCRIPTION	LENGTH
94-1032L-26	Locking Handle	10 ¼ in 26 cm
94-1032L	Locking Handle	13 in 33 cm
94-1032L-45	Locking Handle	17 ¾ in 45 cm
94-1032-26	Non-Locking Handle	10 ¼ in 26 cm
94-1032	Non-Locking Handle	13 in 33 cm
94-1032-45	Non-Locking Handle	17 ¾ in 45 cm



Bullet Nose Dissector

Monopolar Laparoscopy Forceps, 5.0 mm, Insulated

CODE	DESCRIPTION	LENGTH
94-1033L-26	Locking Handle	10 ¼ in 26 cm
94-1033L	Locking Handle	13 in 33 cm
94-1033L-45	Locking Handle	17 ¾ in 45 cm
94-1033-26	Non-Locking Handle	10 ¼ in 26 cm
94-1033	Non-Locking Handle	13 in 33 cm
94-1033-45	Non-Locking Handle	17 ¾ in 45 cm



Dolphin Nose Grasper

Monopolar Laparoscopy Forceps, 5.0 mm, Insulated

CODE	DESCRIPTION	LENGTH
94-1034-26	Non-Locking Handle	10 ¼ in 26 cm
94-1034	Non-Locking Handle	13 in 33 cm
94-1034-45	Non-Locking Handle	17 ¾ in 45 cm



Single-use products are sold sterile.
Kits Include Forceps, 12 ft (3.6 m) Cord and Surgi-Wipe.

Questions? Contact us at customerservice@symmetrysurgical.com
or learn more at symmetrysurgical.com

Ultra™ Monopolar Laparoscopic Forceps

Cup Tooth Grasper

Direct Drive, Monopolar Laparoscopy Forceps, 5.0 mm, Insulated

CODE	DESCRIPTION	LENGTH	
94-1035L-26	Locking Handle	10 ¼ in	26 cm
94-1035L	Locking Handle	13 in	33 cm
94-1035L-45	Locking Handle	17 ¾ in	45 cm
94-1035-26	Non-Locking Handle	10 ¼ in	26 cm
94-1035	Non-Locking Handle	13 in	33 cm
94-1035-45	Non-Locking Handle	17 ¾ in	45 cm
94-1035SL-26	Slidelock Handle	10 ¼ in	26 cm
94-1035SL	Slidelock Handle	13 in	33 cm
94-1035SL-45	Slidelock Handle	17 ¾ in	45 cm



Micro Fenestrated Grasper

Direct Drive, Monopolar Laparoscopy Forceps, 5.0 mm, Insulated

CODE	DESCRIPTION	LENGTH	
94-1044-26	Non-Locking Handle	10 ¼ in	26 cm
94-1044	Non-Locking Handle	13 in	33 cm
94-1044-45	Non-Locking Handle	17 ¾ in	45 cm



Needle Nose Dissector

Direct Drive, Monopolar Laparoscopy Forceps, 5.0 mm, Insulated

CODE	DESCRIPTION	LENGTH	
94-1045-26	Non-Locking Handle	10 ¼ in	26 cm
94-1045	Non-Locking Handle	13 in	33 cm
94-1045-45	Non-Locking Handle	17 ¾ in	45 cm

Ultra™ Monopolar Laparoscopic Forceps

Fundus Tri-Max Grasper

Monopolar Laparoscopy Forceps, 5.0 mm, Insulated

CODE	DESCRIPTION	LENGTH
94-1050L-26	Locking Handle	10 ¼ in 26 cm
94-1050L	Locking Handle	13 in 33 cm
94-1050L-45	Locking Handle	17 ¾ in 45 cm

Wave Grasper

Direct Drive, Monopolar Laparoscopy Forceps, 5.0 mm, Insulated

CODE	DESCRIPTION	LENGTH
94-1051L-26	Locking Handle	10 ¼ in 26 cm
94-1051L	Locking Handle	13 in 33 cm
94-1051L-45	Locking Handle	17 ¾ in 45 cm
94-1051-26	Non-Locking Handle	10 ¼ in 26 cm
94-1051	Non-Locking Handle	13 in 33 cm
94-1051-45	Non-Locking Handle	17 ¾ in 45 cm



Alligator Grasper

Direct Drive, Monopolar Laparoscopy Forceps, 5.0 mm, Insulated

CODE	DESCRIPTION	LENGTH
94-1052L-26	Locking Handle	10 ¼ in 26 cm
94-1052L	Locking Handle	13 in 33 cm
94-1052L-45	Locking Handle	17 ¾ in 45 cm

Single-use products are sold sterile.
Kits Include Forceps, 12 ft (3.6 m) Cord and Surgi-Wipe.

Questions? Contact us at customerservice@symmetrysurgical.com
or learn more at symmetrysurgical.com

Ultra™ Monopolar Laparoscopic Forceps

Atraumatic Grasper

Direct Drive, Monopolar Laparoscopy Forceps, 5.0 mm, Insulated

CODE	DESCRIPTION	LENGTH
94-1053L-26	Locking Handle	10 ¼ in 26 cm
94-1053L	Locking Handle	13 in 33 cm
94-1053L-45	Locking Handle	17 ¾ in 45 cm
94-1053SL-26	Slidelock Handle	10 ¼ in 26 cm
94-1053SL	Slidelock Handle	13 in 33 cm
94-1053SL-45	Slidelock Handle	17 ¾ in 45 cm



Retraction Grasper

Monopolar Laparoscopy Forceps, 5.0 mm, Insulated

CODE	DESCRIPTION	LENGTH
94-1056L-26	Locking Handle	10 ¼ in 26 cm
94-1056L	Locking Handle	13 in 33 cm
94-1056L-45	Locking Handle	17 ¾ in 45 cm
94-1056-26	Non-Locking Handle	10 ¼ in 26 cm
94-1056	Non-Locking Handle	13 in 33 cm
94-1056-45	Non-Locking Handle	17 ¾ in 45 cm



Right Angle Dissector

Direct Drive, Monopolar Laparoscopy Forceps, 5.0 mm, Insulated

CODE	DESCRIPTION	LENGTH
94-1063-26	Non-Locking Handle	10 ¼ in 26 cm
94-1063	Non-Locking Handle	13 in 33 cm
94-1063-45	Non-Locking Handle	17 ¾ in 45 cm



Ultra™ Monopolar Laparoscopic Forceps

Soft Curve Dissector

Direct Drive, Monopolar Laparoscopy Forceps, 5.0 mm, Insulated

CODE	DESCRIPTION	LENGTH	
94-1064-26	Non-Locking Handle	10 ¼ in	26 cm
94-1064	Non-Locking Handle	13 in	33 cm
94-1064-45	Non-Locking Handle	17 ¾ in	45 cm



Mixer Right Angle Dissector

Direct Drive, Monopolar Laparoscopy Forceps, 5.0 mm, Insulated

CODE	DESCRIPTION	LENGTH	
94-1066-26	Non-Locking Handle	10 ¼ in	26 cm
94-1066	Non-Locking Handle	13 in	33 cm
94-1066-45	Non-Locking Handle	17 ¾ in	45 cm



Maryland Dissector

Direct Drive, Monopolar Laparoscopy Forceps, 5.0 mm, Insulated

CODE	DESCRIPTION	LENGTH	
94-1067L-26	Locking Handle	10 ¼ in	26 cm
94-1067L	Locking Handle	13 in	33 cm
94-1067L-45	Locking Handle	17 ¾ in	45 cm
94-1067-26	Non-Locking Handle	10 ¼ in	26 cm
94-1067	Non-Locking Handle	13 in	33 cm
94-1067-45	Non-Locking Handle	17 ¾ in	45 cm

Single-use products are sold sterile.
Kits Include Forceps, 12 ft (3.6 m) Cord and Surgi-Wipe.

Questions? Contact us at customerservice@symmetrysurgical.com
or learn more at symmetrysurgical.com

Ultra™ Monopolar Laparoscopic Forceps

Cobra Tooth Grasper

Direct Drive, Monopolar Laparoscopy Forceps, 5.0 mm, Insulated

CODE	DESCRIPTION	LENGTH
94-1069L-26	Locking Handle	10 ¼ in 26 cm
94-1069L	Locking Handle	13 in 33 cm
94-1069L-45	Locking Handle	17 ¾ in 45 cm



Allis Grasper

Direct Drive, Monopolar Laparoscopy Forceps, 5.0 mm, Insulated

CODE	DESCRIPTION	LENGTH
94-1078L-26	Locking Handle	10 ¼ in 26 cm
94-1078L	Locking Handle	13 in 33 cm
94-1078L-45	Locking Handle	17 ¾ in 45 cm
94-1078SL-26	Slidelock Handle	10 ¼ in 26 cm
94-1078SL	Slidelock Handle	13 in 33 cm
94-1078SL-45	Slidelock Handle	17 ¾ in 45 cm



Allis Single Grasper

Direct Drive, Monopolar Laparoscopy Forceps, 5.0 mm, Insulated

CODE	DESCRIPTION	LENGTH
94-1088L-26	Locking Handle	10 ¼ in 26 cm
94-1088L	Locking Handle	13 in 33 cm
94-1088L-45	Locking Handle	17 ¾ in 45 cm

Babcock Grasper

Direct Drive, Monopolar Laparoscopy Forceps, 5.0 mm, Insulated

CODE	DESCRIPTION	LENGTH
94-1089L-26	Locking Handle	10 ¼ in 26 cm
94-1089L	Locking Handle	13 in 33 cm
94-1089L-45	Locking Handle	17 ¾ in 45 cm
94-1089SL-26	Slidelock Handle	10 ¼ in 26 cm
94-1089SL	Slidelock Handle	13 in 33 cm
94-1089SL-45	Slidelock Handle	17 ¾ in 45 cm



Single-use products are sold sterile.
Kits Include Forceps, 12 ft (3.6 m) Cord and Surgi-Wipe.

Ultra™ Monopolar Laparoscopic Forceps

DeBakey Grasper

Direct Drive, Monopolar Laparoscopy Forceps, 5.0 mm, Insulated

CODE	DESCRIPTION	LENGTH	
94-1091-26	Locking Handle	10 ¼ in	26 cm
94-1091	Locking Handle	13 in	33 cm
94-1091-45	Locking Handle	17 ¾ in	45 cm



DeBakey Grasper

Direct Drive, Monopolar Laparoscopy Forceps, 5.0 mm, Insulated

CODE	DESCRIPTION	LENGTH	
94-1091L-26	Locking Handle	10 ¼ in	26 cm
94-1091L	Locking Handle	13 in	33 cm
94-1091L-45	Locking Handle	17 ¾ in	45 cm
94-1091SL-26	Slidelock Handle	10 ¼ in	26 cm
94-1091SL	Slidelock Handle	13 in	33 cm
94-1091SL-45	Slidelock Handle	17 ¾ in	45 cm



DeBakey Long Grasper

Direct Drive, Monopolar Laparoscopy Forceps, 5.0 mm, Insulated

CODE	DESCRIPTION	LENGTH	
94-1092L-26	Locking Handle	10 ¼ in	26 cm
94-1092L	Locking Handle	13 in	33 cm
94-1092L-45	Locking Handle	17 ¾ in	45 cm
94-1092SL-26	Slidelock Handle	10 ¼ in	26 cm
94-1092SL	Slidelock Handle	13 in	33 cm
94-1092SL-45	Slidelock Handle	17 ¾ in	45 cm



Single-use products are sold sterile.
Kits Include Forceps, 12 ft (3.6 m) Cord and Surgi-Wipe.

Questions? Contact us at customerservice@symmetrysurgical.com
or learn more at symmetrysurgical.com

Ultra™ Monopolar Laparoscopic Forceps

"Hunter Type" Bowel Grasper

Direct Drive, Monopolar Laparoscopy Forceps, 5.0 mm, Insulated

CODE	DESCRIPTION	LENGTH
94-1093L-26	Locking Handle	10 ¼ in 26 cm
94-1093L	Locking Handle	13 in 33 cm
94-1093L-45	Locking Handle	17 ¾ in 45 cm
94-1093-26	Non-Locking Handle	10 ¼ in 26 cm
94-1093	Non-Locking Handle	13 in 33 cm
94-1093-45	Non-Locking Handle	17 ¾ in 45 cm



Metzenbaum Scissors

Direct Drive, Monopolar Laparoscopy Scissors, 5.0 mm, Insulated

CODE	DESCRIPTION	LENGTH
94-2006-26	Non-Locking Handle	10 ¼ in 26 cm
94-2006	Non-Locking Handle	13 in 33 cm
94-2006-45	Non-Locking Handle	17 ¾ in 45 cm



Straight Scissors

Direct Drive, Monopolar Laparoscopy Scissors, 5.0 mm, Insulated

CODE	DESCRIPTION	LENGTH
94-2009-26	Non-Locking Handle	10 ¼ in 26 cm
94-2009	Non-Locking Handle	13 in 33 cm
94-2009-45	Non-Locking Handle	17 ¾ in 45 cm

Hook Scissors

Direct Drive, Monopolar Laparoscopy Scissors, 5.0 mm, Insulated

CODE	DESCRIPTION	LENGTH
94-2020-26	Non-Locking Handle	10 ¼ in 26 cm
94-2020	Non-Locking Handle	13 in 33 cm
94-2020-45	Non-Locking Handle	17 ¾ in 45 cm

Single-use products are sold sterile.
Kits Include Forceps, 12 ft (3.6 m) Cord and Surgi-Wipe.

Ultra™ Monopolar Laparoscopic Forceps

Accu-SPIN™ Bottle Nose Dissecting Forcep

Monopolar Laparoscopy Forceps, 5.0 mm, Insulated

CODE	DESCRIPTION	LENGTH
97-1030RC	Non-Locking Handle	12 ½ in 32 cm



Accu-SPIN™ Cup Tooth Dissecting Forceps

Monopolar Laparoscopy Forceps, 5.0 mm, Insulated

CODE	DESCRIPTION	LENGTH
97-1035RC	Non-Locking Handle	12 ½ in 32 cm

Accu-SPIN™ Micro Fenestrated Forceps

Monopolar Laparoscopy Forceps, 5.0 mm, Insulated

CODE	DESCRIPTION	LENGTH
97-1042RC	Non-Locking Handle	12 ½ in 32 cm



Accu-SPIN™ Maryland Dissecting Forceps

Monopolar Laparoscopy Forceps, 5.0 mm, Insulated

CODE	DESCRIPTION	LENGTH
97-1067RC	Non-Locking Handle	12 ½ in 32 cm
97-1067RCL	Non-Locking Handle	17 ¾ in 45 cm

Single-use products are sold sterile.
Kits Include Forceps, 12 ft (3.6 m) Cord and Surgi-Wipe.

Questions? Contact us at customerservice@symmetrysurgical.com
or learn more at symmetrysurgical.com

Ultra™ Monopolar Laparoscopic Forceps

Accu-CUT™ Endo Scissors

Monopolar Laparoscopy Scissors, 5.0 mm, Insulated

CODE	DESCRIPTION	LENGTH
97-2006RC	Non-Locking Handle	12½ in 32 cm
97-2006RCL	Non-Locking Handle	17¾ in 45 cm



Accu-CUT™ Hook Scissors

Monopolar Laparoscopy Scissors, 5.0 mm, Insulated

CODE	DESCRIPTION	LENGTH
97-2020RC	Non-Locking Handle	12½ in 32 cm

Insulated Lap Handle

CODE	DESCRIPTION	LENGTH
90-0030	Take Apart	12½ in 32 cm



Bottle Nose Dissecting Forceps

Replacement Jaw Assembly, 5.0 mm, Insulated

CODE	DESCRIPTION	LENGTH
90-1030B	Take Apart	12½ in 32 cm



Ultra™ Monopolar Laparoscopic Forceps

Bottle Nose Dissecting Forceps

Monopolar Laparoscopy Forceps, with Insulated Handle, 5.0 mm

CODE	DESCRIPTION	LENGTH
90-1030I	Take Apart	12 ½ in 32 cm
90-1030IL	Take Apart	17 ¾ in 45 cm



Duckbill Grasping Forceps

Replacement Jaw Assembly, 5.0 mm, Insulated

CODE	DESCRIPTION	LENGTH
90-1040I	Take Apart	12 ½ in 32 cm
90-1040BL	Take Apart	17 ¾ in 45 cm



Fenestrated Dissecting Forceps

Replacement Jaw Assembly, 5.0 mm, Insulated

CODE	DESCRIPTION	LENGTH
90-1042BL	Take Apart	12 ½ in 32 cm
90-1042B	Take Apart	17 ¾ in 45 cm



Fenestrated Dissecting Forceps

Monopolar Laparoscopy Forceps, with Insulated Handle, 5.0 mm

CODE	DESCRIPTION	LENGTH
90-1042I	Take Apart	12 ½ in 32 cm
90-1042IL	Take Apart	17 ¾ in 45 cm



Single-use products are sold sterile.
Kits Include Forceps, 12 ft (3.6 m) Cord and Surgi-Wipe.

Questions? Contact us at customerservice@symmetrysurgical.com
or learn more at symmetrysurgical.com

Ultra™ Monopolar Laparoscopic Forceps

Maryland Dissecting Forceps

Replacement Jaw Assembly, 5.0 mm, Insulated

CODE	DESCRIPTION	LENGTH
90-1065B	Take Apart	12 ½ in 32 cm
90-1065I	Take Apart	12 ½ in 32 cm



Mixer Right Angle Dissector

Replacement Jaw Assembly, 5.0 mm, Insulated

CODE	DESCRIPTION	LENGTH
90-1066B	Take Apart	12 ½ in 32 cm
90-1066I	Take Apart	12 ½ in 32 cm



Maryland Dissector

Replacement Jaw Assembly, 5.0 mm, Insulated

CODE	DESCRIPTION	LENGTH
90-1067B	Take Apart	12 ½ in 32 cm
90-1067BL	Take Apart	12 ½ in 32 cm



Maryland Dissecting Forceps

Monopolar Laparoscopy Forceps, with Insulated Handle, 5.0 mm

CODE	DESCRIPTION	LENGTH
90-1067I	Take Apart	12 ½ in 32 cm
90-1067IL	Take Apart	17 ¾ in 45 cm



Ultra™ Monopolar Laparoscopic Forceps

Mixer Right Angle Dissector

Monopolar Laparoscopy Forceps, with Insulated Handle, 10 mm


CODE	DESCRIPTION	LENGTH
90-1068I	Take Apart	12½ in 32 cm



Insulated Outer Sheath

Take Apart, Insulated

CODE	DESCRIPTION	LENGTH
90-1115	5.0 mm, With Flush Port	12½ in 32 cm
90-1115L	5.0 mm, With Flush Port	17¾ in 45 cm
90-1116	5.0 mm, With Flush Port	12½ in 32 cm
90-1116L	5.0 mm, With Flush Port	17¾ in 45 cm
90-1117	10 mm, Without Flush Port	12½ in 32 cm



Bottle Nose Dissecting Forceps

Monopolar Laparoscopy Forceps, 5.0 mm, Insulated

CODE	DESCRIPTION	LENGTH
90-1030	One-Piece	12½ in 32 cm



Bullet Nose Dissecting Forceps

Monopolar Laparoscopy Forceps, 5.0 mm, Insulated

CODE	DESCRIPTION	LENGTH
90-1033	One-Piece	12½ in 32 cm



Single-use products are sold sterile.
Kits Include Forceps, 12 ft (3.6 m) Cord and Surgi-Wipe.

Questions? Contact us at customerservice@symmetrysurgical.com
or learn more at symmetrysurgical.com

Ultra™ Monopolar Laparoscopic Forceps

Cup Tooth Dissecting Forceps

Monopolar Laparoscopy Forceps, 5.0 mm, Insulated

CODE	DESCRIPTION	LENGTH
90-1035	One-Piece	12½ in 32 cm



Duckbill Grasping Forceps

Monopolar Laparoscopy Forceps, 5.0 mm, Insulated

CODE	DESCRIPTION	LENGTH
90-1040	One-Piece	12½ in 32 cm



Micro Fenestrated Dissecting Forceps

Monopolar Laparoscopy Forceps, 5.0 mm, Insulated

CODE	DESCRIPTION	LENGTH
90-1044	One-Piece	12½ in 32 cm



Maryland Dissecting Forceps

One-Piece, Monopolar Laparoscopy Forceps, 5.0 mm, Insulated

CODE	DESCRIPTION	LENGTH
90-1067	5.0 mm	12½ in 32 cm
90-1065	10 mm	12½ in 32 cm

Ultra™ Monopolar Laparoscopic Forceps

Hook Biopsy Forceps

Monopolar Laparoscopy Forceps, 5.0 mm, Insulated

CODE	DESCRIPTION	LENGTH
90-1073	One-Piece	12½ in 32 cm



Biopsy Forceps

Monopolar Laparoscopy Forceps, 5.0 mm, Insulated

CODE	DESCRIPTION	LENGTH
90-1075	One-Piece	12½ in 32 cm



Hook Scissors

Monopolar Laparoscopy Forceps, 5.0 mm, Insulated

CODE	DESCRIPTION	LENGTH
90-2020	One-Piece	12½ in 32 cm



Single-use products are sold sterile.
Kits Include Forceps, 12 ft (3.6 m) Cord and Surgi-Wipe.

Questions? Contact us at customerservice@symmetrysurgical.com
or learn more at symmetrysurgical.com

Olsen® Monopolar Blade Electrodes

Blade Electrodes

Single-Use, Straight


CODE	DESCRIPTION	LENGTH	QTY
97000	Straight Electrode	2 ³ / ₄ in 70 mm	50
97001	Straight Electrode	6 in 152 mm	25
97002	Straight Electrode	8 in 203 mm	25



Blade Electrodes, with Midas Touch® Non-Stick, Straight

Single-Use, Standard Exposure, 10 per Box

CODE	DESCRIPTION	LENGTH
30-0001	Straight Blade Electrode	2 ³ / ₄ in 70 mm
30-0002	Straight Blade Electrode	4 in 102 mm
30-0003	Straight Blade Electrode	5 in 127 mm
30-0004	Straight Blade Electrode	6 in 152 mm



Blade Electrodes, with Midas Touch® Non-Stick, Straight

Single-Use, Modified 2.0 mm Exposure, 10 per Box

CODE	DESCRIPTION	LENGTH
30-1011	Straight Blade Electrode	2 ³ / ₄ in 70 mm
30-1012	Straight Blade Electrode	4 in 102 mm
30-1013	Straight Blade Electrode	5 in 127 mm
30-1014	Straight Blade Electrode	6 in 152 mm

Single-use products are sold sterile.
Kits Include Forceps, 12 ft (3.6 m) Cord and Surgi-Wipe.

Olsen® Electrode Extenders

Electrode Extenders

Reusable, 1 Each

CODE	DESCRIPTION	LENGTH
67510	Electrode Extender	3 in 76 mm
67511	Electrode Extender	4 in 102 mm
67512	Electrode Extender	5 in 127 mm
67513	Electrode Extender	6 in 152 mm
67516	Electrode Extender	7 in 178 mm
67514	Electrode Extender	8 in 203 mm
67515	Electrode Extender	10 in 254 mm



ENERGY

Ultra™ Laparoscopic Electrodes

Laparoscopic Monopolar Electrodes

Single-Use, 10 per Box

CODE	DESCRIPTION	LENGTH
90-9525-S	J-Hook	13in 33 cm
90-9525-45S	J-Hook	17 ³ / ₄ 45 cm
90-9522-S	L-Hook	13in 33 cm
90-9522-45S	L-Hook	17 ³ / ₄ 45 cm
90-9523-S	Spatula	13in 33 cm
90-9523-45S	Spatula	17 ³ / ₄ 45 cm



Single-use products are sold sterile.
Kits Include Forceps, 12 ft (3.6 m) Cord and Surgi-Wipe.

Questions? Contact us at customerservice@symmetrysurgical.com
or learn more at symmetrysurgical.com

Olsen® Monopolar Loop Electrodes

LLETZ Loop Electrodes, 2.2 in (5.6 cm) Length

Reusable, 1 Each

CODE	LOOP SIZE
69504	15 x 10 mm
69503	20 x 10 mm
69507	20 x 15 mm

LLETZ Loop Electrodes, 5 in (12.7 cm) Length

Reusable, 1 Each

CODE	LOOP SIZE
69100	15 x 8.0 mm
69300	20 x 8.0 mm
69000	10 x 10 mm
69505	10 x 12 mm
69500	15 x 12 mm
69400	20 x 12 mm
69200	15 x 15 mm
69506	20 x 15 mm
69502	20 x 20 mm



LLETZ Square Loop Electrodes, 5 in (12.7 cm) Length

Reusable, 1 Each

CODE	LOOP SIZE
69501	10 x 7.0 mm
69001	10 x 10 mm



Loop Electrodes

Reusable, 1 Each

CODE	LOOP SIZE		LENGTH	
66300	1/8 in	3.2 mm	1 5/8 in	41 mm
66301	1/4 in	6.4 mm	1 3/4 in	44 mm
66302	5/16 in	7.9 mm	1 13/16 in	59 mm
66306	5/16 in	7.9 mm	4 5/16 in	110 mm
66303	3/8 in	9.5 mm	1 7/8 in	48 mm
66307	3/8 in	9.5 mm	4 3/8 in	111 mm
66304	1/2 in	12.7 mm	2.0 in	51 mm
66308	1/2 in	12.7 mm	4 1/2 in	114 mm
66305	3/4 in	19.1 mm	2 1/4 in	57 mm



Single-use products are sold sterile.
Kits Include Forceps, 12 ft (3.6 m) Cord and Surgi-Wipe.

Questions? Contact us at customerservice@symmetrysurgical.com
or learn more at symmetrysurgical.com

Olsen® Monopolar Needle Electrodes

Needle Electrodes, Blunt

Reusable, Straight, Standard 10 mm Exposure, 1 Each

CODE	DESCRIPTION	LENGTH
66114	Blunt Needle-Tip Electrode	2 ¾ in 70 mm
66115	Blunt Needle-Tip Electrode	4 in 102 mm
66116	Blunt Needle-Tip Electrode	5 in 127 mm
66117	Blunt Needle-Tip Electrode	6 in 152 mm
66118	Blunt Needle-Tip Electrode	7 in 178 mm
66119	Blunt Needle-Tip Electrode	8 in 203 mm



Needle Electrodes, Standard

Reusable, Straight, Standard 10 mm Exposure, 1 Each


CODE	DESCRIPTION	LENGTH
66101	Standard Needle-Tip Electrode	2 ¾ in 70 mm
66102	Standard Needle-Tip Electrode	4 in 102 mm
66103	Standard Needle-Tip Electrode	5 in 127 mm
66104	Standard Needle-Tip Electrode	6 in 152 mm
66105	Standard Needle-Tip Electrode	7 in 178 mm
66106	Standard Needle-Tip Electrode	8 in 203 mm



Needle Electrodes, Sharp

Reusable, Straight, 1 Each

CODE	DESCRIPTION	LENGTH
66107	Sharp Needle-Tip Electrode	2 ¾ in 70 mm
66108	Sharp Needle-Tip Electrode	4 in 10.2 cm
66109	Sharp Needle-Tip Electrode	5 in 127 mm
66110	Sharp Needle-Tip Electrode	6 in 152 mm
66111	Sharp Needle-Tip Electrode	7 in 178 mm
66112	Sharp Needle-Tip Electrode	8 in 203 mm



Single-use products are sold sterile.
Kits Include Forceps, 12 ft (3.6 m) Cord and Surgi-Wipe.

Questions? Contact us at customerservice@symmetrysurgical.com
or learn more at symmetrysurgical.com

Olsen® Monopolar Needle Electrodes

Needle Electrode, Short/Close-Up

Reusable, Straight, Standard 10 mm Exposure, 1 Each

CODE	DESCRIPTION	LENGTH
66800	Needle-Tip Electrode	2 in 51 mm



Needle Electrodes, Blunt

Reusable, Straight, Modified 2.0 mm Exposure, 1 Each

CODE	DESCRIPTION	LENGTH
67003	Blunt Needle-Tip Electrode	2 3/4 in 70 mm
67004	Blunt Needle-Tip Electrode	6 in 152 mm



Needle Electrodes, Standard

Reusable, Straight, Modified 2.0 mm Exposure, 1 Each

CODE	DESCRIPTION	LENGTH
67000	Standard Needle-Tip Electrode	2 3/4 in 70 mm
67001	Standard Needle-Tip Electrode	4 in 102 mm
67002	Standard Needle-Tip Electrode	6 in 152 mm
67008	Standard Needle-Tip Electrode	8 3/4 in 222 mm



Needle Electrodes, Standard

Reusable, Straight, Modified 2.0 mm Exposure, 1 Each

CODE	DESCRIPTION	LENGTH
67008	Standard Needle-Tip Electrode, Bayonet	8 3/4 in 222 mm

Single-use products are sold sterile.
Kits Include Forceps, 12 ft (3.6 m) Cord and Surgi-Wipe.

Olsen® Monopolar Needle Electrodes

Needle Electrodes

Single-Use, Straight

CODE	DESCRIPTION	LENGTH	QTY
97100	Needle Electrode	2 ³ / ₄ in 70 mm	50
97101	Needle Electrode	6 in 152 mm	10



Needle Electrodes, Straight

Single-Use, Tungsten Wire Micro Dissection, 1 Each

CODE	DESCRIPTION	LENGTH
97110	Standard Sharp Needle Electrode	1 ¹ / ₂ in 38 mm
97111	Standard Sharp Needle Electrode	2 in 51 mm
97112	Standard Sharp Needle Electrode	2 ³ / ₈ in 60 mm
97116	Standard Sharp Needle Electrode	3 ¹ / ₂ in 89 mm
97117	Standard Sharp Needle Electrode	4 ¹ / ₂ in 114 mm
97114	Standard Sharp Needle Electrode	6 in 152 mm



Needle Electrodes, 45° Angled Shaft

Single-Use, Tungsten Wire Micro Dissection, 10 mm from Tips, 1 Each


CODE	DESCRIPTION	LENGTH
97115	Standard Sharp Needle Electrode	2 ³ / ₈ in 60 mm



Needle Electrode, Bayonet Shaft

Single-Use, Tungsten Wire Micro Dissection, 1 Each

CODE	DESCRIPTION	LENGTH
97125	Ultra Sharp Needle Electrode	6 in 152 mm



Single-use products are sold sterile.
Kits Include Forceps, 12 ft (3.6 m) Cord and Surgi-Wipe.


Questions? Contact us at customerservice@symmetrysurgical.com
or learn more at symmetrysurgical.com

Olsen® Monopolar Needle Electrodes

Needle Electrodes, Straight

Single-Use, Ultra Sharp, Tungsten Wire Micro Dissection, 1 Each


CODE	DESCRIPTION	LENGTH
97130	Ultra Sharp Needle Electrode	1¼ in 32 mm
97131	Ultra Sharp Needle Electrode	1⅝ in 41 mm
97132	Ultra Sharp Needle Electrode	2 in 51 mm
97133	Ultra Sharp Needle Electrode	2½ in 64 mm
97134	Ultra Sharp Needle Electrode	4⅝ in 118 mm
97135	Ultra Sharp Needle Electrode	5½ in 140 mm
97136	Ultra Sharp Needle Electrode	6⅝ in 168 mm



Needle Electrodes, Angled Tip, 3.0 mm

Single-Use, Ultra Sharp, Tungsten Wire Micro Dissection, 1 Each


CODE	DESCRIPTION	LENGTH
97141	Ultra Sharp Needle Electrode	1⅝ in 41 mm
97143	Ultra Sharp Needle Electrode	2 in 51 mm
97144	Ultra Sharp Needle Electrode	4 in 102 mm
97145	Ultra Sharp Needle Electrode	4⅝ in 118 mm
97147	Ultra Sharp Needle Electrode	6⅝ in 168 mm
97148	Ultra Sharp Needle Electrode	7 in 178 mm



Needle Electrodes, Angled Shaft, 10 mm

Single-Use, Ultra Sharp, Tungsten Wire Micro Dissection, 1 Each

CODE	DESCRIPTION	LENGTH
97150	Ultra Sharp Needle Electrode	2 in 51 mm
97152	Ultra Sharp Needle Electrode	4⅞ in 124 mm
97154	Ultra Sharp Needle Electrode	6⅞ in 175 mm



Needle Electrodes, Curved Shaft

Single-Use, Ultra Sharp, Tungsten Wire Micro Dissection, 1 Each

CODE	DESCRIPTION	LENGTH
97170	Ultra Sharp Needle Electrode	4 in 102 mm
97171	Ultra Sharp Needle Electrode	5⅞ in 149 mm




Single-use products are sold sterile.
Kits Include Forceps, 12 ft (3.6 m) Cord and Surgi-Wipe.

Olsen® Monopolar Needle Electrodes

Needle Electrodes, with Midas Touch® Non-Stick

Single-Use, Standard Exposure, Straight, 10 per Box


CODE	DESCRIPTION	LENGTH
30-1001	Needle Electrode	2¾ in 70 mm
30-1002	Needle Electrode	4 in 102 mm
30-1003	Needle Electrode	5 in 127 mm
30-1004	Needle Electrode	6 in 152 mm



Needle Electrode, with Midas Touch® Non-Stick

Single-Use, Modified 2.0 mm Exposure, Bayonet, 10 per Box

CODE	DESCRIPTION	LENGTH
30-1015	Needle Electrode	5¾ in 146 mm



Olsen® Suction Coagulators

Suction Coagulators

Single-Use, 6 in (152 mm), 10 French, 10 per Box

CODE	DESCRIPTION
20-8020	Foot-Controlled, Requires 960xx Series Adapter
20-8040	Hand-Controlled



Single-use products are sold sterile.
Kits Include Forceps, 12 ft (3.6 m) Cord and Surgi-Wipe.

Questions? Contact us at customerservice@symmetrysurgical.com
or learn more at symmetrysurgical.com

Olsen® Electrosurgical Cords

Bipolar Cords

Reusable, 1 Each

CODE	DESCRIPTION	PLUG	LENGTH
87006	Shrouded Instrument Connector	Fixed-Pin Generator	12 ft 3.6 m
87110	European Flat Instrument Connector	Fixed-Pin Generator	12 ft 3.6 m
50-1106	45° Angled Instrument Connector	Fixed-Pin Generator	12 ft 3.6 m
87600	Instrument Connector	Wolf-Style Generator	12 ft 3.6 m



Bipolar Cords

Single-Use, 1 Each

CODE	DESCRIPTION	PLUG	LENGTH
40-1105	Straight Instrument Connector	Fixed-Pin Generator	12 ft 3.6 m
40-1106	45° Angled Instrument Connector	Fixed-Pin Generator	12 ft 3.6 m
40-1107	Straight Instrument Connector	Fixed-Pin Generator	12 ft 3.6 m
40-1114	Instrument Connector	Wolf-Style Generator	12 ft 3.6 m
40-1115	Wolf-Style Instrument Connector	Valleylab-Style Generator	12 ft 3.6 m



Integrated Irrigation Tubing and Bipolar Cord Set for Rotary Pumps

Integrated Irrigation Tubing and Bipolar Cord Set for Rotary Pumps, 10 per Box

CODE	DESCRIPTION	PLUG	LENGTH
40-3003R	On/Off Clamp	Fixed-Pin Generator	12 ft 3.6 m



Monopolar Cords

Reusable, 1 Each

CODE	DESCRIPTION	PLUG	LENGTH
20002	2.0 mm European-Style Recessed Pin Female	5/32 in (4.0 mm) Male Generator	10 ft 3.0 m
20003	TUR/Resectoscope Female	Valleylab-Style Male Generator	10 ft 3.0 m
20006	2.0 mm European-Style Recessed Pin Female	Valleylab-Style Male Generator	10 ft 3.0 m

Single-use products are sold sterile.
Kits Include Forceps, 12 ft (3.6 m) Cord and Surgi-Wipe.

Olsen® Electrosurgical Cords



Monopolar Cords

Reusable, 1 Each

CODE	DESCRIPTION	PLUG	LENGTH
90018	Olsen® Threaded Instrument Connector	$\frac{5}{32}$ in (4.0 mm) Male Generator	10 ft 3.0 m
90020	Olsen® Threaded Instrument Connector	Valleylab-Style Male Generator	10 ft 3.0 m
90054	Laparoscopic Instrument Connector	$\frac{5}{32}$ in (4.0 mm) Male Generator	10 ft 3.0 m
90056	Laparoscopic Instrument Connector	Valleylab-Style Male Generator	10 ft 3.0 m
90065	3.0 mm Olympus-Style Snare Connector	$\frac{5}{32}$ in (4.0 mm) Male Generator	10 ft 3.0 m
90067	3.0 mm Olympus-Style Snare Connector	Valleylab-Style Male Generator	10 ft 3.0 m
50-2107	Endoscopic Instrument Connector, Straight	Valleylab-Style Male Generator	10 ft 3.0 m
50-2111	Endoscopic Instrument Connector, Angled 90°	Valleylab-Style Male Generator	10 ft 3.0 m



Monopolar Cords

Single-Use, 10 per Box

CODE	DESCRIPTION	PLUG	LENGTH
40-2110	Endoscopic Instrument Connector, Straight	Valleylab-Style Generator	10 ft 3.0 m
40-2106	Endoscopic Instrument Connector, Angled 90°	Valleylab-Style Generator	10 ft 3.0 m



Monopolar Cords

Single-Use, 50 per Box

CODE	DESCRIPTION	PLUG	LENGTH
92002	Resectoscope Instrument Connector, Angled 90°	Crimped Wire Generator	10 ft 3.0 m

Single-use products are sold sterile.
Kits Include Forceps, 12 ft (3.6 m) Cord and Surgi-Wipe.

Questions? Contact us at customerservice@symmetrysurgical.com
or learn more at symmetrysurgical.com

Olsen® Exposed Top Forceps

Adson Forceps

Reusable, Serrated Tips, Straight, 1 Each



CODE	DESCRIPTION	TIP SIZE	LENGTH
42509	Insulated, Exposed Top	1.0 mm	5 in 127 mm

Bayonet Forceps

Reusable, Insulated, Smooth Tips, 1 Each



CODE	DESCRIPTION	TIP SIZE	LENGTH
41207	Insulated, Exposed Top	1.0 mm	7¾ in 197 mm
41206	Insulated, Exposed Top	1.5 mm	7¾ in 197 mm
41406	Insulated, Exposed Top	1.0 mm	8 in 203 mm
41405	Insulated, Exposed Top	2.0 mm	8 in 203 mm
41607	Insulated, Exposed Top	2.0 mm	8½ in 216 mm

DeBakey Forceps

Reusable, 1 Each



CODE	DESCRIPTION	TIP SIZE	LENGTH
44012	Insulated, Exposed Top	1.0 mm	6 in 127 mm
44006	Insulated, Exposed Top	1.5 mm	7¾ in 197 mm
44004	Insulated, Exposed Top	1.5 mm	9½ in 241 mm

Potts-Smith Forceps


Reusable, Exposed Top, Serrated Tips, 1 Each



Length (L)	7½ in 191 mm	8½ in 216 mm	9¾ in 248 mm	10 in 254 mm	12 in 305 mm
TIP SIZE	INSULATED	INSULATED	INSULATED	INSULATED	INSULATED
0.5 mm	42211		42411		
1.0 mm	42209	42309	42410		44005
1.5 mm	42208	42308	42409	44013	44000
2.0 mm	42207	42307	42408		

Single-use products are sold sterile.
Kits Include Forceps, 12 ft (3.6 m) Cord and Surgi-Wipe.

Olsen® Monopolar Adapters



Monopolar Generator Adapter

Electrosurgical Accessory, 1 Each

CODE	DESCRIPTION	PLUG
96000	Crimped Wire $\frac{1}{8}$ in (3.2 mm) Cord Female	$\frac{5}{32}$ in (4.0 mm) Male Generator Plug
96001	Crimped Wire $\frac{1}{8}$ in (3.2 mm) Cord Female	Valleylab-Style Male Generator Plug
96002	European/Wolf-Style Female	Valleylab-Style Male Generator Plug

History of Symmetry Surgical's Energy Portfolio

In 1973, Olsen® Medical introduced the first line of single-use monopolar electrosurgical handles for use as electrosurgical pens and pencils.

Now part of Symmetry Surgical, Olsen®'s innovative energy portfolio has expanded to include our Extendable Electrosurgical Pencil as well as a wide array of electrodes and monopolar forceps, cords and adapters to satisfy your needs.



Symmetry Surgical's Olsen® Bipolar Forceps

In 1988, Olsen® Medical introduced the first line of cost-effective, single-use bipolar forceps, and ten years later we introduced our innovative Midas Touch® non-stick bipolar forceps.

Today, as part of Symmetry Surgical, Olsen® continues to focus on innovative and cost-effective ways to improve patient outcomes.



Single-use products are sold sterile.
Kits Include Forceps, 12 ft (3.6 m) Cord and Surgi-Wipe.

Questions? Contact us at customerservice@symmetrysurgical.com or learn more at symmetrysurgical.com

Olsen® Monopolar Pencils

Electrosurgical Pencils

Reusable, Blade Electrode, Foot-Activated, 1 Each



CODE	DESCRIPTION	PLUG
65000	Attached 10 ft (3.0m) Cord	Valleylab-Style Generator
65006	Attached 10 ft (3.0m) Cord	European/Wolf-Style Generator

Electrosurgical Pencil

Reusable, Blade Electrode, Hand-Activated, Button Switch, 1 Each



CODE	DESCRIPTION	PLUG
75000	Attached 10 ft (3.0 m) Cord	3-Prong Generator

Electrosurgical Pencils

Single-Use, Blade Electrode, Hand-Activated, Button Switch, 50 per Box



CODE	DESCRIPTION	PLUG
75530	Attached 10 ft (3.0 m) Cord	3-Prong Generator
75540	Attached 10 ft (3.0 m) Cord, Holster, Tip Cleaner	3-Prong Generator
75541	Attached 10 ft (3.0 m) Cord, Holster	3-Prong Generator

Olsen® Monopolar Pencils



Electrosurgical Pencils

Single-Use, Blade Electrode, Hand-Activated, Rocker Switch, 50 per Box

CODE	DESCRIPTION	PLUG
75610	Attached 10 ft (3.0 m) Cord	3-Prong Generator
75620	Attached 10 ft (3.0 m) Cord, Holster, Tip Cleaner	3-Prong Generator
75621	Attached 10 ft (3.0 m) Cord, Holster	3-Prong Generator

Electrosurgical Pencils

Single-Use, Blade Electrode, Foot-Activated



CODE	DESCRIPTION	PLUG	QTY
95012	Attached 10 ft (3.0 m) Cord	Crimped Wire Connector (uses 960 Series Adapter)	100/Box
95017	Attached 10 ft (3.0 m) Cord, Packaged with Holster and Tip Cleaner	Crimped Wire Connector (uses 960 Series Adapter)	50/Box



Extendable Electrosurgical Pencils, with Midas Touch® Non-Stick

Single-Use, Midas Touch Non-Stick Blade Electrode, Hand-Activated, Button Switch, 25 per Box

CODE	DESCRIPTION	ADJUSTMENT
75800	Attached 10 ft (3.0 m) Cord with 3-Prong Generator Plug, Standard Blade	Easily Adjusts from 2 ½ to 5 ½ in
75802	Attached 10 ft (3.0 m) Cord with 3-Prong Generator Plug, Holster, Standard Blade	Easily Adjusts from 2 ½ to 5 ½ in
75810	Attached 10 ft (3.0 m) Cord with 3-Prong Generator Plug, Modified Blade	Easily Adjusts from 2 ½ to 5 ½ in
75820	Attached 10 ft (3.0 m) Cord with 3-Prong Generator Plug, Modified Needle	Easily Adjusts from 2 ½ to 5 ½ in

Single-use products are sold sterile.

Kits Include Forceps, 12 ft (3.6 m) Cord and Surgi-Wipe.

Questions? Contact us at customerservice@symmetrysurgical.com or learn more at symmetrysurgical.com

Olsen® Energy Accessories

Anesthesia Circuit Shield

Laser Safety Product, 6 per Box

CODE	DESCRIPTION	SIZE
60-7100	Laser Safety Product	12 x 48 in 305 mm x 1.2 m

Holsters

Single-Use, Sterile, 50 per Box

CODE	DESCRIPTION
77000	Holster Packaged with Tip Cleaner
77001	Holster Only
77005	Holster with Self-Retaining Metal Clip
77006	Holster with Self-Retaining Metal Clip, Packaged with Tip Cleaner



77005

Graves Vaginal Speculums, with Smoke Tube

CODE	DESCRIPTION	SIZE
50001	Insulated, with Smoke Tube	3 x ¾ in 76 x 20 mm
50101	Insulated, with Smoke Tube	4 x 1 ⅜ in 102 x 35 mm
50201	Insulated, with Smoke Tube	4 ½ x 1 ½ in 114 x 35 mm
50002	Non-Insulated, Non-Reflective	3 x ¾ in 76 x 20 mm
50102	Non-Insulated, Non-Reflective	4 x 1 ⅜ in 102 x 35 mm
50202	Non-Insulated, Non-Reflective	4 ½ x 1 ½ in 114 x 35 mm



Open-Sided Vaginal Speculum

CODE	DESCRIPTION	SIZE
50301	Insulated, with Smoke Tube	4 ½ x 1 ⅜ in 114 x 35 mm



Single-use products are sold sterile.
Kits Include Forceps, 12 ft (3.6 m) Cord and Surgi-Wipe.

Olsen® Energy Accessories

Binocular Eye Shield

Single-Use, Laser Safety Product

CODE	DESCRIPTION	QTY
60-7000-OL	Eye Shield	12/Box



Personnel Shield

Single-Use, Drapes Between Sterile Field & Anesthesiologist, 6 per Box

CODE	DESCRIPTION	SIZE
60-7200	Sterile Safety Drape	30 x 60 in 762 mm x 1.5 m

ProClean™ Tip Cleaning Pads with Adhesive Back


Single-Use, Sterile

CODE	DESCRIPTION	QTY
97500	Electrosurgical Accessory	50/Box
97501	Electrosurgical Accessory	144/Box



Replacement Stylets for Irrigating Bipolar Forceps

CODE	DESCRIPTION	LENGTH	QTY
08-1130	For Irrigating Bipolar Forceps		12/Box
60-0300	For Irrigating Bipolar Forceps	10 ¼ in 260 mm	5/Box



Single-use products are sold sterile.
Kits Include Forceps, 12 ft (3.6 m) Cord and Surgi-Wipe.

Questions? Contact us at customerservice@symmetrysurgical.com
or learn more at symmetrysurgical.com