




# Reusable Tips for the Symmetry Sharp Kerrison® Handle



- Convenient – No more sharpening! When the tips dull, dispose of them and insert new tips
- Innovative – Patent-protected Symmetry Surgical® technology
- Versatile – Three tip sizes and three handle sizes allow for nine possible configurations
- Reliable – Manufactured in Germany from high-quality stainless steel

# Introducing the Symmetry Sharp Kerrison®



Visit [symmetrysurgical.com](https://symmetrysurgical.com) or contact your local Sales Representative to learn more about the new, innovative Symmetry Sharp Kerrison®.

CODE	DESCRIPTION	QTY
53-1662R	Reusable Tip for Symmetry Sharp Kerrison® Handle, 2mm	Pack of 2 tips
53-1663R	Reusable Tip for Symmetry Sharp Kerrison® Handle, 3mm	Pack of 2 tips
53-1664R	Reusable Tip for Symmetry Sharp Kerrison® Handle, 4mm	Pack of 2 tips
53-1674	Symmetry Sharp Kerrison® 5" Handle for 7" Kerrison	Single Unit
53-1675	Symmetry Sharp Kerrison® 6" Handle for 8" Kerrison	Single Unit
53-1676	Symmetry Sharp Kerrison® 7" Handle for 9" Kerrison	Single Unit
53-1674RC	RapidClean® Symmetry Sharp Kerrison Handle 5" for 7" Kerrison	Single Unit
53-1675RC	RapidClean® Symmetry Sharp Kerrison Handle 6" for 8" Kerrison	Single Unit
53-1676RC	RapidClean® Symmetry Sharp Kerrison Handle 7" for 9" Kerrison	Single Unit

## Optional tray for holding tips

CODE	DESCRIPTION
1-4800	Riley™ Base, Lid, & Mat for 4 X 8 Minitainer

