


Introducing Reddick® Cholangiogram Catheter

- Scoop tip curved introducer and balloon design for stable catheter placement – without clamping
- Depth marks to visualize catheter location
- Dual lumen, non-radiopaque catheter
- Needle introducer option

Simplifies dye-enhanced imaging for more complete stone detection and anatomy identification for gallbladder removal



Reddick® Cholangiogram Catheter

Reddick® Scoop Tip Cholangiogram Catheter



Introducer Options:

Scoop Tip Introducer and Needle Introducer



Complementary Laparoscopic Instrument

Reddick-Saye Screw Retractor Kit

Ideal when grasping instruments fail. Screws into distended gall bladder or cystic ovary wall for positive fixation and excellent retraction.



Visit [symmetrysurgical.com](https://www.symmetrysurgical.com) or contact your local Sales Representative for more information.

Reddick® Scoop Tip Cholangiogram Catheter

| CODE | DESCRIPTION | DIAMETER | LENGTH |
|----------|---|----------|--------|
| e2401-50 | Reddick® Cholangiogram Catheter | 4F | 50 cm |
| e2401-51 | Reddick® Cholangiogram Stiffer Catheter | 4F | 50 cm |
| e2401-52 | Reddick® Cholangiogram Catheter with Needle Introducer (<i>Needle Introducer: 13 gauge by 8 inches</i>) | 4F | 50 cm |
| e2401-53 | Reddick® Cholangiogram Stiffer Catheter with Needle Introducer (<i>Needle Introducer: 13 gauge by 8 inches</i>) | 4F | 50 cm |

Reddick® Cholangiogram Catheter

| CODE | DESCRIPTION | DIAMETER | LENGTH |
|----------|---|----------|--------|
| e2400-50 | Reddick® Cholangiogram Catheter | 4F | 50 cm |
| e2400-51 | Reddick® Cholangiogram Stiffer Catheter | 4F | 50 cm |
| e2400-52 | Reddick® Cholangiogram Catheter with Needle Introducer (<i>Needle Introducer: 13 gauge by 8 inches</i>) | 4F | 50 cm |
| e2400-53 | Reddick® Cholangiogram Stiffer Catheter with Needle Introducer (<i>Needle Introducer: 13 gauge by 8 inches</i>) | 4F | 50 cm |

Laparoscopic Instruments (4 per box)

| CODE | DESCRIPTION | DIAMETER | LENGTH |
|----------|----------------------------------|----------|--------|
| e2504-01 | Reddick-Saye Screw Retractor Kit | 5 mm | 30 cm |

1 Reddick®-Saye Screw Retractor Kit and 1 13G Aspiration Needle



3034 Owen Drive | Antioch, TN 37013 | 615.883.9090 | 800.251.3000 | customerservice@symmetrysurgical.com | [symmetrysurgical.com](https://www.symmetrysurgical.com)

©2018 Symmetry Surgical, Inc. All rights reserved. SM 0076