


## Introducing PretzelFlex™

### *Reusable Surgical Retraction System*

The PretzelFlex™ retraction system is a unique reusable pretzel-shaped organ and tissue retractor for use in minimally invasive surgery. It provides improved access and visualization of surgical sites. Available in 3mm and 5mm sizes.

- Precision-machined to maintain rigidity and shape
- Smooth slim-line shaft designed to ensure tissue is not penetrated or damaged
- Innovative segment technology ensures minimal space is required inside the abdomen to form into pretzel shape




# PretzelFlex™


# PretzelFlex™




Fully reusable system



Impressive strength and stability to support large organs




Creates shape in smaller space envelope compared to triangular retraction



Large surface area ensures increased contact

## Available in 3 mm and 5 mm sizes



Visit [symmetrysurgical.com](https://www.symmetrysurgical.com) or contact your local Sales Representative for more information.

### PretzelFlex™ Reusable Surgical Retraction System

CODE	DESCRIPTION	PLATFORM WIDTH	WORKING LENGTH
SHR7001	PretzelFlex™ 5mm	100 mm	350 mm
SHR7101	Replacement sterilization sleeve	(n/a)	(n/a)
SHR7301	PretzelFlex™ 3 mm	70 mm	275 mm
SHR7331	Replacement sterilization sleeve	(n/a)	(n/a)

PretzelFlex™ is proudly distributed by Symmetry Surgical® and is a product of Surgical Innovations Ltd. All trademarks and registered trademarks are the property of their respective owners.

Distributed by  
**symmetry surgical**®

3034 Owen Drive | Antioch, TN 37013 | 615.883.9090 | 800.251.3000 | [customerservice@symmetrysurgical.com](mailto:customerservice@symmetrysurgical.com) | [symmetrysurgical.com](https://www.symmetrysurgical.com)