


## FastClamp™

### Laparoscopic Clamping System

FastClamp™ is a high specification laparoscopic clamping solution. It ensures quick, stable, unobtrusive instrument and scope fixation.

#### FastClamp™ clamping arm:

- Has multiple degrees of freedom yet can be securely locked to allow instruments to be positioned accurately and quickly
- Securely grips either 5mm or 10mm instruments by way of the 'soft jaws' which also minimize the risk of instrument damage
- Uniquely designed to allow one or two arms to be used with each table clamp
- Two arm lengths available: standard and bariatric, ideally complementing PretzelFlex™ and other retraction instruments



#### FastClamp™ table clamp:

- Easily positioned with a wide range of options of adjustment for secure yet rapid clamping to any standard table rail
- The clamp handle's toggle action ensures a rigid connection to the operating table
- Whole system can be released in seconds

Visit [symmetrysurgical.com](http://symmetrysurgical.com) or contact your local Sales Representative for more information.

#### FastClamp™ Laparoscopic Clamping System

##### FastClamp™ Clamping Arm Options

CODE	DESCRIPTION
A-1000	FastClamp™ System Includes: 1 Table Clamp - TC-1000-01 1 Instrument Clamp - IC-1000-01 2 Vertical Bars Standard Length - BC-1030-01 & BC-1050-01
A-1001	Extended FastClamp™ System Includes: 1 Table Clamp - TC-1000-01 1 Extended Instrument Clamp - IC-1001-01 2 Vertical Bars - BC-1030-01 & BC-1090-01

##### FastClamp™ Table Clamp Options

CODE	DESCRIPTION
TC-1000-01	Table Clamp
IC-1000-01	Instrument Clamp
IC-1001-01	Extended Instrument Clamp

##### Bar Options

CODE	DESCRIPTION
BC-1030-01	30mm Bar
BC-1050-01	50mm Bar
BC-1070-01	70 mm Bar
BC-1090-01	90mm Bar