



Introducing the Symmetry Secto® Gauze Dissectors

- Enhanced design and durability for sturdy dissection and retraction
- Mix and match with a variety of reusable handles
- Efficiency through versatility of dissectors



Symmetry Secto®
Gauze Dissectors

Symmetry Secto[®] Gauze Dissectors

Ideal for dissection, retraction, port cleaning and probing



Durable

- Up to seven layers of 100% surgical cotton gauze for high-capacity fluid absorption

Only Pay for What You Use

- Each is packaged individually sterile and ready for use

Perfect for Port Cleaning

- Both the 5mm and 10mm one-piece single-use gauze dissectors fit flush inside the port lumen, giving a clean wipe with every pass

Efficiency Through Versatility

- Streamline inventory by consolidating to one dissector that can be used for both open and laparoscopic procedures

Waste Less, Spend Less

- Innovative reusable handle saves 80% of waste

Controlled Leverage

- All dissector handles are made of either stainless steel (reusable) or anodized aluminum (single-use) for sturdy dissection and retraction

Mix and Match

- The same disposable tips thread onto a variety of reusable handles, including: laparoscopic and open types, short and long lengths, and 5mm and 10mm widths

Universal Disposable Tips

Each disposable tip is packaged individually sterile and threads securely onto any corresponding reusable handle.



5mm Disposable Tip

Aluminum frame wrapped in 5 layers of 100% surgical cotton gauze

CODE: SG5



10mm Disposable Tip

Aluminum frame wrapped in 7 layers of 100% surgical cotton gauze

CODE: SG10

Symmetry Secto® Gauze Dissectors

Reusable Handles

Each reusable handle is made with stainless steel that can be steam sterilized for use with corresponding disposable tips.

Short Reusable Handles



5mm (For use with SG5 Disposable Tips)
CODE: GDH5

10mm (For use with SG10 Disposable Tips)
CODE: GDH10



5mm (For use with SG5 Disposable Tips)
CODE: SGH532

10mm (For use with SG10 Disposable Tips)
CODE: SGH1032

Long Reusable Handles



5mm (For use with SG5 Disposable Tips)
CODE: LDH5

10mm (For use with SG10 Disposable Tips)
CODE: LDH10



5mm (For use with SG5 Disposable Tips)
CODE: SGH5

10mm (For use with SG10 Disposable Tips)
CODE: SGH10

Secto® One-Piece Single-Use Gauze Dissectors













Disposable gauze dissector tips fused onto lightweight aluminum shafts that are packaged individually sterile and ready for immediate use.



5mm
CODE: TDSG5

10mm
CODE: TDSG10

Symmetry Secto® Gauze Dissectors

Product	Description	Code	Reusable	Single Use
	Secto® 5mm Gauze Dissector Single-Use Tip, Working Length 2.5cm, Box/10	SG5		X
	Secto® 10mm Gauze Dissector Single-Use Tip, Working Length 3cm, Box/10	SG10		X
	Secto® 5mm Reusable Handle, Short, Working Length 20.5cm	SGH532	X	
	Secto® 10mm Reusable Handle, Short, Working Length 20.5cm	SGH1032	X	
	Secto® 5mm Laparoscopic Reusable Handle, Long, Working Length 33cm	SGH5	X	
	Secto® 10mm Laparoscopic Reusable Handle, Long, Working Length 33cm	SGH10	X	
	Secto® 5mm Reusable Handle, Short, Angle Tip, Working Length 8.5cm	GDH5	X	
	Secto® 10mm Reusable Handle, Short, Angle Tip, Working Length 8.5cm	GDH10	X	
	Secto® 5mm Reusable Handle, Long, Angle Tip, Working Length 17.5cm	LDH5	X	
	Secto® 10mm Reusable Handle, Long, Angle Tip, Working Length 17.5cm	LDH10	X	
	Secto® 5mm One-Piece Single-Use Gauze Dissector, Overall Length 35cm, Box/10	TDSG5		X
	Secto® 10mm One-Piece Single-Use Gauze Dissector, Overall Length 35cm, Box/10	TDSG10		X

Visit symmetrysurgical.com or contact your local Sales Representative for more information.



3034 Owen Drive | Antioch, TN 37013 | 615.883.9090 | 800.251.3000 | customerservice@symmetrysurgical.com | symmetrysurgical.com

©2016 Symmetry Surgical, Inc. All rights reserved. SM0065