





Introducing the Olsen® Precision Single-Use Non-Stick Bipolar Forceps


- Engineered for balance, precision, visibility and thermal control to provide premium performance
- Available in irrigating and non-irrigating designs
- Non-stick Cermet coating technology



Olsen® Precision Single-Use Non-Stick Bipolar Forceps

Olsen® Precision Single-Use Non-Stick Bipolar Forceps

Designed for balance, precision, visibility and thermal control



Visit symmetrysurgical.com or contact your local Sales Representative for more information.

Olsen® Precision Non-Stick Bipolar Forceps

CODE	OVERALL LENGTH	WORKING LENGTH	TIP SIZE	QTY
20-1571IMT	7 in (18 cm)	2 1/2 in (6.5 cm)	0.5 mm	5/box
20-1572IMT	7 in (18 cm)	2 1/2 in (6.5 cm)	1.0 mm	5/box
20-1573IMT	7 in (18 cm)	2 1/2 in (6.5 cm)	1.5 mm	5/box
20-1581IMT	8 in (20 cm)	3 1/2 in (9.0 cm)	0.5 mm	5/box
20-1582IMT	8 in (20 cm)	3 1/2 in (9.0 cm)	1.0 mm	5/box
20-1583IMT	8 in (20 cm)	3 1/2 in (9.0 cm)	1.5 mm	5/box
20-1591IMT	9 in (23 cm)	4 1/2 in (11.5 cm)	0.5 mm	5/box
20-1592IMT	9 in (23 cm)	4 1/2 in (11.5 cm)	1.0 mm	5/box
20-1593IMT	9 in (23 cm)	4 1/2 in (11.5 cm)	1.5 mm	5/box

Olsen® Precision Non-Stick Bipolar Forceps, Irrigating Tip

CODE	OVERALL LENGTH	WORKING LENGTH	TIP SIZE	QTY
20-0571IMT	7 in (18 cm)	2 1/2 in (6.5 cm)	0.5 mm	5/box
20-0572IMT	7 in (18 cm)	2 1/2 in (6.5 cm)	1.0 mm	5/box
20-0573IMT	7 in (18 cm)	2 1/2 in (6.5 cm)	1.5 mm	5/box
20-0581IMT	8 in (20 cm)	3 1/2 in (9.0 cm)	0.5 mm	5/box
20-0582IMT	8 in (20 cm)	3 1/2 in (9.0 cm)	1.0 mm	5/box
20-0583IMT	8 in (20 cm)	3 1/2 in (9.0 cm)	1.5 mm	5/box
20-0591IMT	9 in (23 cm)	4 1/2 in (11.5 cm)	0.5 mm	5/box
20-0592IMT	9 in (23 cm)	4 1/2 in (11.5 cm)	1.0 mm	5/box
20-0593IMT	9 in (23 cm)	4 1/2 in (11.5 cm)	1.5 mm	5/box

Bipolar Cords & Tubing

CODE	DESCRIPTION	QTY
40-1102	12 ft. (3.6M) Bipolar Cord, Single-Use, Twin Pin Generator Plug, Sterile	10/box
40-1103	12 ft. (3.6M) Bipolar Cord, Single-Use, Twin Pin Generator Plug, Sterile	50/box
40-1105-28	12 ft. (3.6M) Bipolar Cord, Single-Use, 28mm Fixed Pin (Unitized) Generator Plug, Sterile	10/box
40-1105-39	12 ft. (3.6M) Bipolar Cord, Single-Use, 39 mm Fixed Pin (Unitized) Generator Plug, Sterile	10/box
40-3001R	12 ft. (3.6M) Integrated Irrigation Tubing & Bipolar Cord Set, Single-Use, On/Off Clamp, Rotary Pump, Sterile	10/box
40-3003R-39	12 ft. (3.6M) Integrated Irrigation Tubing & Bipolar Cord Set, Single-Use, On/Off Clamp, Rotary Pump, 39mm Fixed Pin (Unitized) Generator Plug, Sterile	10/box



3034 Owen Drive | Antioch, TN 37013 | 615.883.9090 | 800.251.3000 | customerservice@symmetrysurgical.com | symmetrysurgical.com

©2016 Symmetry Surgical, Inc. All rights reserved. SM 0061