

The Genuine Bookwalter® Magrina Retractor System

- The latest Bookwalter® innovations are designed to enhance function and stability
- Optimized for vaginal hysterectomy



Bookwalter® Magrina
Retractor System



“Good exposure is the key to good surgery.”

-John Bookwalter, M.D.

- Retractor blades remain on ring without tissue resistance
- Better visualization because ring conforms to patient anatomy
- Lightweight components reduce set weight
- Compatible with all genuine Bookwalter® components
- More than 30 years of innovation, performance and reliability
- Patented portfolio is surgeon designed and developed



CODE	DESCRIPTION	QTY
50-4692	Bookwalter® Magrina Kit includes the following	
50-4688	Bookwalter® 12" Lightweight Table Post	1
50-4554	Bookwalter® Post Coupling	1
50-4582	Bookwalter® Lightweight Horizontal Flex Bar	1
50-4689	Bookwalter® Magrina Segmented Ring w/ Suture Slots	1
50-4686	Bookwalter® Self-Retaining Tilt Ratchet	4
50-4591	Bookwalter® Magrina Lateral Retractor Blade	2
50-4592	Bookwalter® Magrina Small Deaver Retractor Blade	1
50-4593	Bookwalter® Magrina Large Deaver Retractor Blade	1
50-4594	Bookwalter® Magrina Posterior Retractor Blade	1

