

# Laparoscopy Catalog



# Table of Contents

Laparoscopic Portfolio Overview	4
Symmetry VersaGrip™ Instruments	5-9
• Handles	6
• Scissors	6
• Dissectors	7
• Graspers	8-9
RapidClean® Modular Instruments	10-20
• Scissors	12
• Dissectors	13
• Graspers	14-16
• Replacement	17
• Ordering Additional Inserts	18-20
Symmetry Surgical® One-Piece Instruments	21-33
• Scissors	25
• Dissectors	26-27
• Graspers	28-33
Energy Instruments	34-38
• Graspers/Kleppinger	35
• Cords	36
• Electrodes	37-38
• Accessories	38
Needle Holders & Knot Pushers	39-41
• Needle Holders	40
• Knot Pushers	41
• Accessories	41
Symmetry Surgical® Laparoscopic Accessories	42-48
• Access Accessories	43-44
• Intervention Accessories	46-47
• Single-Use Accessories	48
• Miscellaneous	48
Customer Support	49

# Laparoscopic Portfolio

## Overview

We designed our laparoscopic solutions specifically to meet patient and surgeon needs - focusing on high quality, versatile and durable instrumentation to inspire the highest levels of confidence. Our broad portfolio of ergonomically designed, single-piece and modular laparoscopic instruments and accessories are designed to enable minimally invasive surgery in laparoscopic specialties, including colorectal, gynecologic, bariatric, general and urologic surgery.






# Laparoscopic Portfolio




## Introducing Symmetry VersaGrip™

Symmetry VersaGrip™ is the only custom fit laparoscopic instrument designed to move with you. Symmetry VersaGrip™ includes a two-piece modular design and compatibility to fit dissectors, scissors and graspers — a versatile laparoscopic handle for every laparoscopic surgery.







# Symmetry VersaGrip™

## Handles

CODE	SHAFT Diameter: 5mm	Handle	
VG-1000-S	33cm	Symmetry VersaGrip™ Handle, Small, non-ratcheted	
VG-1000-45-S	45cm		
VG-1000-M	33cm	Symmetry VersaGrip™ Handle, Medium, non-ratcheted	
VG-1000-45-M	45cm		
VG-1000-L	33cm	Symmetry VersaGrip™ Handle, Large, non-ratcheted	
VG-1000-45-L	45cm		

## Scissors

CODE Insert Only	CODE Small	CODE Medium	CODE Large	SHAFT Diameter: 5mm	TIP	
VG-2020Q	VG-2020-S	VG-2020-M	VG-2020-L	33cm	Hook Scissor	
VG-2020Q-45	VG-2020-45-S	VG-2020-45-M	VG-2020-45-L	45cm		
VG-2045Q	VG-2045-S	VG-2045-M	VG-2045-L	33cm	Metzenbaum Scissors, Straight	
VG-2045Q-45	VG-2045-45-S	VG-2045-45-M	VG-2045-45-L	45cm		
VG-2044Q	VG-2044-S	VG-2044-M	VG-2044-L	33cm	Mini Metzenbaum Scissors, Curved	
VG-2044Q-45	VG-2044-45-S	VG-2044-45-M	VG-2044-45-L	45cm		
VG-2006Q	VG-2006-S	VG-2006-M	VG-2006-L	33cm	Endo-Metzenbaum Scissors, Curved	
VG-2006Q-45	VG-2006-45-S	VG-2006-45-M	VG-2006-45-L	45cm		










# Symmetry VersaGrip™

## Dissectors

CODE Insert Only	CODE Small	CODE Medium	CODE Large	SHAFT Diameter: 5mm	TIP	
VG-1030Q	VG-1030-S	VG-1030-M	VG-1030-L	33cm	Dolphin (Bottle Nose) Dissector	
VG-1030Q-45	VG-1030-45-S	VG-1030-45-M	VG-1030-45-L	45cm		
VG-1067Q	VG-1067-S	VG-1067-M	VG-1067-L	33cm	Maryland Dissector, Short, Bent Left	
VG-1067Q-45	VG-1067-45-S	VG-1067-45-M	VG-1067-45-L	45cm		
VG-1064Q	VG-1064-S	VG-1064-M	VG-1064-L	33cm	Kelly Dissector, Left Bent	
VG-1064Q-45	VG-1064-45-S	VG-1064-45-M	VG-1064-45-L	45cm		
VG-1066Q	VG-1066-S	VG-1066-M	VG-1066-L	33cm	Mixer Dissector	
VG-1066Q-45	VG-1066-45-S	VG-1066-45-M	VG-1066-45-L	45cm		
VG-1068Q	VG-1068-S	VG-1068-M	VG-1068-L	33cm	Mixer, Angle Combo, Long	
VG-1068Q-45	VG-1068-45-S	VG-1068-45-M	VG-1068-45-L	45cm		
VG-1060Q	VG-1060-S	VG-1060-M	VG-1060-L	33cm	Tube Dissector	
VG-1060Q-45	VG-1060-45-S	VG-1060-45-M	VG-1060-45-L	45cm		
VG-1164Q	VG-1164-S	VG-1164-M	VG-1164-L	33cm	Tissue Forceps	
VG-1164Q-45	VG-1164-45-S	VG-1164-45-M	VG-1164-45-L	45cm		




# Symmetry VersaGrip™

## Graspers

CODE Insert Only	CODE Small	CODE Medium	CODE Large	SHAFT Diameter: 5mm	TIP	
VG-1035Q	VG-1035-S	VG-1035-M	VG-1035-L	33cm	Grasping Forcep with Reservoir	
VG-1035Q-45	VG-1035-45-S	VG-1035-45-M	VG-1035-45-L	45cm		
VG-1092Q	VG-1092-S	VG-1092-M	VG-1092-L	33cm	Debaquey Grasper	
VG-1092Q-45	VG-1092-45-S	VG-1092-45-M	VG-1092-45-L	45cm		
VG-1033Q	VG-1033-S	VG-1033-M	VG-1033-L	33cm	Grasper With Cross-Serrations	
VG-1033Q-45	VG-1033-45-S	VG-1033-45-M	VG-1033-45-L	45cm		
VG-1032Q	VG-1032-S	VG-1032-M	VG-1032-L	33cm	Aggressive Grasper	
VG-1032Q-45	VG-1032-45-S	VG-1032-45-M	VG-1032-45-L	45cm		
VG-1022Q	VG-1022-S	VG-1022-M	VG-1022-L	33cm	Duckbill Grasper	
VG-1022Q-45	VG-1022-45-S	VG-1022-45-M	VG-1022-45-L	45cm		
VG-1008Q	VG-1008-S	VG-1008-M	VG-1008-L	33cm	Tenaculum Grasper	
VG-1008Q-45	VG-1008-45-S	VG-1008-45-M	VG-1008-45-L	45cm		
VG-1049Q	VG-1049-S	VG-1049-M	VG-1049-L	33cm	Claw Retracting Forceps	
VG-1049Q-45	VG-1049-45-S	VG-1049-45-M	VG-1049-45-L	45cm		
VG-2070Q	VG-2070-S	VG-2070-M	VG-2070-L	33cm	Park Grasper	
VG-2070Q-45	VG-2070-45-S	VG-2070-45-M	VG-2070-45-L	45cm		
VG-1042Q	VG-1042-S	VG-1042-M	VG-1042-L	33cm	Schertl Fenestrated Grasper	
VG-1042Q-45	VG-1042-45-S	VG-1042-45-M	VG-1042-45-L	45cm		

# Symmetry VersaGrip™

## Graspers


CODE Insert Only	CODE Small	CODE Medium	CODE Large	SHAFT Diameter: 5mm	TIP	
VG-1050Q	VG-1050-S	VG-1050-M	VG-1050-L	33cm	Fundus Grasper	
VG-1050Q-45	VG-1050-45-S	VG-1050-45-M	VG-1050-45-L	45cm		
VG-1016Q	VG-1016-S	VG-1016-M	VG-1016-L	33cm	Alligator Grasper	
VG-1016Q-45	VG-1016-45-S	VG-1016-45-M	VG-1016-45-L	45cm		
VG-1069Q	VG-1069-S	VG-1069-M	VG-1069-L	33cm	Cobra Grasper	
VG-1069Q-45	VG-1069-45-S	VG-1069-45-M	VG-1069-45-L	45cm		
VG-1031Q	VG-1031-S	VG-1031-M	VG-1031-L	33cm	Dorsey Fenestrated Grasper	
VG-1031Q-45	VG-1031-45-S	VG-1031-45-M	VG-1031-45-L	45cm		
VG-1043Q	VG-1043-S	VG-1043-M	VG-1043-L	33cm	Single Action Wave Grasper	
VG-1043Q-45	VG-1043-45-S	VG-1043-45-M	VG-1043-45-L	45cm		
VG-1079Q	VG-1079-S	VG-1079-M	VG-1079-L	33cm	Babcock Grasper	
VG-1079Q-45	VG-1079-45-S	VG-1079-45-M	VG-1079-45-L	45cm		
VG-1041Q	VG-1041-S	VG-1041-M	VG-1041-L	33cm	Double Action Wave Grasper	
VG-1041Q-45	VG-1041-45-S	VG-1041-45-M	VG-1041-45-L	45cm		
VG-1087Q	VG-1087-S	VG-1087-M	VG-1087-L	33cm	Allis Grasper	
VG-1087Q-45	VG-1087-45-S	VG-1087-45-M	VG-1087-45-L	45cm		
VG-1093Q	VG-1093-S	VG-1093-M	VG-1093-L	33cm	Hunter Grasper	
VG-1093Q-45	VG-1093-45-S	VG-1093-45-M	VG-1093-45-L	45cm		
VG-1089Q	VG-1089-S	VG-1089-M	VG-1089-L	33cm	Teardrop Babcock Grasper	
VG-1089Q-45	VG-1089-45-S	VG-1089-45-M	VG-1089-45-L	45cm		



# Laparoscopic Portfolio

## Modular Instruments

Symmetry Surgical's RapidClean® Modular instrument line has more than 50 patterns of graspers, dissectors and scissors. All instruments are reusable and detachable, able to disassemble into three components: Carbon Fiber handle + Insulated shaft with flush port and 360° rotation knob + interchangeable insert/tip.



# Features & Benefits

## RapidClean®




### INTERCHANGEABLE TIPS

Cost effective replacement and versatility. More than 50 patterns of jaw inserts.



### 360° ROTATION KNOB

Controlled precision. Easily rotate the shaft for optimum jaw position.



### FLUSH PORT

Allows for optimal cleaning of the instrument.



### INSULATED SHAFT

Durable and high dielectric strength polymer shaft. Added protection against thermal damage.



### 3-PIECE DETACHABLE

Intuitive assembly. Quick and easy disassembly.



### 45° CAUTERY POST

Reduces electrocautery cord interference. 100% monopolar product line.



### ERGONOMIC DESIGN CARBON FIBER HANDLE

To increase comfort during an operation and reduce hand fatigue. Interchangeable based on procedure/function. Non-conductive for safety.


# Modular Instruments

## Scissors

CODE	SHAFT Diameter: 5mm	TIP	HANDLE
96-2003NL	33cm	Metzenbaum, Curved, Double Action	 Non-locking handle
96-2003NL-45	45cm		
96-2006NL	33cm	Metzenbaum, Curved, Slimline Scissors, Double Action	
96-2006NL-45	45cm		
96-2045NL	33cm	Metzenbaum Straight Scissors, Double Action	
96-2020NL	33cm	Hook Scissors, Single Action	
96-2020NL-45	45cm		
96-2044NL-45	45cm	Micro Scissors, Extremely Delicate Tips, Single Action	
96-1145NL	33cm	Metzenbaum, Curved Serrated Scissors, Double Action	
96-11137	33cm	Micro Scissors, Extremely Delicate Tips, Single Action	
96-11117	45cm	Metzenbaum, Curved, Serrated Scissors, Flat Dissecting Blades, Double Action	

# Modular Instruments



## Dissectors

CODE	SHAFT Diameter: 5mm	TIP	HANDLE
96-1022NL	33cm	Duckbill/Tapered Forceps, Cross Serrated, Serrated, Double Action	 Non-locking handle
96-1030NL	33cm	Bottle Nose (Dolphin) Dissector, Short Jaw, Double Action	
96-1030NL-45	45cm	Bottle Nose (Dolphin) Dissector, Short Jaw, Double Action	
96-1033NL	33cm	Bullet Nose Dissector, With Cup, Double Action	
96-1066NL	33cm	Mixer Dissector, 90° Angle, Double Action	
96-1066NL-45	45cm	Mixer Dissector, 90° Angle, Double Action	
96-1067NL	33cm	Maryland Dissector, Double Action	
96-1067NL-45	45cm	Maryland Dissector, Double Action	
96-1167NL	33cm	Maryland Dissector, Diamond Serrations, Double Action	

\*RapidClean® Instruments also available with Traditional Ratchet Locking Handle and Locking Handle with Ratchet-On-Demand

# Modular Instruments

## Graspers

CODE	SHAFT Diameter: 5mm	TIP	HANDLE
96-1008R	33cm	Tenaculum Forceps, Traumatic, 1x1 Sharp Teeth,	 Locking Handle with Ratchet on Demand
96-1008R-45	45cm	Double Action	
96-1015R	33cm	Claw Grasper With 2x3 Teeth, Traumatic, Single Action	
96-1016R	33cm	Alligator Grasper, Traumatic, Double Action	
96-1016R-45	45cm		
96-1018R	33cm	Fenestrated Grasper, Micro Jaw, Double Action	
96-1020R	33cm	Strong, Traumatic Grasper, Double Action	
96-1020R-45	45cm		
96-1024R	33cm	Serrated Grasper, Rounded Jaws, Double Action	
96-1035R	33cm	Cup Tooth Grasper, Double Action	
95-1035R-45	45cm		
96-1042R	33cm	Fenestrated Grasper, Slightly, Flat Jaws, Atraumatic, Double Action	
96-1042NL	33cm	Fenestrated Grasper, Slightly, Flat Jaws, Atraumatic, Double Action, Non-locking Handle	

\*RapidClean® Instruments also available with Non-Locking Handle and Traditional Ratchet Locking Handle



# Modular Instruments


## Graspers

CODE	SHAFT Diameter: 5mm	TIP	HANDLE
96-1050R	33cm	Fundus Grasper, Slightly Flat Jaws, Atraumatic, Double Action	 Locking Handle with Ratchet on Demand
96-1050R-45	45cm		
96-1051R	33cm	Wave Grasper with Atraumatic Serrations, Single Action	
96-1051R-45	45cm		
96-1053R	33cm	Wave Grasper with Atraumatic Serration, Double Action	
96-1053R-45	45cm		
96-1056R	33cm	Retraction Grasper with Atraumatic Teeth, Single Action	
96-1069R	33cm	Snake Tooth (Cobra) Traumatic Grasper, Double Action, 2x2 Teeth	
96-1069R-45	45cm		
96-1070R	33cm	Biopsy Spoon Grasper, Single Action	
96-1070R-45	45cm		
96-1071R	33cm	Biopsy Spoon, with Spike, Double Action	
96-1078R	33cm	Allis Atraumatic Grasper, 3x4 Teeth, Double Action	
96-1079R	33cm	Babcock Grasper, Atraumatic, Double Action	

\*RapidClean® Instruments also available with Non-Locking Handle and Traditional Ratchet Locking Handle

# Modular Instruments






## Graspers

CODE	SHAFT Diameter: 5mm	TIP	HANDLE
96-1089R	33cm	Babcock Grasper, Atraumatic, Double Action	 Locking Handle with Ratchet on Demand
96-1089R-45	45cm		
96-1092R	33cm	DeBakey Straight Grasper, Atraumatic, 40mm Jaw Length, Double Action	
96-1092R-45	45cm		
96-1093R	33cm	Hunter Bowel Grasper, Double Action	

\*RapidClean® Instruments also available with Non-Locking Handle and Traditional Ratchet Locking Handle

# Modular Instruments

## Replacement Components

CODE	SHAFT	DESCRIPTION	
96-0526	33cm	Sheath with flush port (Diameter: 5mm)	
96-0526-45	45cm		
96-0126	33cm	Sheath with flush port (Diameter 10mm)	
96-0126-45	45cm		
96-0002NL		Non-Locking Handle	
96-0003R		Locking Handle with Ratchet on Demand	
96-0004TR		Traditional Ratchet Locking Handle	

# Modular Instruments

## Ordering Additional Inserts

Use the codes listed below to order additional inserts for our modular instrument line. Use the full instrument code in the left-hand column as a reference to help find your needed insert code.

FULL INSTRUMENT CODE	CODE FOR INSERT ONLY	SHAFT	DIAM	TYPE	DESCRIPTION
96-2003NL	96-2003Q	33cm	5mm	Scissors	Metzenbaum Curved Scissors, Double Action
96-2003NL-45	96-2003Q-45	45cm			
96-2006NL	96-2006Q	33cm	5mm	Scissors	Metzenbaum Curved Slim-Line Scissors, Double Action
96-2006NL-45	96-2006Q-45	45cm			
96-2045NL	96-2045Q	33cm	5mm	Scissors	Metzenbaum Straight Scissors, Double Action
96-2020NL	96-2020Q	33cm	5mm	Scissors	Hook Scissors, Single Action
96-2020NL-45	96-2020Q-45	45cm			
96-2044NL-45	96-2044Q-45	45cm	5mm	Scissors	Micro Scissors, Extremely Delicate Tips, Single-Action
96-1145NL	96-1145Q	33cm	5mm	Scissors	Metzenbaum Curved Serrated Scissors, Double Action
96-11137	96-11137Q	33cm	5mm	Scissors	Micro Scissors, Extremely Delicate Tips, Single Action
96-11117	96-11117Q	45cm	5mm	Scissors	Metzenbaum Curved Serrated Scissors, Flat Dissecting Blades, Double Action
96-1022NL	96-1022Q	33cm	5mm	Dissectors	Duckbill/Tapered Forceps, Cross Serrated, Double Action
96-1030NL	96-1030Q	33cm	5mm	Dissectors	Bottle Nose (Dolphin) Dissector, With Cup, Short Jaw, Double Action
96-1030NL-45	96-1030Q-45	45cm			
96-1033NL	96-1033Q	33cm	5mm	Dissectors	Bullet Nose Dissector, With Cup, Double Action
96-1066NL	96-1066Q	33cm	5mm	Dissectors	Mixer Dissector, 90° Angle, Double Action
96-1066NL-45	96-1066Q-45	45cm			
96-1167NL	96-1167Q	33cm	10mm	Dissectors	Maryland Dissector, Diamond Serrations, Double Action
96-1067NL	96-1067Q	33cm	5mm	Dissectors	Maryland Dissector, Double Action
96-1067NL-45	96-1067Q-45	45cm			

# Modular Instruments

## Ordering Additional Inserts

Use the codes listed below to order additional inserts for our modular instrument line. Use the full instrument code in the left-hand column as a reference to help find your needed insert code.

FULL INSTRUMENT CODE	CODE FOR INSERT ONLY	SHAFT	DIAM	TYPE	DESCRIPTION
96-1008R	96-1008Q	33cm	5mm	Graspers	Tenaculum Forceps, Traumatic, 1x1 Sharp Teeth, Double Action
96-1008R-45	96-1008Q-45	45cm			
96-1015R	96-1015Q	33cm	10mm	Graspers	Claw Grasper With 2x3 Teeth, Traumatic, Single Action
96-1016R	96-1016Q	33cm	5mm	Graspers	Alligator Grasper, Traumatic, Double Action
96-1016R-45	96-1016Q-45	45cm			
96-1018R	96-1018Q	33cm	5mm	Graspers	Fenestrated Grasper, Micro Jaw, Double Action
96-1020R	96-1020Q	33cm	5mm	Graspers	Strong, Traumatic Grasper, Double Action
96-1020R-45	96-1020Q-45	45cm			
96-1024R	96-1024Q	33cm	5mm	Graspers	Serrated Grasper, Rounded Jaws, Double Action
96-1035R	96-1035Q	33cm	5mm	Graspers	Cup Tooth Grasper, Double Action
96-1035R-45	96-1035Q-45	45cm			
96-1042NL	96-1042Q	33cm	5mm	Graspers	Fenestrated Grasper, Slightly Flat Jaws, Atraumatic, Double Action
96-1042R					
96-1050R	96-1050Q	33cm	5mm	Graspers	Fundus Grasper, Slightly Flat Jaws, Atraumatic, Double Action
96-1050R-45	96-1050Q-45	45cm			
96-1051R	96-1051Q	33cm	5mm	Graspers	Wave Grasper With Atraumatic Serrations, Double Action
96-1051R-45	96-1051Q-45	45cm			
96-1053R	96-1053Q	33cm	5mm	Graspers	Wave Grasper With Atraumatic Serrations, Double Action
96-1053R-45	96-1053Q-45	45cm			
96-1056R	96-1056Q	33cm	5mm	Graspers	Retraction Grasper With Atraumatic Teeth, Single Action
96-1069R	96-1069Q	33cm	5mm	Graspers	Snake Tooth (Cobra) Traumatic Grasper, Double Action, 2x2 Teeth
96-1069R-45	96-1069Q-45	45cm			



# Modular Instruments

## Ordering Additional Inserts

Use the codes listed below to order additional inserts for our modular instrument line. Use the full instrument code in the left-hand column as a reference to help find your needed insert code.

FULL INSTRUMENT CODE	CODE FOR INSERT ONLY	SHAFT	DIAM	TYPE	DESCRIPTION
96-1070R	96-1070Q	33cm	5mm	Graspers	Biopsy Spoon Grasper, Single Action
96-1070R-45	96-1070Q-45	45cm			
96-1071R	96-1071Q	33cm	5mm	Graspers	Biopsy Spoon, With Spike, Double Action
96-1078R	96-1078Q	33cm	5mm	Graspers	Allis Atraumatic Grasper, 3x4 Teeth, Double Action
96-1079R	96-1079Q	33cm	10mm	Graspers	Babcock Grasper, Atraumatic, Double Action
96-1089R	96-1089Q	33cm	5mm	Graspers	Babcock Grasper, Atraumatic, Double Action
96-1089R-45	96-1089Q-45	45cm			
96-1092R	96-1092Q	33cm	5mm	Graspers	DeBakey Straight Grasper, Atraumatic, 40mm Jaw Length, Double Action
96-1092R-45	96-1092Q-45	45cm			
96-1093R	96-1093Q	33cm	5mm	Graspers	Hunter Bowel Grasper, Double Action

# Laparoscopic Portfolio

---


## One-Piece Instruments

With multiple handle configuration options and a wide range of tips, our one-piece laparoscopic instrument line continues our focus on quality, variety, durability and confidence to fully meet surgeon's handheld instrument needs. Made by manufacturers in Germany and the USA.



# Features & Benefits

## One-Piece Instruments: Slide Lock 95-Series



### ONE-PIECE REUSABLE

Ease of use, no need to assemble. Rigid, stable and solid. Even pressure distribution.



### INSULATED SHAFT

Durable and high dielectric strength polymer shaft. Added protection against thermal damage.



### ROTATION KNOB OPTION

Ergonomic feature to easily rotate the shaft for optimum jaw position.



### SLIDE LOCK MECHANISM

Precise and ergonomic ratcheting mechanism. Secure grasping and retraction. Locking instruments are available with the "SLIDE LOCK" locking feature. This unique locking device provides dual action for the surgeon's grasping needs in terms of durability, comfort and strength



### ADDITIONAL FEATURES

All instruments are designed with minimal moving parts and are made in the USA.



### STAINLESS STEEL HANDLE

Made with surgical grade stainless steel for maximum precision durability.



### FLUSH PORT

Allows for optimal cleaning of the instrument.


# Features & Benefits

## One-Piece Instruments: Direct Drive 95-Series



### ONE-PIECE REUSABLE

Ease of use, no need to assemble. Rigid, stable and solid. Even pressure distribution.



### FLUSH PORT

Allows for optimal cleaning of the instrument.



### SPECIAL "NON-STICK" COATING

Special coating on the jaws provides "non-stick" properties and facilitates wiping and cleaning. Additionally, it provides non-reflective properties and increased mechanical wear properties



### RATCHET ON DEMAND LOCKING SYSTEM

Locking instruments are available with the "RATCHET ON DEMAND" locking feature. This unique locking device provides dual action for the surgeon's dissection and grasping needs.

Handle Variety: Slide Lock, Ratchet on Demand, Non-locking insulated, Non-locking Non-cautery options available.



### ADDITIONAL FEATURES

All instruments are designed with minimal moving parts and are made with surgical grade stainless steel for maximum durability. Made in the USA.



### ROTATION KNOB OPTION

12 stop location at 360° rotation option



### PATENTED DIRECT DRIVE SYSTEM

All instruments include the patented "HINGELESS" direct drive system. This system is dependable and has no small hinges or links. Tests validated it to be stronger in grasping and dissection mode than conventional hinged systems.

# Features & Benefits

## One-Piece Instruments: Accu 97-Series



### ONE-PIECE REUSABLE

Ease of use, no need to assemble. Rigid, stable and solid. Even pressure distribution.



### INSULATED SHAFT

Durable and high dielectric strength polymer shaft. Added protection against thermal damage.



### ONE-HANDED ACCURATE ROTATION

Controlled precision. Easily rotate the shaft for optimum jaw position. Add ratchet lock for graspers.



### 90° CAUTERY POST

Insulated and cautery post handles for dissectors and scissors and non-insulated handles for graspers.







### FLUSH PORT

Allows for optimal cleaning of the instrument.



# One-Piece Instruments

## Scissors

CODE	SHAFT Diameter: 5mm	TIP	HANDLE
95-2020-45	45cm	Hook Straight Scissors, Direct Drive, Single Action	 Non-locking handle w/ Rotation
97-2020RC	33cm	UltraCut Hook Scissors, Single Action, Insulated, Monopolar	 Non-locking handle w/ Rotation
97-2020RCL	45cm		
95-2006	33cm	Metzenbaum UltraCut Curved Scissors, Direct Drive, Double Action, Insulated	 Non-locking handle w/ Rotation
95-2006-45	45cm		
97-2006RC	33cm	UltraCut Metzenbaum Endo Scissors, 13mm Jaw, Double Action, Insulated, Monopolar	 Non-locking handle w/ Rotation
97-2006RCL	45cm		

# One-Piece Instruments


## Dissectors

CODE	SHAFT Diameter: 5mm	TIP	HANDLE
90-1030	36cm	Bottle Nose/Dolphin Nose Dissector, 16 mm Jaw, Double Action, Insulated, Monopolar	  Non-locking handle
90-1030RC	36cm	Bottle Nose/ Dolphin Nose Dissector, Double Action, Insulated, Monopolar	  Non-locking handle w/ Rotation
90-1030RCL	45cm		
90-1035R	32cm	Cup Tooth Dissector, 12 mm Jaw, Double Action, Insulated Shaft	  Ratchet handle
90-1035RL	45cm		
90-1055	32cm	Cup Tooth Dissector, 12 mm Jaw, Double Action, Insulated Shaft	  Ratchet handle
90-1035	35cm	Cup Tooth Dissector, 12mm Jaw, Double Action, Insulated, Monopolar	  Non-locking handle
90-1035L	45cm		
97-1030RC	36cm	Bottle Nose/Dolphin Nose Dissector, Double Action, Insulated, Monopolar	  Non-locking handle w/ Rotation
97-1030RCL	45cm		
97-1035RC	32cm	Cup Tooth Dissector, Double Action, Insulated, Monopolar	  Non-locking handle w/ Rotation
95-1067	33cm	Maryland Dissector Curved, Direct Drive, Double Action	  Non-locking handle w/ Rotation
95-1067-45	45cm		

# One-Piece Instruments












## Dissectors

CODE	SHAFT Diameter: 5mm	TIP	HANDLE
90-1067	34cm	Maryland Dissector, 23mm Jaw, Double Action, Insulated, Monopolar	 Non-locking handle
97-1067RC	33cm	Maryland Dissector, Double Action, Insulated, Monopolar	 Non-locking handle w/ Rotation
97-1067RCL	45cm		
97-1063RC	35cm	Maryland Tapered Dissector, Double Action, Insulated, Monopolar	 Non-locking handle w/ Rotation
90-1044	32cm	Micro Fenestrated Dissector, 13mm Jaw, Double Action, Insulated, Monopolar	 Monopolar, Non-Locking handle
98-1143	36cm	Mixer Dissector, Right Angle, 19mm Jaw, Double Action, Insulated Monopolar	 Non-locking handle w/ Rotation

CODE	SHAFT Diameter: 10mm	TIP	HANDLE
90-1068	32cm	Mixer Dissector, 90° Right Angle, Double Action, Insulated, Monopolar	 Non-locking handle

# One-Piece Instruments

## Graspers

CODE	SHAFT Diameter: 5mm	TIP	HANDLE
95-1052L	33cm	Alligator Grasper, Traumatic, Direct Drive, Double Action, Insulated, Monopolar	  Locking Ratchet on Demand w/Rotation
94-1078L-45	33cm	Allis Grasper, Atraumatic, Direct Drive, Double Action, Insulated, Monopolar	  Locking Ratchet on Demand w/Rotation
95-1078L	33cm	Allis Grasper, Atraumatic, Direct Drive, Double Action, Insulated, Monopolar	  Locking Ratchet on Demand w/Rotation
98-1115	37cm	Long Allis Grasper, Atraumatic, 28mm Jaw, Double Action, Insulated, Monopolar	  Non-locking handle w/Rotation
98-1188	36cm	Atraumatic Grasper, 22mm Jaw, Double Action	  Ratchet handle w/Rotation
98-1188L	45cm		
95-1049SL-DA	35cm	Atraumatic Retraction Grasper, 22mm Jaw, Double Action	  Slide Locking
95-1049SL-DA-ROT	35cm	Atraumatic Retraction Grasper, 22mm Jaw, Double Action	  Slide Locking w/Rotation
95-1049SL-DA-45	45cm	Atraumatic Retraction Grasper, 22mm Jaw, Double Action, Insulated Shaft	  Slide Locking

# One-Piece Instruments

## Graspers

CODE	SHAFT Diameter: 5mm	TIP	HANDLE
95-1049SL-DA-ROT-45	45cm	Atraumatic Retraction Grasper, 22mm Jaw, Double Action, Insulated Shaft	 Slide Locking w/Rotation
95-1048SL-ROT-45	45cm	Atraumatic Retraction Grasper, 22mm Jaw, Single Action	 Slide Locking w/Rotation
95-1049SL	33cm	Atraumatic Retraction Grasper, 22mm Jaw, Single Action	 Slide Locking
95-1049SL-45	45cm	Atraumatic Retraction Grasper, 22mm Jaw, Single Action	 Slide Locking w/Rotation
95-1049SL-ROT	35cm	Atraumatic Retraction Grasper, 22mm Jaw, Single Action	 Slide Locking w/Rotation
95-1049SL-ROT-45	45cm	Atraumatic Retraction Grasper, 22mm Jaw, Single Action, Insulated Shaft	 Slide Locking w/Rotation
95-1053L	33cm	Atraumatic Wave Grasper, Direct Drive, Double Action, Insulated, Monopolar	 Locking Ratchet on Demand w/Rotation
95-1048SL-DA-ROT	33cm	Traumatic Retraction Grasper, Serrated, 22mm Jaw, Double Action	 Slide Locking w/Rotation
95-1048SL-ROT	35cm	Traumatic Retraction Grasper, Serrated, 22mm Jaw, Single Action	 Slide Locking w/Rotation
95-1048SL-45	45cm	Traumatic Retraction Grasper, Serrated, 22mm Jaw, Single Action, Insulated Shaft	 Slide Locking



















# One-Piece Instruments

## Graspers

CODE	SHAFT Diameter: 5mm	TIP	HANDLE
97-1056RC	37cm	Traumatic Retraction Grasper, Single Action	  Ratchet Handle w/Rotation
95-1089L	33cm	Babcock Grasper, Atraumatic, Direct Drive, Double Action, Insulated, Monopolar	  Locking Ratchet on Demand w/Rotation
95-1089L-45	45cm		
90-1005	32cm	Claw Extractor Grasper, Traumatic, Heavy 2x3 Teeth, Single Action	  Spring Handle
95-1069L	33cm	Cobra Tooth Grasper, Traumatic, Direct Drive, Double Action	  Locking Ratchet on Demand w/Rotation
95-1069L-45	45cm		
95-1092L	33cm	DeBakey Straight Grasper, 40mm Jaw, Direct Drive, Double Action, Insulated, Monopolar	  Locking Ratchet on Demand w/Rotation
95-1092L-45	45cm		
97-1031RC	32cm	Dorsey Grasper, Fenestrated, 30mm Jaw, Double Action, Insulated Shaft	  Ratchet Handle w/Rotation
97-1031RCL	46cm		
95-1031SL	34.5cm	Dorsey Grasper, Fenestrated, 30mm Jaw, Single Action, Insulated Shaft	  Slide Locking
95-1042L	33cm	Fenestrated Grasper, Direct Drive, Double Action	  Locking Ratchet on Demand w/Rotation














# One-Piece Instruments

## Graspers

CODE	SHAFT Diameter: 5mm	TIP	HANDLE
90-1050	32cm	Fundus Trimax Grasper, 16mm Jaw, Double Action, Insulated Shaft	  Ratchet Handle
95-1050SL-DA-ROT	33.4cm	Fundus Trimax Grasper, 16mm Jaw, Double Action, Insulated Shaft	  Slide Locking w/ Rotation
98-1125	38cm	Glassman Grasper, 40mm Jaw, Double Action	  Ratchet Handle w/ Rotation
98-1125L	45cm		
95-1093L	33cm	Hunter Bowel Grasper, 20mm Jaw, Direct Drive, Double Action, Insulated, Monopolar	  Locking Ratchet on Demand w/ Rotation
97-1090RC	37cm	Hunter Bowel Grasper, 20mm Jaw, Double Action, Insulated Shaft	  Non-locking handle w/ Rotation
97-1090RCL	45cm		
98-1257	36cm	Spoon Grasper 13mm Jaw, Double Action	  Standard Handle
90-1008	32cm	Tenaculum Forceps, Traumatic, 1x1 Sharp Teeth, Double Action, Insulated Shaft	  Ratchet Handle
90-1008L	45cm		
97-1008RC	33cm	Tenaculum Grasper, Traumatic, Double Action	  Ratchet Handle



# One-Piece Instruments

## Graspers

CODE	SHAFT Diameter: 10mm	TIP	HANDLE	
90-1015	33cm	Claw Extractor Grasper, Traumatic, Heavy 2x3		 Spring Handle
90-1015L	45cm	Teeth, 29mm Jaw, Single Action		
98-1270	35.1cm	Claw Grasper, 29mm Jaw, Single Action		 Ratchet Handle w/Rotation
90-1083	34.5cm	Duval Grasper , Double Action		 Ratchet Handle
98-1274	35cm	Spoon Grasper, 25mm Jaw, Single Action		 Auto-Closure Handle w/Rotation
98-1277L	45cm	Spoon Grasper, 25mm Jaw, Single Action		 Standard Handle w/Rotation
90-1070	32cm	Spoon Grasper, Stone Retrieval Forceps, 25mm Jaw, Single Action		 Spring Handle
95-1070	33cm	Spoon Grasper, Stone Retrieval Forceps, Direct Drive, Double Action		 Non-locking handle w/Rotation
90-1007	32cm	Tenaculum Forceps, Traumatic, 1x1 Sharp Teeth, Double Action		 Ratchet Handle




# One-Piece Instruments

## Graspers

CODE	SHAFT Diameter: 10mm	TIP	HANDLE
90-1009-H	32cm	Tenaculum Forceps, Traumatic, 1x1 Sharp Teeth, Double Action, Insulated Shaft	 Ratchet Handle
90-1009-HX	32cm	Tenaculum Forceps, Traumatic, 1x1 Sharp Teeth, Double Action, Insulated Shaft	 Spring Handle
90-1009L	45cm	Tenaculum Forceps, Traumatic, 1x1 Sharp Teeth, Double Action, Insulated Shaft	

# NEW

## Graspers

CODE	TOTAL LENGTH	WORKING LENGTH	JAW LENGTH	JAW WIDTH	TIP	HANDLE
DG520	55cm	35cm	20mm	5mm	20 mm stainless steel jaws and PEEK handle	  

# Laparoscopic Portfolio

---

## Energy Instruments









Symmetry Surgical's Laparoscopic Energy offerings include monopolar and bipolar electrosurgery devices, reusable and single-use monopolar electrodes, bipolar forceps and accessories.



# Energy Instruments

## Graspers

CODE	LENGTH Diameter: 4.5mm	TIP	HANDLE
BC3310	34cm	Mectra Labs® Single Use Kleppinger-Style Bipolar Serrated Forceps w/ Cable	
BC3320	45cm	Mectra Labs® Single Use Kleppinger-Style Bipolar Serrated Forceps w/ Cable	

CODE	LENGTH	DESCRIPTION
90-1011		Symmetry Ultra™ Kleppinger Type Bipolar Forceps, Ring Handle 
90-1012		Symmetry Ultra™ Kleppinger Type Bipolar Forceps, Ring Handle 
90-1013		Symmetry Ultra™ Kleppinger Type Bipolar Forceps, Thumb Ring Handle 
90-1014		Symmetry Ultra™ Kleppinger Type Bipolar Forceps, Thumb Ring Handle 
90-1016	32cm	Insert, Kleppinger Type Waffle Tong 
90-1017	45cm	Insert, Kleppinger Type Waffle Tong Grasper 
90-1018	32cm	Insert, Kleppinger Type Paddle Tong 
90-1019	45cm	Insert, Kleppinger Type Paddle Tong 



# Energy Instruments

## Cords

CODE	LENGTH	DESCRIPTION
------	--------	-------------

90-1020	10ft	Bipolar Cable for SSI/Wolf Kleppinger Type Forceps & Wolf Generator, Reusable	
90-1021	10ft	Bipolar Cable for SSI/Wolf Kleppinger Type Forceps & Birtcher, Cabot & Valley Lab Type Generators, Reusable	
90018	10ft	Olsen® Monopolar Cord, Generator to Standard Olsen Instrument, Reusable	
90020	10ft	Olsen® Monopolar Cord, Valleylab Generator to Standard 8-32 UNC Olsen® Instrument, Reusable	
90056	10ft	Olsen® Monopolar Cord, Valleylab Style Generator Plug to 4mm Laparoscopic Instrument, Reusable	
40-2106	10ft	Olsen® Monopolar Cord, Bovie Plug to 4mm Socket Angled 90° (Valleylab® Style), Single Use, 10/Box	
40-2110	10ft	Olsen® Monopolar Cord, Bovie Plug to Straight 4.0mm Socket (Valleylab® Style), Single Use, 10/Box	
50-2107	10ft	Olsen® Monopolar Cord, Endoscopic Instrument, 90° Plug, Reusable	
50-2111	10ft	Olsen® Monopolar Cord, Endoscopic Instrument, Straight Plug, Reusable	




# Energy Instruments

## Electrodes

CODE	LENGTH	DESCRIPTION	
90-9522	33cm	L-Hook Electrode, Monopolar, Reusable	
90-9522L	45cm		
90-9523	32cm	Spatula Electrode, Monopolar, Reusable	
90-9523L	45cm		
90-9525	32cm	J-Hook Electrode, Monopolar, Reusable	
90-9525L	45cm		
90-9526	32cm	Needle Electrode, Monopolar, Reusable	
90-9526L	45cm		
90-9527	32cm	Button Tip Electrode, Monopolar, Reusable	
90-9527L	45cm		
90-9522-H	32cm	L-Hook Electrode with Taper	
90-4005	32cm	J-Hook Suction/Irrigation Electrode, with Valve, Insulated, Monopolar, Reusable	
90-4010	32cm	Button Tip Suction/Irrigation Electrode, with Valve, Insulated, Monopolar, Reusable	
90-4010L	45cm		
90-4012	32cm	L-Hook Suction/Irrigation Electrode, with Valve, Insulated, Monopolar, Reusable	
90-4012L	45cm		
90-4015	32cm	Spatula Suction/Irrigation Electrode, with Valve, Insulated, Monopolar, Reusable	
90-4015L	45cm		

# Energy Instruments

## Single-Use Electrodes

CODE	LENGTH	DESCRIPTION	
90-9522-S	33cm	L-Hook Electrode, Monopolar, Single-Use (10/box)	
90-9523-S	33cm	Spatula Electrode, Monopolar, Single-Use (10/box)	
90-9525-S	33cm	J-Hook Electrode, Monopolar, Single-Use (10/box)	

## Accessories

CODE	DESCRIPTION	IMAGE
90-1011SC	Sealing Cap for Kleppinger	

# Laparoscopic Portfolio

---


## Needle Holders & Knot Pushers

Symmetry Surgical's stainless steel and tungsten carbide needle holders and knot pushers provide durability, comfort and dexterity for suturing and knot tying.



# Needle Holders & Knot Pushers

## Needle Holders

CODE	LENGTH Diameter:5mm	TIP	HANDLE
68-1564	33cm	Tungsten Carbide Symmetry® Needle Holder, 30° Angle	  In-line Handle
90-2008	35cm	Tungsten Carbide Needle Holder, Right Curved	  In-line Handle
90-2010	35cm	Tungsten Carbide Needle Holder, Straight	  In-line Handle
90-2010L	45cm		
90-2011	35cm	Tungsten Carbide Needle Holder, Left Curved	  In-line Handle
90-2011L	45cm		
90-3035	32cm	Tungsten Carbide Ergonomic Needle	  Squeeze Handle with Ratchet
90-3038	32cm	Tungsten Carbide Needle Holder, Curved	  Squeeze Handle with Ratchet
90-3038L	45cm		
90-3039	32cm	Tungsten Carbide Needle Holder, Straight	  Squeeze Handle with Ratchet

# Needle Holders & Knot Pushers

## Knot Pushers

CODE	LENGTH	DESCRIPTION
90-5097	32cm	Knot Pusher, Close
90-5099	29.85cm	Knot Pusher, Open
90-5099L	45cm	



## Accessories

CODE	LENGTH Diameter: 5mm	DESCRIPTION
91-0630	15cm	Roeder Loop Applicator





# Laparoscopic Portfolio

---










## Accessories

Symmetry Surgical's laparoscopic line offers both reusable and single use devices for Access, Intervention and Closure.







# Accessories

## Access Accessories

CODE	LENGTH	DESCRIPTION	
52420-01	89 x 178 x 114 mm	Nathanson Style Liver Retractor, Small Size, 6.5mm diameter	
52420-02	114 x 216 x 203mm	Nathanson Style Liver Retractor, Medium Size, 6.5mm diameter	
52420-03	165 x 279 x 254mm	Nathanson Style Liver Retractor, Large Size, 6.5mm diameter	
52420-00		Nathanson Style Liver Retractor Kit, Includes One Small, One Medium, One Large	
90-6005	18cm	Hasson Retractor "S"	
90-9800		Fan Retractor, Curved, 5mm diameter	
68-1590		Symmetry® Fan Liver Retractor, 5mm diameter	
90-5054	10cm	Veress Needle with Luer Lock, Reusable, 2mm diameter	
90-5055	12.7cm	Veress Needle with Luer Lock, Reusable, 2mm diameter	

# Accessories

## Access Accessories

CODE	SHAFT Diameter: 5mm	DESCRIPTION	
90-5057	18cm	Veress Needles with Luer Lock, Reusable, 2mm diameter	
90-5058	8cm	Veress Needle with Luer Lock, Reusable, 2mm diameter	
90-5059	20cm	Veress Needle with Luer Lock, Reusable, 2mm diameter	
CD7212		Symmetry® CO2 Insufflation Needle, Single-Use, 14 Gauge, (10/box)	

# Accessories

## Access Accessories

CODE	QTY	DESCRIPTION	
50-5730		Bookwalter® Laparoscopic Support System *includes the following:	
50-4688	1	Bookwalter® 12" light-weight table post	
50-4554	1	Bookwalter® Table Post Coupling	
50-4805	1	Bookwalter® Universal Flex Arm	
50-4806	1	Bookwalter® Flex Arm Ring Adapter	
50-4807	1	Bookwalter® Flex Arm Laparoscopic Instrument Adapter	
50-4808	1	Bookwalter® Flex Arm Camera Lens Adapter	
50-4809	1	Bookwalter® Flex Arm Camera Housing Adapter	
68-1593	1	Fan Retractor, 5 mm, Articulating, Three Fingers	
68-1594	1	Fan Retractor, 10 mm, Articulating, Five Fingers	

# Accessories

## Access Accessories

NEW

CODE	WORKING LENGTH	PLATFORM WIDTH	DESCRIPTION
------	----------------	----------------	-------------

SHR7001	350mm	100mm	PretzelFlex™ 5mm
---------	-------	-------	------------------

SHR7101	(n/a)	(n/a)	Replacement sterilization sleeve
---------	-------	-------	----------------------------------

SHR7301	275mm	70mm	PretzelFlex™ 3mm
---------	-------	------	------------------

SHR7331	(n/a)	(n/a)	Replacement Sterilization Sleeve
---------	-------	-------	----------------------------------



PretzelFlex™ is proudly distributed by Symmetry Surgical® and is a product of Surgical Innovations Ltd. All trademarks and registered trademarks are the property of their respective owners.

# Accessories

## Access Accessories

NEW

CODE	DESCRIPTION
A-1000	FastClamp™ System Includes: 1 Table Clamp - TC-1000-01 1 Instrument Clamp - IC-1000-01 2 Vertical Bars Standard Length - BC-1030-01 & BC-1050-01
A-1001	Extended FastClamp™ System Includes: 1 Table Clamp - TC-1000-01 1 Extended Instrument Clamp - IC-1001-01 2 Vertical Bars - BC-1030-01 & BC-1090-01
TC-1000-01	FastClamp™ Table Clamp
IC-1000-01	FastClamp™ Instrument Clamp
IC-1001-01	FastClamp™ Extended Instrument Clamp
BC-1030-01	FastClamp™ 30mm Bar
BC-1050-01	FastClamp™ 50mm Bar
BC-1070-01	FastClamp™ 70mm Bar
BC-1090-01	FastClamp™ 90mm Bar
A-1000	FastClamp™ System



\*PretzelFlex™ retractor  
not included

PretzelFlex™ is proudly distributed by Symmetry Surgical® and is a product of Surgical Innovations Ltd. All trademarks and registered trademarks are the property of their respective owners.







# Accessories

## Intervention Accessories

CODE	LENGTH	DESCRIPTION	
90-5048	32cm	Ultra™ Puncture Needle with Aspiration, Reusable, 10 Gauge	
90-5049	32cm	Ultra™ Puncture Needle with Aspiration, Reusable, 12 Gauge	
90-5050	32cm	Ultra™ Puncture Needle with Aspiration, Reusable, 19 Gauge	
90-5050L	45cm		
90-5051	32cm	Ultra™ Puncture Needle with Aspiration, Reusable, 15 Gauge	
90-5051L	45cm		
90-5052	32cm	Ultra™ Puncture Needle with Aspiration, Reusable, 22 Gauge	
90-5052L	45cm		
90-5070	32cm	Symmetry® Cannula, Irrigating/Aspirating Dual Valve	
90-5070L	45cm		
90-5053	32cm	Ultra™ Palpation Probe	
90-5053L	45cm		
90-5600	32cm	Myoma Screw	
90-5605	45cm		
90-1076	32cm	Cholangiogram Guide Curved	
91-5035	39cm	Cholangiogram Clamp, 5mm diameter	

# Accessories

## Intervention Accessories

CODE	LENGTH	DESCRIPTION	
SG5		Secto® Dissector, Gauze, 5 mm diameter, 2.5 cm working length, single-use tip, 10/box	
		*For use with SGH5 handle	
SG10		Secto® Dissector, Gauze, 10 mm diameter, 3.4 cm working length, single-use tip, 10/box	
		*For use with SGH10 handle	
SGH5	33cm	Secto® Handle, Laparoscopic, Long, 5mm diameter	
SGH10	33cm	Secto® Handle, Laparoscopic, Long, 10mm diameter	









# Accessories

## Single-Use Accessories

CODE	DESCRIPTION	
PS500	JosNoe Medical Scope Warmer, Single-Use, 10/box	
TDSG5	Secto® 5mm One-Piece Single-Use Gauze Dissector, Length 35cm, 10/box	
TDSG10	Secto® 10mm One-Piece Single-Use Gauze Dissector, Length 35cm, 10/box	
DC5000	Fog Reduction Device, Sponge + Anti-Fog, 20/box	

# Accessories


## Single-Use Accessories

CODE	LENGTH Diameter: 4F	DESCRIPTION	
e2401-50	50 cm	Reddick® Scoop Tip Cholangiogram Catheter	
e2401-51	50 cm	Reddick® Scoop Tip Cholangiogram Stiffer Catheter	
e2401-52	50 cm	Reddick® Scoop Tip Cholangiogram Catheter with Needle Introducer ( <i>Needle Introducer: 13 gauge by 8 inches</i> )	
e2401-53	50 cm	Reddick® Scoop Tip Cholangiogram Stiffer Catheter with Needle Introducer ( <i>Needle Introducer: 13 gauge by 8 inches</i> )	
e2400-50	50 cm	Reddick® Cholangiogram Catheter	
e2400-51	50 cm	Reddick® Cholangiogram Stiffer Catheter	
e2400-52	50 cm	Reddick® Cholangiogram Catheter with Needle Introducer ( <i>Needle Introducer: 13 gauge by 8 inches</i> )	
e2400-53	50 cm	Reddick® Cholangiogram Stiffer Catheter with Needle Introducer ( <i>Needle Introducer: 13 gauge by 8 inches</i> )	

CODE	LENGTH Diameter: 5mm	DESCRIPTION	
e2504-01	30 cm	Reddick®-Saw Screw Retractor Kit 1 Reddick®-Saw Screw Retractor Kit and 113G Aspiration Needle	

# Accessories

## Miscellaneous

CODE	DESCRIPTION	
90-0000	Caps, Choli Flush Port, (10/pk)	

# Customer Support

---

Our highly trained customer support team is available to take your calls from 7:00 am to 5:30 pm (CST), Monday - Friday. This group is dedicated to providing a full range of support services, including:

- Order Placement
- Order Inquiries
- Delivery Information
- Product Information
- Pricing Inquiries

If we can be of assistance, please contact us using the methods below:

## PHONE

US & Canada: +1 800-251-3000

International: +1 615-883-9090

## EMAIL

customerservice@symmetrysurgical.com

internationalorders@symmetrysurgical.com

## WEB

symmetrysurgical.com

-Live Chat on [symmetrysurgical.com](https://www.symmetrysurgical.com):

Click the “Can We Help You” window in the bottom corner of any page on [symmetrysurgical.com](https://www.symmetrysurgical.com) to be connected directly to a member of our sales or customer service team. We can chat in real time and answer your questions about our company and surgical instrument portfolio.

(Available Monday through Friday, 7:30 am - 5:30 pm CST)

## MAILING ADDRESS

Symmetry Surgical  
Attn: Customer Service  
3034 Owen Drive  
Antioch, TN 37013



# Notes

---

